

Popis ispitnih metoda u fleksibilnom području akreditacije, status 2023-03-22
The valid list of the test methods in the flexible scope of accreditation, status 2023-03-22
Ovaj popis odnosi se na Prilog potvrdi o akreditaciji broj: 1073 koji je dostupan na www.akreditacija.hr
This list is related to Annex to Accreditation Certificate Number: 1073 which is available at www.akreditacija.hr
Izmjene su vidljive na kraju dokumenta / Changes visible at the end of the document.
FLEKSIBILNO PODRUČJE AKREDITACIJE / FLEXIBLE SCOPE OF ACCREDITATION

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
A1	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	03/01/2012	02/08/2019
		mg/l							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		Ag	0,001	Mo	0,004				
		Al	0,005	Na	0,005				
		As	0,005	Ni	0,005				
		B	0,013	P	0,050				
		Ba	0,001	Pb	0,005				
		Be	0,001	Sb	0,023				
		Ca	0,005	Se	0,011				
		Cd	0,001	Si	0,011				
		Co	0,001	Sn	0,015				
		Cr	0,001	Sr	0,001				
		Cu	0,001	Te	0,018				
A2	Eluat otpada Waste eluat	Fe	0,005	Ti	0,001				
		K	0,013	Tl	0,051				
		Li	0,001	U	0,010				
		Mg	0,003	V	0,012				
		Mn	0,001	Zn	0,005				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
A3	Mulj, sediment Sludge, sediment	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN 16170:2016 (EN 16170:2016)	13/09/2013	20/03/2023
		mg/kg							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		Al	0,83	Ni	0,44				
		As	0,83	P	8,33				
		Cd	0,16	Pb	2,00				
		Cr	0,17	Zn	0,47				
	Cu	0,87							
A4	A4-1 Otpad Waste	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	13/09/2013	20/03/2023
		mg/kg							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		As	2,62	Ni	0,44				
		Cd	0,16	Pb	2,00				
		Co	0,74	Sb	3,24				
		Cr	0,17	Tl	2,30				
		Cu	0,43	V	0,93				
		Gd	0,16	Zn	0,48				
		Mn	0,12						
	A4-2 Kruta oporabljena goriva Solid recovered fuels	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	05/05/2014	21/03/2023
		mg/kg							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		As	7,09	Mo	0,68				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Ba	0,60	Ni	0,78				
		Be	0,07	Pb	4,19				
		Cd	3,92	Sb	4,74				
		Co	0,24	Tl	3,39				
		Cr	0,26	V	1,13				
		Cu	5,34	Zn	5,27				
		Mn	0,05						
A5	Tlo Soil	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN 16170:2016 (EN 16170:2016)	23/07/2014	20/02/2023
		mg/kg							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		Al	0,40	Ni	0,44				
		Cd	0,16	Pb	2,0				
		Cr	0,17	Zn	0,47				
		Cu	0,87	As	2,62				
		Fe	0,93	Mo	1,97				
		Mn	0,03	Co	0,74				
A6	A6-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode, eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more)	Određivanje hlapljivih organskih spojeva metodom analize para iznad otopine (headspace) plinskom kromatografijom s masenom spektrometrijom (GC-MS/MS) Determination of volatile organic compounds in water by headspace method gas chromatography / mass spectrometry (GC-MS/MS)				GC-MS/MS	HRN EN ISO 10301:2002 (ISO 10301:1997; EN ISO 10301:1997)	23/07/2014	10/01/2023
		ug/L							
		Spoj /Compound		Granica kvantifikacija / Limit of quantification					
		kloroform / chloroform		0,09					
		tetraklorugljik / carbon tetrachloride		0,06					
		1,2-dikloretan / 1,2-dichloroethane		0,05					
		1,1,1-trikloretan / 1,1,1-trichloroethane		0,03					
		trikloreten (1,1,2-trikloreten) / trichloroethene		0,02					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification		Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)	(1,1,2- trichloroethene)					
		tetrakloreten (1,1,2,2-tetrakloreten) / tetrachloroethene (1,1,2,2- tetrachloroethene)	0,04				
		benzen / benzene	0,04				
		toluen / toluene	0,02				
		ksileni / xylene	0,04				
		etilbenzen / ethylbenzene	0,02				
		1,2,4- triklorbenzen / 1,2,4-trichlorobenzene	0,02				
		1,2,3-triklorbenzen / 1,2,3-trichlorobenzene	0,02				
		1,3,5-triklorbenzen / 1,3,5-trichlorobenzene	0,02				
		bromdiklormetan / bromodichloromethane	0,04				
		dibromklormetan / dibromochloromethane	0,05				
		bromoform / bromoform	0,04				
		diklormetan / dichloromethane	0,07				
		1,1-dikloreten / 1,1-dichloroethene	0,04				
		stiren / styrene	0,05				
		heksaklorbutadien / hexachlorobutadiene	0,04				
		vinil-klorid / vinil-chloride	0,03				
		1,2,4-trimetilbenzen / 1,2,4-trimethylbenzene	0,04				
		1,3,5-trimetilbenzen / 1,3,5-trimethylbenzene	0,02				
		tetrahidrofuran / tetrahydrofuran	0,02				
		tetrahidrotiofen / tetrahydrothiophen	0,04				
		1,1-dikloreten / 1,1-dichloroethane	0,03				
		cis-1,2-dikloreten / cis-1,2-dichloroethene	0,03				
		trans-1,2-dikloreten / trans-1,2-dichloroethene	0,05				
		1,1,2-trikloreten / 1,1,2-trichloroethane	0,04				
		1,1,2,2-tetrakloreten / 1,1,2,2-tetrachloroethane	0,04				
		heksakloreten / hexachloroethane	0,03				
		epiklorhidrin / epichlorohydrin	0,2				
		n-heksan / n-hexane	0,02				
		1,2-dibrometan / 1,2-dibromethane	0,05				
	A6-2 Tlo, otpad, mulj,	Određivanje BTEXa i heksaklorbutadiena metodom analize para iznad otopine (Headspace) plinskom kromatografijom sa		GC-MS/MS	Vlastita metoda In-house method	16/6/2014	14/6/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change												
	sediment Soil, waste, sludge, sediment	masenom spektrometrijom (GC-MS/MS) Determination of BTEX and hexachlorobutadiene method Headspace gas chromatography mass spectrometry (GC-MS/MS) mg/kg <table><tr><td>Spoj /Compound</td><td>Granica određivanja/ Limit of detection</td></tr><tr><td>benzen / benzene</td><td>0,0002</td></tr><tr><td>toluen / toluene</td><td>0,0001</td></tr><tr><td>ksileni / xylene</td><td>0,0002</td></tr><tr><td>etilbenzen / ethylbenzene</td><td>0,0001</td></tr><tr><td>heksaklorbutadien / hexachlorobutadiene</td><td>0,0002</td></tr></table>	Spoj /Compound	Granica određivanja/ Limit of detection	benzen / benzene	0,0002	toluen / toluene	0,0001	ksileni / xylene	0,0002	etilbenzen / ethylbenzene	0,0001	heksaklorbutadien / hexachlorobutadiene	0,0002		SOP-LEK-38,39/142 V. izdanje / edition (2022-05-20)		
Spoj /Compound	Granica određivanja/ Limit of detection																	
benzen / benzene	0,0002																	
toluen / toluene	0,0001																	
ksileni / xylene	0,0002																	
etilbenzen / ethylbenzene	0,0001																	
heksaklorbutadien / hexachlorobutadiene	0,0002																	
A7	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode, eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more) Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea	Određivanje pesticida u vodama metodom plinske kromatografije s masenom spektrometrijom Determination of pesticide residues in water – gas chromatographic method with mass spectrometry 0,0001 µg/l Cibutrin (irgarol), diklorvos / cibutrin (irgarol), dichlorvos 0,0006 µg/l Heptaklor i heptaklorepoksid / hepachlor and heptachlor epoxide 0,0007 µg/l Dikofol / Dicofol 0,001 µg/l Kloroneb, alfa-BHC (alfa-HCH), heksaklorbenzen (HCB), pentakloroanisol, beta-BHC (beta-HCH), gama-BHC (gama-HCH ili Lindan), delta-BHC (delta-HCH), klorbenzid, endosulfan eter, pentaklorotioanisol, aldrin, 4,4'-diklorobenzofenon, fenson, izodrin, cis-klordan, trans-klordan, o,p'-DDE (2,4-DDE), alfa-endosulfan (endosulfan I), klorfenson, p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), cis-nonaklor, trans-nonaklor, endrin aldehyd, 4,4'-metoksiklor olefin, endosulfan sulfat, p,p'-DDT (4,4-DDT), 2,4'-	GC-MS/MS	Vlastita metoda In-house method SOP-LEK-31-33,37/181 V. izdanje / edition (2020-05-28)	15/11/2017	31/05/2022												

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	water (transitional and coastal water; territorial sea)	metoksiklor (metoksiklor), tetradifon, mireks, pentaklorbenzen, endrin keton, etilan (pertan) / <i>Chloroneb, alph-BHC (alpha-HCH), hexachlorobenzene (HCB), pentachloroanisole, beta-BHC (beta-HCH), gamma-BHC (gamma-HCH or Lindane), chlorbenside, endosulfan ether; pentachlorothioanisole, aldrin, 4,4'-dichlorobenzophenone, fenison, isodrin, cis-chlordane, trans-chlordane, o,p'-DDE (2,4-DDE), alpha-endosulfan (endosulfan I), chlofenson (Ovex), p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), cis-nonachlor, trans-nonachlor, endrin aldehyde, 4,4'-metoxychlor olefin, endosulfan sulfate, p,p'-DDT (4,4-DDT), 2,4'-metoxychlor (metoxychlor), tetradifon, mirex, pentachlorobenzene, endrin ketone, ethylan (perthane)</i> 0,002 µg/l Bifenoks / Bifenox 0,008 µg/l Aklonifen / Aclonifen 0,01 µg/l Klorpirifos metil, pirimifos metil, klorpirifos etil (klorpirifos), kvinalfos, piridafention, EPN, fosalon, azinfos metil, azinfos etil, isazofos, pirimifos etil, difenilamin, kvintozen, diklofluanid, tolilfluanid, pendimetalin, nitrofen, nitalin, terbufos, fenklorfos, fention, bromofos etil, protiofos, etion, klortiofos, sulprofos, karbofenotion, leptofos, forat, fonofos, paration, triazofos, piperonil butoksid, alidoklor, etridiazol, pebulat, propaklor, cikloat, di-alat I i II (cis i trans), klomazon, terbutilazin, propizamid, pirimetanil, teflutrin, tri-alat, dimetaklor, acetoklor, transflutrin, metolaklor, linuron, triadimefon, MGK-264 I i II, difenamid, ciprodinil, penkonazol, triadimenol, procimidon, triflumizol, paklobutrazol, flutriafol, flutolanil, fludioksonil, pretilaklor, oksadiazon, miklobutanil, bupirimat, klorfenapir, heksazinon, tebukonazol, resmetrin I i II, fenpropatrin, tebufenpirad, fenotrin I i II (cis i trans), piriproksifen, akrinatri I i II, fenarimol, trans-permetrin, flukvinkonazol, ciflutrin I – IV, piridaben, cipermetrin I i II, flucitrinat I i II, fluvalinat I i II, deltametrin, tributilfosfat (TBP), trikloroetilfosfat (TCEP), trikloropropilfosfat (TCPP) / <i>Chlorpyrifos methyl, pirimiphos methyl, chlorpyrifos ethyl (chlorpyrifos), quinalphos, pyridaphenthion, EPN, phosalone, azinphos methyl, azinphos</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		ethyl, isazophos, pirimiphos ethyl, diphenylamine, quintozene, dichlofluanid, tolylfluanid, pendimethalin, nitrofen, nitralin, terbufos, fenchlorphos, fenthion, bromophos ethyl, prothiofos, ethion, chlorthiophos, sulprofos, carbophenothion, leptophos, phorate, fonofos, parathion, triazophos, piperonyl butoxide, allidochlor, etridiazole, pebulate, propachlor, cycloate, diallate I i II (cis i trans), clomazone, terbuthylazine, propyzamide, pyrimethanil, tefluthrin, triallate, dimethachlor, acetochlor, transfluthrin, metolachlor, linuron, triadimefon, MGK-264 I i II, diphenamid, cyprodinil, penconazole, triadimenol, procymidone, triflumizole, paclobutrazol, flutriafol, flutolanil, fludioxonil, pretilachlor, oxadiazon, myclobutanil, bupirimate, chlorfenapyr, hexazinone, tebuconazole, resmethrin I i II, fenpropathrin, tebufenpyrad, phenothrin I i II (cis i trans), pyriproxyfen, acrinathrin I i II, fenarimol, transpermethrin, fluquinconazole, cyfluthrin I – IV, pyridaben, cypermethrin I i II, flucythrinate I i II, fluvalinate I i II, deltamethrin, tributyl phosphate (TBP), trichloroethyl phosphate (TCEP), tris chloropropyl phosphate (TCPP) 0,02 µg/l Prometrin / Prometryn 0,021 µg/l Folpet / Folpet 0,025 µg/l Mankozeb / Mancozeb Ukupno / In total 147				
A8	A8-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode, eluat otpada,	Određivanje odabranih sredstava za zaštitu bilja u vodi metodom tekućinske kromatografije s masenom spektrometrijom <i>Determination of selected plant treatment agents in water – liquid chromatographic method with mass spectrometry</i> 0,002 µg/l Bromacil, trifluralin, kinoksifen, cianazin, alaklor, simazin, diuron, heksazinon, dimetenamid, metazaklor, metolaklor, sekbumeton, N,N-dietil-m-toluamid, sebutilazin, mevinfos, atrazin desetil, diazinon, terbutilazin,	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31-33,37/183 V. izdanje / <i>edition</i> (2020-05-29)	14/12/2017	01/06/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	morska voda (prijelazne i priobalne vode, teritorijalno more) Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)	terbutrin, malation, klorfenvinfos, diklofenak, propazin, klorpirifos, atrazin, izoproturon, dimetoat, klorbromuron, metomitron, metoksuron, terbumeton, metalaksil, klorotoluron, fluometuron I, fluometuron II, imidaklopid, linuron, monuron, neburon, triadimefon, buturon, desmetil klorotoluron, metribuzin, azoksistrobin, ciprokonazol, propikonazol, acetamiprid / <i>Bromacil, trifluralin, quinoxifen, cyanazine, alachlor, simazine, diuron, hexazinon, dimethenamide, metazachlor, metolachlor, secbumeton, N,N- diethyl-m-toluamide, sebutylazine, mevinphos, atrazine desethyl, diazinon, terbuthylazine, terbutryn, malathion, chlorfenvinphos, dichlofenac, propazine, chlorpyrifos (Dursban), atrazine, isoproturon, dimethoate, chlorbromuron, metomitron, metoxuron, terbumeton, metalaxyl, chlorotoluron, fluometuron I, fluometuron II, imidaklopid, linuron, monuron, neburon, triadimefon, buturon, chlorotoluron-desmethyl, metribuzin, azoxystrobin, cyproconazole, propiconazole, acetamiprid</i> 0,006 µg/l Metobromuron / metobromuron 0,008 µg/l Atrazin-desizopropil /Atrazine-desisopropyl 0,012 µg/l Monolinuron / monolinuron 0,02 µg/l Piretrin, ometoat, flazasulfuron / pyrethrin, omethoat, flazasulfuron Ukupno / In total 52				
	A8-2 Voda za ljudsku potrošnju, prirodne mineralne, prirodne izvorske, stolne vode, podzemne, površinske vode, procjedne i otpadne	Određivanje kiselih pesticida u vodi metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of acidic pesticides in water – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,002 µg/L	LC-MS/MS	<i>Vlastita metoda In-house method SOP-LEK-31- 33,37/182 III. izdanje / edition (2020-02-03) modificirana</i>	20/04/2016	06/06/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	vode i eluati otpada <i>Water for human consumption, natural spring water, natural mineral water, table water, ground water, surface water, leachate water, waste water and waste eluate</i>	2,4-DP, bentazon, dikamba, MCPB, 2,4-DB, MCPA, fenoprop (Silveks), 2,4-D, ioksinil / <i>2,4-DP, bentazon, dicamba, MCPB, 2,4-DB, MCPA, fenoprop (Silvex), 2,4-D, ioxynil</i>		<i>/modified DIN 38407-35:2010</i>		
A9	A9-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,005 mg/kg 1-naftalenacetamid, 2,3,5-Trimetakarb, 2,4,5-T, 2,4-D, 2,4-dimetilanilin, 3,4,5-Trimetakarb, 3-Decen-2-on, 3-hidroksikarbofuran, 4-klorofenoksioktenu kiselina, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Akrinatri, Alaklor, Alanikarb, Aldikarb, Aldikarb-sulfoksid, Aldikarb-sulfon (Aldoksikarb), Alidoklor, Ametotradin, Ametrin, Amidosulfuron, Aminokarb, Amisulbrom, Ancimidol, Anilazin, Atraton, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a, Azadirachtin, Azakonazol, Azametifos, Azinfos-etil, Azinfos-metil, Aziprotrin, Azoksistrobin, Benalaksil, Bendiokarb, Benfurakarb, Benodanil, Benomil, Bensulfuron-metil, Bensulid, Bentazon, Bentiavalikarb-izopropil, Bentivalikarb, Benziladenin, Benzoilprop-etil, Benzoksimat, Benzovindflupir, Bifenazat-diazen, Bifentrin, Biksafen, Bispiribak-natrij, Bitertanol, Boskalid, Bromacil, Bromofos-etil, Bromoksinil, Bromkonazol, Bupirimat, Buprofezin, Butafenacil, Butaklor, Butilat, Butoksikarboks, Butralin, Buturon, Centrimonium klorid, Cianazin,	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	21/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Cianfenfos, Cianofos, Ciantraniliprol, Ciazofamid, Ciflufenamid, Ciflumetofen, Cihalofof-butil, Ciklanilid, Cikloat, Cikloksidim, Cikluron, Cimiazol, Cimoksanil, Cinidon-etil, Ciprodinil, Ciprokonazol, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfoksid, Demeton-S-metil-sulfon, Desmedifam, Diafentiuron, Diallate, Diazinon, Dietofenkarb, Difenamid, Difenokonazole, Difenoksuron, Difetialon, Diflubenzuron, Diflufenikan, Diklobutrazol, Diklofenak, Diklofention, Diklofop, Diklofop-metil, Diklorprop, Diklorvos, Dikrotofos, Dimefuron, Dimetaklor, Dimetenamid, Dimetilan, Dimetoat, Dimetomorf (E,Z), Dimoksisstrobin, Dinikonazol, Dinoseb, Dinotefuran, Dinoterb, Dioksakarb, Dioksation, Dipropetrin, Disulfoton-sulfoksid, Disulfoton-sulfon, Diuron (DCMU), DMST, DNOC, Dodemorf, Dodin, Edifenfos, Enamektin-benzoat, Endosulfan sulfat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etiofenkarb-sulfon, Etion, Etrimol, Etofenproks, Etofumezat, Etoksakonazol, Etoprofos, Etrimfos, Famoksadon, Fampur, Fenamidon, Fenamifos, Fenamifos-sulfoksid, Fenamifos-sulfon, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenmedifam, Fenobukarb, Fenoksaprop-etil, Fenoksikarb, Fenoprop, Fenotiokarb, Fenotrin, Fenpikoksamid, Fenpirazamin, Fenpiroksimat, Fenpropatrin, Fenpropidin, Fenpropimorf, Fensulfotion, Fensulfotion-sulfon, Fention, Fention-okson, Fention-okson-sulfoksid, Fention-okson-sulfon, Fention-sulfoksid, Fention-sulfon, Fentoat, Fenuron, Fipronil, Fipronil-sulfid, Fipronil-sulfon, Flamprop-izopropil, Flamprop-metil, Flamprop-M-izopropil, Flazasulfuron, Flazuron, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazifop-P-butyl, Fluazinam, Flucipirim, Flucitrat, Fludioksonil, Flufenacet, Flufenoksuron, Fluksapiroksad, Flukvinkonazol, Flumetralin, Fluoksastrobin, Fluometuron, Fluopikolid, Fluopiram, Fluridon, Flurokloridon, Fluroksipir, Flusilazol, Flutiacet-metil, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Foksim, Fonofos, Forat, Forat-okson, Forat-okson-sulfon, Forat-sulfoksid, Forat-sulfon, Forklorfenuron, Formetanat hidroklorid, Fosalon, Fosfamidon, Fosfolan, Fosmet, Fosmet-okson, Fostiazat, Fuberidazol, Furalaksil, Furatiokarb, Furmecikloks, Haloksifop, Haloksifop-2-etoksietil, Haloksifop-metil, Haloksifop-P-metil, Heksaflumuron, Heksakonazol, Heksazinon, Heksitiazoks, Heptenofos, Hidrametilnon, Imazakvin, Imazalil, Imazamoks, Imazapir, Imazetapir, Imibenkonazol, Imidaklopid, Indoksakarb, Ioksonil, Ipkonazol, Iprobenfos, Iprovalikarb, Ivermektin B1a, Izazofos, Izofenfos-metil, Izofetamid,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Izokarbamid, Izokarbofos, Izoksaben, Izoksadifen-etil, Izoksaflutol, Izoksaflutol-dikenotril, Izoksation, Izometiozin, Izopirazam, Izoprokarb, Izopropalin, Izoprotiolan, Izoproturon, Jodosulfuron-metil, Kadusafos, Karanjin, Karbaril (NAC), Karbendazim, Karbetamid, Karbofenotion, Karbofuran, Karboksini, Karboksini sulfoksidi, Karbosulfan, Karfentrazon-etil, Kinosulfuron, Kletodim, Klimbazol, Klodinafop (slobodna kiselina), Klodinafop-propargil, Klorfentezin, Klokventocet-meksil, Klorazon, Klorbromuron, Klorimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Kloridazon-desfenil, Kloroantraniliprol, Kloroksuron, Klorotoluron, Klorpirifos, Klorpirifos-metil, Klorprofam, Klorosulfuron, Klorotiamid, Klorotianidin, Kresoksime-metil, Klorfumat, Klorotoksifos, Kumafos, Kumatetralil, Klorinalfos, Klorinklorak, Klorinmerak, Klorinoklamin, Klorinoksifen, Klorizalofop (slobodna kiselina), Klorizalofop-etil, Klorizalofop-metil, Klorizalofop-P, Klorizalofop-P-etil, Klorizalofop-P-tefuril, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), MCPB, Mefenacet, Mefenpir-dietil, Mefentriklonazol, Mefosfolan, Mekarbam, Meksakarbat, Mepanipirim, Mepronil, Meptildinokap, Metabenzotiazuron, Metaflumizon, Metakrifos, Metalaksil i Metalaksil-M, Metamidofos, Metamitron, Metazaklor, Metazaklor ESA, Metazaklor OA, Metidation, Metiokarb, Metiokarb-sulfoksidi, Metiokarb-sulfon, Metkonazol, Metobromuron, Metoksifenozi, Metoksuron, Metolaklor, Metolcarb (MTMC), Metomil, Metominostrobin (E, Z), Metoprotin, Metrafenon, Metribuzin, Metsulfuron-metil, Mevinfos, Miklobutanil, Milbemektin A3, Molinat, Monokrotofos, Monolinuron, Monuron, N-(2,4-dimetilfenil)-N-metilformamidin, N-(2,4-Dimetilfenil)formamid, N,N-dietil-m-toluamid, Napropamid, N-desmetil-propamokarb, Neburon, Nereistoksin oksalat, Nikosulfuron, Nitenpiram, Norflurazon, Novaluron, Nuarimol, Ofurak, Oksadiargil, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Oksikarboksini, Oksimatin, Ometoat, Orizalin, Paklobutrazol, Paraokson-etil, Paraokson-metil, Pebulat, Pencikuron, Pendimetalin, Penflufen, Penkonazol, Pentiopirad, Petoksamid, Picaridin (Icaridin), Pikoksistrobin, Pikolinafen, Pimetrozin, Piperofos, Piperonil-butoksidi, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Pirazokson, Piretrini, Piridaben, Piridafention, Piridafol, Piridilil, Piridat, Pirifeonks, Pirimetanil, Pirimifos-metil, Pirimikarb, Pirimikarb-desmetil, Piriofenon, Piriproksifen, Potasan, Praletrin, Pretilaklor, Profam, Profenofos, Profoksidi, Prokloraz,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Prokloraz met. BTS 44595, Prokloraz met. BTS 44596, Prokvinazid, Promekarb, Prometon, Prometrin, Propaklor, Propakvizafop, Propamokarb, Propamokarb-N-oksidi, Propanil, Propargit, Propazin, Propetamfos, Propikonazol, Propizamid, Propoksur, Prosulfokarb, Prosulfuron, Protiokonazol-destio, Rabenzazol, Resmetrin, Rimsulfuron, Rotenon, Saflufenacil, Sebutilazin, Sekbumeton, Setoksidim, Siduron, Silafluofen, Siltiofam, Simazin, Simetrin, Spinetoram, Spinosad (suma Spinosin A i Spinsin D), Spirodiklofen, Spiroksamin, Spiromezifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glukozid, Spirotetramat-keto-hidroksi, Spirotetramat-mono-hidroksi, Sulfoksaflo, Sulfosulfuron, Sulfotep, Sulkotri, Sulprofos, TCMTB, Tebufenozid, Tebufenpirad, Tebukonazol, Tebupirimifos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrion, Temefos, Tepraloksidim, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetraconazole, Tetraklorinfos (CVMP), Tetrametrin, Tiabendazol, Tiakloprid, Tiametoksam, Tiazuron, Tifensulfuron-metil, Tiobenkarb, Tiociklam, Tiociklam hidrogen oksalat, Tiodikarb, Tiofanat-metil, Tiofanoks, Tiofanoks-sulfoksidi, Tiofanoks-sulfoni, Tiometon, Tionazin, Tiram, Tolfenpirad, Tolilfluorid, Tolklofos-metil, Tralometrin, Triadimefon, Triadimenol, Trialat, Triazamat, Triazofos, Tribenuron-metil, Triciklazol, Trietazin, Trifloksistobin, Triflumizol, Triflumizol metabolit FM-6-1, Triflururon, Triflururon-sulfoni, Triflorin, Triklorpir, Triklorfon, Trisulfuron, Tritikonazol, Tritosulfuron, Valifenalat, Vamidotio, Vernolat, Vinklozolin, XMC, Zoksamid / <i>1-naphthaleneacetamide, 2,3,5-Trimethacarb, 2,4,5-T, 2,4-D, 2,4-dimethylaniline, 3,4,5-Trimethacarb, 3-Decen-2-one, 3-Hydroxycarbofuran, 4-Chlorophenoxyacetic acid, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Acrinathrin, Alachlor, Alanycarb, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Amidosulfuron, Aminocarb, Amisulbrom, Ancymidol, Anilazine, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a, Azaconazole, Azadirachtin, Azamethiphos, Azinphos-ethyl, Azinphos-methyl, Aziprotryne, Azoxystrobin, Benalaxyl, Bendiocarb, Benfuracarb, Benodanil, Benomyl, Bensulfuron-methyl, Bensulide, Bentazone, Benthiavalicarb, Benthiavalicarb-isopropyl, Benzovindflupyr, Benzoximate, Benzoylprop-ethyl, Benzyladenine, Bifenazate-diazene, Bifenthrin, Bispyribac-sodium, Bitertanol, Bixafen,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Boscalid, Bromacil, Bromophos-ethyl, Bromoxynil, Bromuconazole, Bupirimate, Buprofezin, Butachlor, Butafenacil, Butoxycarboxim, Butralin, Buturon, Butylate, Cadusafos, Carbaryl (NAC), Carbendazim, Carbetamide, Carbofuran, Carbophenothion, Carbosulfan, Carboxin, Carboxin sulfoxide, Carfentrazone-ethyl, Centrimonium chloride, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chloridazon-desphenyl, Chlorotoluron, Chloroxuron, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorsulfuron, Chlorthiamid, Cinidon-ethyl, Cinosulfuron, Clethodim, Climbazole, Clodinafop (free acid), Clodinafop-propargyl, Clofentezine, Clomazone, Cloquintocet-mexyl, Clothianidin, Coumaphos, Coumatetralyl, Crotoxyphos, Crufomate, Cyanazine, Cyanofenphos, Cyanophos, Cyantraniliprole, Cyazofamid, Cyclanilide, Cycloate, Cycloxydim, Cycluron, Cyflufenamid, Cyflumetofen, Cyhalofop-butyl, Cymiazole, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Di-allate, Diazinon, Dichlofenac, Dichlofenthion, Dichlorprop, Dichlorvos, Diclobutrazol, Diclofop, Diclofop-methyl, Dicrotophos, Diethofencarb, Difenoconazole, Difenoxuron, Difethialone, Diflubenzuron, Diflufenican, Dimefuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimetilan, Dimoxystrobin, Diniconazole, Dinoseb, Dinotefuran, Dinoterb, Dioxacarb, Dioxathion, Diphenamid, Dipropetryn, Disulfoton-sulfone, Disulfoton-sulfoxide, Diuron (DCMU), DMST, DNOC, Dodemorph, Dodine, Edifenphos, Emamectin-benzoate, Endosulfan-sulfate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethiofencarb-sulfone, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrinfos, Famoxadone, Famphur, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenothiocarb, Fenoxaprop-ethyl, Fenoxycarb, Fenpicoxamid, Fenpropathrin, Fenpropidin, Fenpropimorph, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fensulfothion-sulfone, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfide, Fipronil-sulfone, Flamprop-isopropyl, Flamprop-methyl, Flamprop-M-isopropyl, Flazasulfuron, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluacrypyrim, Fluazifop, Fluazifop-P-butyl,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Fluazinam, Fluazuron, Flucythrinate, Fludioxonil, Flufenacet, Flufenoxuron, Flumetralin, Fluometuron, Fluopicolide, Fluopyram, Fluoxastrobin, Fluquinconazole, Fluridone, Flurochloridone, Fluroxypyr, Flusilazole, Fluthiacet-methyl, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Furmecyclox, Haloxyfop, Haloxyfop-2-ethoxyethyl, Haloxyfop-methyl, Haloxyfop-P-methyl, Heptenophos, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazamox, Imazapyr, Imazaquin, Imazethapyr, Imibenconazole, Imidacloprid, Indoxacarb, Iodosulfuron-methyl, Ioxynil, Ipconazole, Iprobenfos, Iprovalicarb, Isazofos, Isocarbamid, Isocarbophos, Isofenphos-methyl, Isofetamid, Isomethiozin, Isoprocab, Isopropalin, Isoprothiolane, Isoproturon, Isoprazam, Isoxaben, Isoxadifen-ethyl, Isoxaflutole, Isoxaflutole dikenotriple, Isoxathion, Ivermectin B1a, Karanjin, Kresoxim-methyl, Lenacil, Linuron, Lufenuron, Malaoxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), MCPB, Mecarbam, Mefenacet, Mefenpyr-diethyl, Mefentrifluconazole, Mepanipyrim, Mephosfolan, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metazachlor ESA, Metazachlor OA, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprottryne, Methoxyfenozide, Metobromuron, Metolachlor, Metolcarb (MTMC), Metominostrobin, Metoxuron, Metrafenone, Metribuzin, Metsulfuron-methyl, Mevinphos, Mexacarbate, Milbemectin A3, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N-(2,4-Dimethylphenyl)formamide, N-(2,4-dimethylphenyl)-N-methylformamidine, N,N-diethyl-m-toluamide, Napropamide, N-desmethyl-propamocarb, Neburon, Nereistoxin oxalate, Nicosulfuron, Nitenpyram, Norflurazon, Novaluron, Nuarimol, Ofurace, Omethoate, Oryzalin, Oxadiargyl, Oxadiazon, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapiprolin, Oxycarboxin, Oxydemeton-methyl, Oxytrine, Paclobutrazol, Paraoxon-ethyl, Paraoxon-methyl, Pebulate, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthiopyrad, Pethoxamid, Phenmedipham, Phenothrin, Phenthoate, Phorate, Phorate-oxon, Phorate-oxon-sulfone, Phorate-sulfone, Phorate-sulfoxide, Phosalone, Phosfolan, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picaridin (Icaridin), Picolinafen, Picoxystrobin, Piperonyl-butoxide, Piperophos,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Potasan, Prallethrin, Pretilachlor, Prochloraz, Prochloraz met. BTS 44595, Prochloraz met. BTS 44596, Profenofos, Profoxydim, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propamocarb-N--oxide, Propanil, Propaquizafop, Propargite, Propazine, Propetamphos, Propham, Propiconazole, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prosulfuron, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrazoxon, Pyrethrins, Pyridaben, Pyridafol, Pyridalyl, Pyridaphenthion, Pyridate, Pyrifenox, Pyrimethanil, Pyriofenone, Pyriproxifen, Quinalphos, Quinclorac, Quinmerac, Quinoclamine, Quinoxifen, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rabenzazole, Resmethrin, Rimsulfuron, Rotenone, Saflufenacil, Sebuthylazine, Secbumeton, Sethoxydim, Siduron, Silafluofen, Silthiofam, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spiromesifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glucoside, Spirotetramat-keto-hydroxy, Spirotetramat-mono-hydroxy, Spiroxamine, Sulcotrione, Sulfosulfuron, Sulfotep, Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebupirimfos, Tebutam, Tebuthiuron, Teflubenzuron, Tembotrione, Temephos, Tepraloxymid, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiamethoxam, Thidiazuron, Thifensulfuron-methyl, Thiobencarb, Thiocyclam hydrogene oxalate, Thiodicarb, Thiofanox, Thiofanox-sulfone, Thiofanox-sulfoxide, Thiometon, Thionazin, Thiophanate-methyl, Thiram, Tolclofos-methyl, Tolfenpyrad, Tolyfluanid, Tralomethrin, Triadimefon, Triadimenol, Tri-allate, Triasulfuron, Triazamate, Triazophos, Tribenuron-methyl, Trichlorfon, Triclopyr, Tricyclazole, Trietazine, Trifloxystrobin, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triflusulfuron-methyl, Triforine, Triticonazole, Tritosulfuron, Valifenalate, Vamidothion, Vernolate, Vinclozolin, XMC, Zoxamide Ukupno / In total 553 Određivanje kiselih pesticida nakon alkalne hidrolize u hrani				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-DB; 2,4,5-T; Bentazone				
	A9-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,005 mg/kg 1-naftalenacetamid, 2,3,5-Trimetakarb, 2,4,5-T, 2,4-D, 2,4-DB, 2,4-dimetilanilin, 3,4,5-Trimetakarb, 3-Decen-2-on, 3-hidroksikarbofuran, 4-klorofenoksiocetna kiselina, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Alaklor, Alanikarb, Aldikarb, Aldikarb-sulfoksid, Aldikarb-sulfon (Aldoksikarb), Alidoklor, Ametoktradin, Ametrin, Amidosulfuron, Aminokarb, Amisulbrom, Ancimidol, Anilazin, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a, Azadirahtin, Azakonazol, Azametifos, Azinfos-etil, Azinfos-metil, Aziprotrin, Azoksistrobin, Benalaksil, Benfurakarb, Benodanil, Benomil, Bensulfuron-metil, Bensulid, Bentazon, Bentiavalikarb-izopropil, Bentalikarb, Benziladenin, Benzilprop-etil, Benzoksimat, Benzovindflupir, Bifenazat, Bifenazat-diazin, Bifentrin, Biksafen, Bispiribak-natrij, Bitertanol, Boskalid, Bromacil, Bromofos-etil, Bromoksinil, Bromokonazol, Bupirimat, Buprofezin, Butafenacil, Butaklor, Butilat, Butoksikarboksim, Butralin, Buturon, Centrimonium klorid, Cianazin, Cianfenfos, Cianofos, Ciantraniliprol, Ciazofamid, Ciflufenamid, Ciflumenofen, Cihalofof-butil, Ciklanilid, Cikloat, Cikloksidim, Cikluron, Cimiazol, Cimoksanil, Cinidon-	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	21/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		etil, Ciprodinil, Ciprokonazol, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfoksid, Demeton-S-metil-sulfon, Desmedifam, Diafentiuron, Diallylate, Diazinon, Dietofenkarb, Difenokonazole, Difenoksuron, Difetialon, Diflubenzuron, Diflufenikan, Diklobutrazol, Diklofenak, Diklofention, Diklofop, Diklofop-metil, Diklorprop, Diklorvos, Dikrotofos, Dimefuron, Dimetaklor, Dimetenamid, Dimetilan, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinoseb, Dinotefuran, Dinoterb, Dioksation, Dipropettrin, Disulfoton-sulfoksid, Disulfoton-sulfon, Diuron (DCMU), DMST, DNOC, Dodemorf, Dodin, Edifenfos, Enamektin-benzoat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etiofenkarb-sulfon, Etion, Etirimol, Etofenproks, Etofumezat, Etoksakonazol, Etoprofos, Etrimfos, Famoksadon, Fampur, Fenamidon, Fenamifos, Fenamifos-sulfoksid, Fenamifos-sulfon, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenobukarb, Fenoksaprop-etil, Fenoksikarb, Fenoprop, Fenotiokarb, Fenotrin, Fenpikoksamid, Fenpirazamin, Fenpiroksimat, Fenpropatrin, Fenpropidin, Fensulfotion, Fensulfotion-sulfon, Fention, Fention-okson, Fention-okson-sulfoksid, Fention-okson-sulfon, Fention-sulfoksid, Fention-sulfon, Fentoat, Fenuron, Fipronil, Fipronil-sulfid, Fipronil-sulfon, Flamprop-izopropil, Flamprop-metil, Flamprop-M-izopropil, Flazasulfuron, Flazuron, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazifop-P-butil, Fluazinam, Flucipirim, Flucitrinat, Fludioksonil, Flufenacet, Flufenoksuron, Fluksapiroksad, Flukvinkonazol, Flumetralin, Fluoksastrobin, Fluometuron, Fluopikolid, Fluopiram, Fluridon, Flurokloridon, Fluroksipir, Flusilazol, Flutiacet-metil, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Foksim, Fonofos, Forat, Forat-okson, Forat-okson-sulfon, Forat-sulfoksid, Forat-sulfon, Forklorfenuron, Formetanat hidroklorid, Fosalon, Fosfamidon, Fosfolan, Fosmet, Fosmet-okson, Fostiazat, Fuberidazol, Furalaksil, Furatiokarb, Furmecikloks, Haloksifop, Haloksifop-2-etoksietil, Haloksifop-metil, Haloksifop-P-metil, Heksaflumuron, Heksakonazol, Heksazinon, Heksitiazoks, Heptenofos, Hidrametilnon, Imazakvin, Imazalil, Imazamoks, Imazapir, Imazetapir, Imibenkonazol, Imidakloprid, Indoksakarb, Ioksonil, Ipkonazol, Iprobenfos, Iprovalikarb, Ivermektin B1a, Izazofos, Izofenfos-metil, Izofetamid, Izokarbamid, Izokarbofos, Izoksaben, Izoksadifen-etil, Izoksaflutol, Izoksaflutol dikenotril, Izoksation, Izometiozin, Izopirazam, Izoprokarb, Izopropalin, Izoprotiolan, Izoproturon, Jodosulfuron-metil, Kadusafos,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Karanjin, Karbendazim, Karbetamid, Karbofenotion, Karbofuran, Karboksini, Karboksini sulfoksidi, Karbosulfan, Karfentrazon-etil, Kinosulfuron, Kletodim, Klimbazol, Klodinafop (slobodna kiselina), Klodinafop-propargil, Klofentezin, Klokintocet-meksil, Klomazon, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Kloridazon-desfenil, Kloroantraniliprol, Kloroksuron, Klorotoluron, Klorpirifos, Klorpirifos-metil, Klorprofam, Klorosulfuron, Klortiamid, Klotianidin, Kresoksime-metil, Krofumat, Krotoksifos, Kumafos, Kumatetralil, Kvinalfos, Kvinklorak, Kvinmerak, Kvinoklamin, Kvinoksifen, Kvizalofop (slobodna kiselina), Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), MCPB, Mefenacet, Mefenpir-dietil, Mefentriklonazol, Mefosfolan, Mekarbam, Meksakarbat, Mepanipirim, Mepronil, Meptildinokap, Metabenzotiazuron, Metaflumizon, Metakrifos, Metalaksil i Metalaksil-M, Metamidofos, Metamitron, Metazaklor, Metazaklor ESA, Metazaklor OA, Metidation, Metkonazol, Metobromuron, Metoksifenozi, Metoksuron, Metolaklor, Metolcarb (MTMC), Metomil, Metominostrobin (E, Z), Metoprotin, Metrafenon, Metribuzin, Metsulfuron-metil, Mevinfos, Miklobutanil, Milbemektin A3, Molinat, Monokrotofos, Monolinuron, Monuron, N-(2,4-dimetilfeil)-N-metilformamidin, N-(2,4-Dimetilfenil)formamid, N,N-dietil-m-toluamid, Napropamid, N-desmetil-propamokarb, Neburon, Nereistoksin oksalat, Nikosulfuron, Nikotin, Nitenpiram, Norflurazon, Novaluron, Nuairimol, Ofurak, Oksadiargil, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Oksikarboksini, Oksimatrini, Ometoat, Orizalin, Paklobutrazol, Paraokson-etil, Paraokson-metil, Pebulat, Pencikuron, Pendimetalin, Penflufen, Penkonazol, Pentiopirad, Petoksamid, Picaridin (Icaridin), Pikoksistrobin, Pikolinafen, Pimetrozin, Piperofos, Piperonil-butoksid, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Pirazokson, Piretrini, Piridaben, Piridafention, Piridafol, Piridalil, Pirifeonks, Pirimetanil, Pirimifos-metil, Pirimikarb, Pirimikarb-desmetil, Piriofenon, Piriproksifen, Potasan, Pralettrin, Pretlaklor, Profam, Profenofos, Profoksimid, Prokloraz, Prokloraz met. BTS 44595, Prokloraz met. BTS 44596, Prokvinazid, Promekarb, Prometon, Prometrin, Propaklor, Propakvizafop, Propamokarb, Propamokarb-N-oksidi, Propanil, Propargit, Propazin, Propetamfos, Propikonazol, Propizamid, Propoksuri, Prosulfokarb, Prosulfuron,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Protiokonazol-destio, Rabenzazol, Resmetrin, Rimsulfuron, Rotenon, Saflufenacil, Sebutilazin, Sekbumeton, Setoksidim, Siduron, Silafluofen, Siltiofam, Simazin, Simetrin, Spinetoram, Spinosad (suma Spinosin A i Spinsin D), Spirodiklofen, Spiroksamin, Spiromezifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glukozid, Spirotetramat-keto-hidroksi, Spirotetramat-mono-hidroksi, Sulfoksafloor, Sulfosulfuron, Sulfotep, Sulkotripon, Sulprofos, TCMTB, Tebufenozid, Tebufenpirad, Tebukonazol, Tebupirimifos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrion, Temefos, Tepaloksidim, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetraconazole, Tetraklorvinfos (CVMP), Tetrametrin, Tiabendazol, Tiaklopid, Tiametoksam, Tiazuron, Tifensulfuron-metil, Tiobenkarb, Tiociklam hidrogen oksalat, Tiodikarb, Tiofanat-metil, Tiofanoks, Tiofanoks-sulfoksid, Tiofanoks-sulfon, Tiometon, Tionazin, Tolfenpirad, Tolilfluamid, Tolklofos-metil, Triadimefon, Triadimenol, Trialat, Triazamat, Triazofos, Tribenuron-metil, Triciklazol, Trietazin, Trifloksistrobin, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Triflusulfuron-metil, Triflorin, Triklorpir, Triklorfon, Trisulfuron, Tritikonazol, Tritosulfuron, Valifenalat, Vamidotion, Vernolat, XMC, Zoksamid / 1-naphthaleneacetamide, 2,3,5-Trimethacarb, 2,4,5-T, 2,4-D, 2,4-DB, 2,4-dimethylaniline, 3,4,5-Trimethacarb, 3-Decen-2-one, 3-Hydroxycarbofuran, 4-Chlorophenoxyacetic acid, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Alachlor, Alanycarb, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Amidosulfuron, Aminocarb, Amisulbrom, Ancymidol, Anilazine, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a, Azaconazole, Azadirachtin, Azamethiphos, Azinphos-ethyl, Azinphos-methyl, Aziprotryne, Azoxystrobin, Benalaxyl, Benfuracarb, Benodanil, Benomyl, Bensulfuron-methyl, Bensulide, Bentazone, Benthiavalicarb, Benthiavalicarb-isopropyl, Benzovindflupyr, Benzoximate, Benzoylprop-ethyl, Benzyladenine, Bifenazate, Bifenazate-diazene, Bifenthrin, Bispyribac-sodium, Bitertanol, Bixafen, Boscalid, Bromacil, Bromophos-ethyl, Bromoxynil, Bromuconazole, Bupirimate, Buprofezin, Butachlor, Butafenacil, Butoxycarboxim, Butralin, Buturon, Butylate, Cadusafos, Carbendazim, Carbetamide, Carbofuran, Carbophenothion, Carbosulfan, Carboxin, Carboxin sulfoxide, Carfentrazone-ethyl, Centrimonium chloride,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chloridazon-desphenyl, Chlorotoluron, Chloroxuron, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorsulfuron, Chlorthiamid, Cinidon-ethyl, Cinosulfuron, Clethodim, Climbazole, Clodinafop (free acid), Clodinafop-propargyl, Clofentezine, Clomazone, Cloquintocet-mexyl, Clothianidin, Coumaphos, Coumatetralyl, Crotoxyphos, Crufomate, Cyanazine, Cyanofenphos, Cyanophos, Cyantraniliprole, Cyazofamid, Cyclanilide, Cycloate, Cycloxydim, Cycluron, Cyflufenamid, Cyflumetofen, Cyhalofop-butyl, Cymiazole, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Di-allate, Diazinon, Dichlofenac, Dichlofenthion, Dichlorprop, Dichlorvos, Diclobutrazol, Diclofop, Diclofop-methyl, Dicrotophos, Diethofencarb, Difenoconazole, Difenoxuron, Difethialone, Diflubenzuron, Diflufenican, Dimefuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimetilan, Dimoxystrobin, Diniconazole, Dinoseb, Dinotefuran, Dinoterb, Dioxathion, Dipropetryn, Disulfoton-sulfone, Disulfoton-sulfoxide, Diuron (DCMU), DMST, DNOC, Dodemorph, Dodine, Edifenphos, Emamectin-benzoate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethiofencarb-sulfone, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrinfos, Famoxadone, Famphur, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenothiocarb, Fenoxaprop-ethyl, Fenoxycarb, Fencicoxamid, Fenpropathrin, Fenpropidin, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fensulfothion-sulfone, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfide, Fipronil-sulfone, Flamprop-isopropyl, Flamprop-methyl, Flamprop-M-isopropyl, Flazasulfuron, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluacrypyrim, Fluazifop, Fluazifop-P-butyl, Fluazinam, Fluazuron, Flucythrinate, Fludioxonil, Flufenacet, Flufenoxuron, Flumetralin, Fluometuron, Fluopicolide, Fluopyram, Fluoxastrobin, Fluquinconazole, Fluridone, Flurochloridone, Fluroxypyr, Flusilazole, Fluthiacet-methyl, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Furmecyclox, Haloxyfop, Haloxyfop-2-ethoxyethyl, Haloxyfop-methyl, Haloxyfop-P-methyl, Heptenophos, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazamox, Imazapyr, Imazaquin, Imazethapyr, Imibenconazole, Imidacloprid, Indoxacarb, Iodosulfuron-methyl, Ioxynil, Ipconazole, Iprobenfos, Iprovalicarb, Isazofos, Isocarbamid, Isocarbophos, Isofenphos-methyl, Isofetamid, Isomethiozin, Isopropcarb, Isopropalin, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaben, Isoxadifen-ethyl, Isoxaflutole, Isoxaflutole dikenotril, Isoxathion, Ivermectin B1a, Karanjin, Kresoxim-methyl, Linuron, Lufenuron, Malaaxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), MCPB, Mecarbam, Mefenacet, Mefenpyr-diethyl, Mefentrifluconazole, Mepanipyrim, Mephosfolan, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metazachlor ESA, Metazachlor OA, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methomyl, Methoprotryne, Methoxyfenozide, Metobromuron, Metolachlor, Metolcarb (MTMC), Metominostrobin (E, Z), Metoxuron, Metrafenone, Metribuzin, Metsulfuron-methyl, Mevinphos, Mexacarbate, Milbemectin A3, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N-(2,4-Dimethylphenyl)formamide, N-(2,4-dimethylphenyl)-N-methylformamidine, N,N-diethyl-m-toluamide, Napropamide, N-desmethyl-propamocarb, Neburon, Nereistoxin oxalate, Nicosulfuron, Nicotine, Nitenpyram, Norflurazon, Novaluron, Nuarimol, Ofurace, Omethoate, Oryzalin, Oxadiargyl, Oxadiazon, Oxadixyl, Oxamyl, Oxa-sulfuron, Oxathiapiprolin, Oxycarboxin, Oxydemeton-methyl, Oxymatrine, Paclobutrazol, Paraoxon-ethyl, Paraoxon-methyl, Pebulate, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthiopyrad, Pethoxamid, Phenothrin, Phenthoate, Phorate, Phorate-oxon, Phorate-oxon-sulfone, Phorate-sulfone, Phorate-sulfoxide, Phosalone, Phosfolan, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picaridin (Icaridin), Picolinafen, Picoxystrobin, Piperonyl-butoxide, Piperophos, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Potasan, Prallethrin, Pretilachlor, Prochloraz, Prochloraz met. BTS 44595, Prochloraz met. BTS 44596, Profenofos, Profoxydim, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propamocarb-N-oxide, Propanil, Propaquizafop, Propargite, Propazine, Propetamphos, Propham, Propiconazole, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prosulfuron, Prothioconazole-desthio,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrazoxon, Pyrethrins, Pyridaben, Pyridafol, Pyridalyl, Pyridaphenthion, Pyrifenox (E, Z), Pyrimethanil, Pyriofenone, Pyriproxyfen, Quinalphos, Quinclorac, Quinmerac, Quinoclamine, Quinoxifen, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rabenzazole, Resmethrin, Rimsulfuron, Rotenone, Saflufenacil, Sebuthylazine, Secbumeton, Sethoxydim, Siduron, Silafluofen, Silthiofam, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spiromesifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glucoside, Spirotetramat-keto-hydroxy, Spirotetramat-mono-hydroxy, Spiroxamine, Sulcotrione, Sulfosulfuron, Sulfotep, Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebupirimfos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrione, Temephos, Tepraloxymid, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiamethoxam, Thidiazuron, Thifensulfuron-methyl, Thiobencarb, Thiocyclam hydrogene oxalate, Thiodicarb, Thiofanox, Thiofanox-sulfone, Thiofanox-sulfoxide, Thiometon, Thionazin, Thiophanate-methyl, Tolclofos-methyl, Tolfenpyrad, Tolyfluanid, Triadimefon, Triadimenol, Tri-allate, Triasulfuron, Triazamate, Triazophos, Tribenuron-methyl, Trichlorfon, Triclopyr, Tricyclazole, Trietazine, Trifloxystrobin, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triflusaluron-methyl, Triforine, Triticonazole, Tritosulfuron, Valifenalate, Vamidothion, Vernolate, XMC, Zoxamide</i> Ukupno / In total 539 Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		DB; 2,4,5-T; Bentazone				
	A9-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,005 mg/kg 1-naftalenacetamid, 2,4-D, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Alaklor, Aldikarb, Aldikarb-sulfon (Aldoksikarb), Aldikarb-sulfoksid, Alidoklor, Ametoktradin, Ametrin, Aminokarb, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a, Azinfos-etil, Azinfos-metil, Azoksistrobin, Benalaksil, Bendiokarb, Bentazon, Benzovindflupir, Bifenazat, Bifentrin, Bitertanol, Biksafen, Boskalid, Bromacil, Bromokonazol, Bupirimat, Butafenacil, Butoksikarboksime, Buturon, Kadusafos, Karbaril (NAC), Karbendazim, Karbetamid, Karbosulfan, Karboksini, Karfentrazon-etil, Kloroantraniliprol, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Klorotoluron, Kloroksuron, Klorpirifos-metil, Kletodim, Klotentazin, Klotianidin, Kumafos, Kumatetralil, Ciantraniliprol, Ciazofamid, Cikloksidim, Cikluron, Ciprokonazol, Ciprodinil, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfon, Demeton-S-metil-sulfoksid, Desmedifam, Diafenturon, Diazinon, Diklofenak, Diklorprop, Diklobutrazol, Dikrotofos, Dietofenkarb, Difenokonazole, Diflubenzuron, Dimetenamid, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinotefuran, Dioksakarb, Diuron (DCMU), DMST, Dodin, Edifenfos, Emamektin-benzoat, Epoksikonazol, EPTC, Etakonazol, Etofenkarb, Etion, Etrimerol, Etofumezat, Etoprofos, Etofenproks, Etoksazol, Etrimerfos, Famoksadon, Fenamidon, Fenamifos, Fenamifos-sulfon, Fenamifos-sulfoksid, Fenarimol, Fenbukonazol, Fenheksamid, Fenobukarb, Fenoprop, Fenoksikarb, Fenpikoksamid, Fenpropidin, Fenpirazamin, Fenpiroksimat, Fensulfotion, Fention, Fention-okson, Fention-okson-sulfon, Fention-okson-sulfoksid, Fention-sulfon, Fention-sulfoksid, Fipronil, Fipronil-sulfon, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazinam, Fludioksonil, Flufenacet, Flufenoksuron, Fluometuron,	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	21/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Fluopikolid, Fluopiram, Fluoksastrobina, Flukvinkonazol, Fluridon, Flusilazol, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Fluksapiroksad, Fonofos, Forklorfenuron, Formetanat hidroklorid, Fostiazat, Fuberidazol, Furalaksil, Heksakonazol, Heksaflumuron, Heksazinon, Heksitiazoks, Hidrametilnon, Imazalil, Imazetapir, Imidaklopid, Indoksakarb, Ioksonil, Ipkonazol, Iprovalikarb, Isazofos, Izokarbofos, Izofenfos-metil, Izofetamid, Izoprokarb, Izopropalin, Izoproturon, Izopirazam, Izoksaflutol, Izoksaflutol-dikenotril, Ivermektin B1a, Kresoksime-metil, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), Mefenacet, Mefentriklonazol, Mepronil, Meptildinokap, Metaflumizon, Metalaksil i Metalaksil-M, Metamitron, Metazaklor, Metkonazol, Metabenzotiazuron, Metakrifos, Metamidofos, Metidation, Metiokarb, Metiokarb-sulfon, Metiokarb-sulfoksid, Metomil, Metoprotin, Metoksifenozid, Metobromuron, Metolaklor, Metoksuron, Metrafenon, Metribuzin, Mevinfos, Meksakarb, Molinat, Monokrotofos, Monolinuron, Monuron, Miklobutanil, N,N-dietil-m- toluamid, Napropamid, Neburon, Novaluron, Ometoat, Oksadiargil, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Paklobutrazol, Penkonazol, Pencikuron, Pendimetalin, Penflufen, Pentiopirad, Fenmedifam, Fenotrin, Fentoat, Forat, Fosalon, Fosmet, Fosmet-okson, Fosfamidon, Foksim, Pikoksistobin, Pirimikarb, Pirimikarb-desmetil, Pirimifos-metil, Prokloraz, Profenofos, Promekarb, Prometon, Prometrin, Propamokarb, Propargit, Propazin, Profam, Propikonazol, Propoksur, Prosulfokarb, Protiokonazol-destio, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Piridaben, Piridalil, Piridafention, Piridat, Kvinalfos, Kvinoklamin, Kvizalofop (slobodna kiselina), Kvizalofop-etil, Kvizalofop-metil, Kvizalofop- P, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Rotenon, Sebutilazin, Sekbumeton, Siduron, Simazin, Simetrin, Spinetoram L, Spinosad (suma Spinosin A i Spinosin D), Spiromezifen, Spirotetramat, Spirotetramat-enol-glukozid, Spirotetramat-mono-hidroksi, Spiroksamin, Sulfotep, Sulfoksaflor, Sulprofos, TCMTB, Tebukonazol, Tebufenozid, Tebufenpirad, Tebutiuron, Teflubenzuron, Temefos, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetraklorvinfos (CVMP), Tetraconazole, Tetrametrin, Tiabendazol, Tiaklopid, Tiobenkarb, Tiodikarb, Tolfenpirad, Tolilfluaniid, Triadimefon, Triadimenol, Triazofos, Triklorfon, Triklorpir, Triciklazol, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Tritikonazol, Tritosulfuron, Vamidotion, Vernolat, Vinklozolin, XMC, Zoksamid /				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>1-naphthaleneacetamide, 2,4-D, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Alachlor, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a, Azinphos-ethyl, Azinphos-methyl, Azoxystrobin, Benalaxyl, Bendiocarb, Bentazone, Benzovindflupyr, Bifenazate, Bifenthrin, Bitertanol, Bixafen, Boscalid, Bromacil, Bromuconazole, Bupirimate, Butafenacil, Butoxycarboxim, Buturon, Cadusafos, Carbaryl (NAC), Carbendazim, Carbetamide, Carbosulfan, Carboxin, Carfentrazone-ethyl, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chlorotoluron, Chloroxuron, Chlorpyrifos-methyl, Clethodim, Clofentezine, Clothianidin, Coumaphos, Coumatetralyl, Cyantraniliprole, Cyazofamid, Cycloxydim, Cycluron, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Diazinon, Dichlofenac, Dichlorprop, Diclobutrazol (stereo isomer), Dicrotophos, Diethofencarb, Difenoconazole (isomer), Disflubenzuron, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimoxystrobin, Diniconazole, Dinotefuran, Dioxacarb, Diuron (DCMU), DMST, Dodine, Edifenphos, Enamectin-benzoate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrinfos, Famoxadone, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenoxycarb, Fenpicoxamid, Fenpropidin, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fipronil, Fipronil-sulfone, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluazifop, Fluazinam, Fludioxonil, Flufenacet, Flufenoxuron, Fluometuron, Fluopicolide, Fluopyram, Fluoxastrobin, Fluquinconazole, Fluridone, Flusilazole, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazethapyr, Imidacloprid, Indoxacarb, Ioxynil, Ipconazole, Iprovalicarb, Isazofos, Isocarboxiphos, Isafenphos-methyl, Isofetamid, Isoprocacarb, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaflutole, Isoxaflutole dikenotriple, Ivermectin B1a, Kresoxim-methyl, Lenacil, Linuron, Lufenuron, Malaoxon, Malathion,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Mandipropamid, Matrine, MCPA (MCP), Mefenacet, Mefentrifluconazole, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotetryne, Methoxyfenozide, Metobromuron, Metolachlor, Metoxuron, Metrafenone, Metribuzin, Mevinphos, Mexacarbate, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N,N-diethyl-m-toluamide, Napropamide, Neburon, Novaluron, Omethoate, Oxadiargyl, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapiprolin, Oxydemeton-methyl, Paclobutrazol, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthioopyrad, Phenmedipham, Phenothrin, Phenthoate, Phorate, Phosalone, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picoxystrobin, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Prochloraz, Profenofos, Promecarb, Prometon, Prometryn, Propamocarb, Propargite, Propazine, Propham, Propiconazole, Propoxur, Prosulfocarb, Prothioconazole-desthio, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyridaben, Pyridalyl, Pyridaphenthion, Pyridate, Quinalphos, Quinoclamine, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rotenone, Sebutylazine, Secbumeton, Siduron, Simazine, Simetryn, Spinetoram L, Spinosad (sum of spinosyn A and spinosyn D), Spiromesifen, Spirotetramat, Spirotetramat-enol-glucoside, Spirotetramat-mono-hydroxy, Spiroxamine, Sulfotep, Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebuthiuron, Teflubenzuron, Temephos, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiobencarb, Thiodicarb, Tolfenpyrad, Tolyfluanid, Triadimefon, Triadimenol, Triazophos, Trichlorfon, Triclopyr, Tricyclazole, Triflumizole, Triflumizole metabolite FM-6-1, Triflumurion, Triticonazole, Tritosulfuron, Vamidothion, Vernolate, Vinclozolin, XMC, Zoxamide Ukupno / In total 327 Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS)				

Oznaka/ Identification	Materijali/ Proizvodi <i>Materials/ Products</i>	Vrsta ispitivanja/Svojstvo <i>Type of test / Property</i> Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja <i>Test method</i>	U primjeni od/ <i>In use from</i>	Datum posljednje izmjene/ <i>Date of last change</i>
		<i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-DB; 2,4,5-T; Bentazon				
	A9-4 Žitarice i proizvodi od žitarica – visoki udio škroba i/ili proteina te niski udio vode i masti <i>Cereals and cereals products – high starch and/or protein content and low water and fat content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,01 mg/kg 1-naftalenacetamid, 2,4-D, 2,4-dimetilanilin, 3-hidroksikarbofuran, 4-klorofenoksiocena kiselina , Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Akrinatrin, Alaklor, Aldikarb, Aldikarb-sulfon (Aldoksikarb), Aldikarb-sulfoksid, Alidoklor, Ametoktradin, Ametrin, Aminokarb, Atraton, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a, Azinfos-etil, Azinfos-metil, Azoksistrobin, Benalaksil, Bendiokarb, Benfurakarb, Benomil, Bentazon, Benzovindflupir, Benzoksimat, Bifenazat, Bifentrin, Bitertanol, Biksafen, Boskalid, Bromacil, Bromoksinil, Bromukonazol, Bupirimat, Buprofezin, Butaklor, Butafenacil, Butokarboksim, Butoksikarboksim, Buturon, Butilat, Karbaril (NAC), Karbendazim, Karbetamid, Karbofuran, Karbofenotion, Karbosulfan, Karboksim, Karfentrazon-etil, Kloroantraniliprol, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Klodiazon, Klorotoluron, Kloroksuron, Klorprofam, Klorpirifos, Klorpirifos-metil, Kletodim, Klofentezin, Klomazon, Klotianidin, Kumafos, Kumatetrilil, Cianazin, Ciantraniliprol, Ciazofamid, Cikloat, Cikloksidim, Cikluron, Ciflufenamid, Cimoksamil, Ciprokonzol, Ciprodinil, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfon, Demeton-S-metil-sulfoksid, Desmedifam, Diafentiuron, Diallate, Diazinon, Diklofenak, Diklorprop, Diklorvos, Diklobutrazol, Dikrotfos, Dietofenkarb, Difenokonazole, Difetialon, Diflubenzuron, Dimetaklor, Dimetenamid,	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	21/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinotefuran, Dioksakarb, Difenamid, Diuron (DCMU), DMST, Dodin, Edifenfos, Enamektin-benzoat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etion, Etirimol, Etofumezat, Etoprofos, Etofenproks, Etoksakonazol, Etrimfos, Famoksadon, Fenamidon, Fenamifos, Fenamifos-sulfon, Fenamifos-sulfoksid, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenobukarb, Fenoprop, Fenoksikarb, Fenpikoksamid, Fenpropatrin, Fenpropidin, Fenpirazamin, Fenpiroksimat, Fensulfotion, Fention, Fention-okson, Fention-okson-sulfon, Fention-okson-sulfoksid, Fention-sulfon, Fention-sulfoksid, Fenuron, Fipronil, Fipronil-sulfon, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazifop-P-butil, Fluazinam, Flucitrinat, Fludioksonil, Flufenacet, Flufenoksuron, Fluometuron, Fluopikolid, Fluopiram, Fluoksastrobin, Flukvinkonazol, Fluridon, Flusilazol, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Fluksapiroksad, Fonofos, Forklorfenuron, Formetanat hidroklorid, Fostiazat, Fuberidazol, Furalaksil, Furatiokarb, Haloksifop, Heksakonazol, Heksaflumuron, Heksazinon, Heksitiazoks, Hidrametilnon, Imazalil, Imidaklopid, Indoksakarb, Ioksonil, Ipkonazol, Iprovalikarb, Izazofos, Izofenfos-metil, Izofetamid, Izoprokarb, Izopropalin, Izoprotiolan, Izoproturon, Izopirazam, Izoksaflutol, Izoksaflutol-dikenotril, Ivermektin Bla, Kresoksim-metil, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), Mefenacet, Mefentriflukonazol, Mepanipirim, Mepronil, Meptildinokap, Metaflumizon, Metalaksil i Metalaksil-M, Metamitron, Metazaklor, Metkonazol, Metabenzthiazuron, Metakrifos, Metamidofos, Metidation, Metiokarb, Metiokarb-sulfon, Metiokarb-sulfoksid, Metomil, Metoprotrin, Metoksifenozyd, Metobromuron, Metolaklor, Metoksuron, Metrafenon, Metribuzin, Mevinfos, Meksakarbat, Molinat, Monokrotofos, Monolinuron, Monuron, Miklobutanil, N,N-dietil-m-toluamid, Naled, Napropamid, Neburon, Nereistoksin oksalat, Nikotin, Nitenpiram, Norflurazon, Novaluron, Nuarimol, Ometoat, Oksadiargil, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapiprolin, Oksidemeton-metil, Paklobutrazol, Pebulat, Penkonazol, Pencikuron, Pendimetalin, Penflufen, Pentiopirad, Fenmedifam, Fenotrin, Fentoat, Forat, Fosalon, Fosmet, Fosmet-okson, Fosfamidon, Foksim, Pikolinafen, Pikoksistrobin, Piperonil-butoksid, Pirimikarb, Pirimikarb-desmetil, Pirimifos-metil, Pretilaklor, Prokloraz, Profenofos, Promekarb, Prometon, Prometrin, Propaklor, Propamokarb, Propanil, Propakvizafop, Propargit, Propazin, Profam, Propikonazol, Propoksur, Propizamid,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Prokvinazid, Prosulfokarb, Protiokonazol-destio, Pimetrozin, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Piretrini, Piridaben, Piridalil, Piridafention, Piridat, Pirimetanil, Piriofenon, Piriproksifen, Kvinalfos, Kvinoklamin, Kvinoksifen, Kvizalofop (slobodna kiselina), Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Resmetrin, Rotenon, Sebutilazin, Sekbumeton, Siduron, Silafluofen, Simazin, Simetrin, Spinetoram, Spinosad (suma Spinosin A i Spinsin D), Spirodiklofen, Spiromezifen, Spirotetramat, Spirotetramat-enol-glukozid, Spirotetramat-mono-hidroksi, Spiroksamin, Sulfotep, Sulfoksaflor, Sulprofos, TCMTB, Tebukonazol, Tebufenozid, Tebufenpirad, Tebutiuron, Teflubenzuron, Temefos, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetraklorvinfos (CVMP), Tetraconazole, Tetrametrin, Tiabendazol, Tiaklopid, Tiametoksam, Tiazuron, Tiobenkarb, Tiociklam hidrogen oksalat, Tiodikarb, Tiometon, Tiofanat-metil, Tiram, Tolklofos-metil, Tolfenpirad, Tralometrin, Triadimefon, Triadimenol, Trialat, Triazofos, Triklorfon, Triklorpir Triciklazol, Trifloksistrobin, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Tritikonazol, Tritosulfuron, Vamidotion, Vernolat, Vinklozolin, XMC, Zoksamid / <i>1-naphthaleneacetamide, 2,4-D, 2,4-dimethylaniline, 3-Hydroxycarbofuran, 4-Chlorophenoxyacetic acid, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Acrinathrin, Alachlor, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametocradin, Ametryn, Aminocarb, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a, Azinphos-ethyl, Azinphos-methyl, Azoxystrobin, Benalaxyl, Bendiocarb, Benfuracarb, Benomyl, Bentazone, Benzovindflupyr, Benzoximate, Bifenazate, Bifenthrin, Bitertanol, Bixafen, Boscalid, Bromacil, Bromoxynil, Bromuconazole, Bupirimate, Buprofezin, Butachlor, Butafenacil, Butocarboxim, Butoxycarboxim, Buturon, Butylate, Carbaryl (NAC), Carbendazim, Carbetamide, Carbofuran, Carbophenothion, Carbosulfan, Carboxin, Carfentrazon-ethyl, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chlorotoluron, Chloroxuron, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Clethodim, Clofentezine, Clomazone, Clothianidin, Coumaphos, Coumatetralyl, Cyanazine, Cyantraniliprole, Cyazofamid, Cycloate, Cycloxydim, Cycluron, Cyflufenamid, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Di-allate,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Diazinon, Dichlofenac, Dichlorprop, Dichlorvos, Diclobutrazol, Dicrotophos, Diethofencarb, Difenconazole, Difethialone, Diflubenzuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimoxystrobin, Diniconazole, Dinotefuran, Dioxacarb, Diphenamid, Diuron (DCMU), DMST, Dodine, Edifenphos, Emamectin-benzoate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrinfos, Famoxadone, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenoxycarb, Fenpicoxamid, Fenpropathrin, Fenpropidin, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfone, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Floryprauxifen-benzyl, Fluazifop, Fluazifop-P-butyl, Fluazinam, Flucythrinate, Fludioxonil, Flufenacet, Flufenoxuron, Fluometuron, Fluopicolide, Fluopyram, Fluoxastrobin, Fluquinconazole, Fluridone, Flusilazole, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Haloxyfop, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imidacloprid, Indoxacarb, Ioxynil, Ipconazole, Iprovalicarb, Isazofos, Isofenphos-methyl, Isofetamid, Isoprocab, Isopropalin, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaflutole, Isoxaflutole dikenotriole, Ivermectin B1a, Kresoxim-methyl, Lenacil, Linuron, Lufenuron, Malaoxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), Mefenacet, Mefentrifluconazole, Mepanipyrim, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotetryne, Methoxyfenozide, Metobromuron, Metolachlor, Metoxuron, Metrafenone, Metribuzin, Mevinphos, Mexacarbate, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N,N-diethyl-m-toluamide, Naled, Napropamide, Neburon, Nereistoxin oxalate, Nicotine, Nitenpyram, Norflurazon, Novaluron, Nuarimol, Omethoate, Oxadiargyl, Oxadiazon, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapiprolin, Oxydemeton-methyl, Paclobutrazol, Pebulate, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthiopyrad, Phenmedipham, Phenothrin, Phenthoate, Phorate, Phosalone, Phosmet, Phosmet-oxon,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Phosphamidon, Phoxim, Picolinafen, Picoxystrobin, Piperonyl-butoxide, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Pretilachlor, Prochloraz, Profenofos, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propanil, Propaquizafop, Propargite, Propazine, Propham, Propiconazole, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrethrins, Pyridaben, Pyridalyl, Pyridaphenthion, Pyridate, Pyrimethanil, Pyriofenone, Pyriproxyfen, Quinalphos, Quinoclamine, Quinoxifen, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, - Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Resmethrin, Rotenone, Sebuthylazine, Secbumeton, Siduron, Silafluofen, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spiromesifen, Spirotetramat, Spirotetramat-enol-glucoside, Spirotetramat-mono-hydroxy, Spiroxamine, Sulfotep, Sulfoxaflo, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebuthiuron, Teflubenzuron, Temephos, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiachloprid, Thiamethoxam, Thidiazuron, Thiobencarb, Thiocyclam hydrogene oxalate, Thiodicarb, Thiometon, Thiophanate-methyl, Thiram, Tolclofos-methyl, Tolfenpyrad, Tralomethrin, Triadimefon, Triadimenol, Tri-allate, Triazophos, Trichlorfon, Triclopyr, Tricyclazole, Trifloxystrobin, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triticonazole, Tritosulfuron, Vamidothion, Vernolate, Vinclozolin, XMC, Zoxamide</i> Ukupno/ In Total: 392 Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		DB; 2,4,5-T; Bentazone				
	A9-5 "Teški matriksi" "Difficult or unique commodities"	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,01 mg/kg 1-naftalacetamid, 2,3,5-trimetakarb, 2,4,5-T, 3,4,5-trimetakarb, 3-Decen-2-on, Acefat, Acibenzolar-S-metil, Alaklor, Alanikarb, Aldikarb, Aldikarb-sulfoksid, Aldikarb-sulfon (Aldoksikarb), Alidoklor, Ametoktradin, Ametrin, Amidosulfuron, Amisulbrom, Atraton, Atrazin, Atrazin-desetil, Atrazin-desizopropil, Avermektin B1a, Azakonazol, Azametifos, Azinfos-metil, Aziprotrin, Azoksistrobin, Benalaksil, Bendiokarb, Benodanil, Bensulfuron-metil, Bensulid, Bentazon, Benthiavalikarb, Benthiavalikarb-izopropil, Benziladenin, Benzoilprop-etil, Benzoksimat, Benzovindflupir, Bifenazat, Bifenazat-diazen, Bifentrin, Biksafen, Bitertanol, Boscalid, Bromacil, Bromukonazol, Bupirimat, Butafenacil, Butilat, Butoksikarboksime, Butralin, Buturon, Centrimonijev klorid, Ciflufenamid, Ciflumetofen, Cijanazin, Cijanofenfos, Cijantraniliprol, Ciklanilid, Cikloksidim, Cikluron, Cimiazol, Cimoksanil, Cinosulfuron, Ciprodinil, Ciprokonazol, Ciromazin, Cyazofamid, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfoksid, Demeton-S-metil-sulfon, Diazinon, Dietofenkarb, Difenamid, Difenokonazol, Difenoksuron, Diflubenzuron, Diflufenikan, Dijafentiuron, Diklobutrazol, Diklofenak, Dikrotofos, Dimefuron, Dimetaklor, Dimetenamid, Dimetilan, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinoseb, Dinotefuran, Dinoterb, Dioksakarb, Dipropetrin, Disulfoton-sulfoksid, Disulfoton-sulfon, Diuron (DCMU), DMST, Dodemorf, Dodin, Edifenfos, Eamektin-benzoat, Epoksikonazol, Etakonazol, Etiofenkarb, Etion, Etirimol, Etofenproks, Etoksazol, Etoprofos, Etrimfos, Famoksadon, Famphur, Fenamidon, Fenamifos, Fenamifos-sulfoksid, Fenamifos-sulfon, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenmedifam, Fenobukarb, Fenoksaprop-etil,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Fenoksikarb, Fenotiokarb, Fenotrin, Fenpikoksamid, Fenpirazamin, Fenpiroksimat, Fenpropimorf, Fensulfotion, Fensulfotion-sulfon, Fention-okson-sulfoksid, Fention-okson-sulfon, Fention-sulfoksid, Fenuron, Fipronil, Fipronil-sulfid, Fipronil-sulfon, Flamprop-izopropil, Flamprop-metil, Flamprop-M-izopropil, Flazasulfuron, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluakripirim, Fluazifop-P-butil, Fluazinam, Fluazuron, Fludioksonil, Flufenacet, Flufenoksuron, Fluksapiroksad, Flukvinkonazol, Fluoksastobin, Fluopiram, Fluridon, Flusilazol, Flutiacet-metil, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Foksim, Forat, Forat-okson, Forat-okson-sulfon, Forat-sulfoksid, Forat-sulfon, Forklorfenuron, Formetanat hidroklorid, Fosalon, Fosfamidon, Fosfolan, Fosmet, Fosmet-okson, Fostiazat, Fuberidazol, Furalaksil, Furathiokarb, Furmecikloks, Haloksifop-2-etoksietil, Heksafumuron, Heksakonazol, Heksazinon, Heksitiazoks, Heptenofos, Hidrametilnon, Imazakvin, Imazalil, Imibenkonazol, Imidaklopid, Indoksakarb, Ipkonazol, Iprobenfos, Iprovalikarb, Isazofos, Isofetamid, Isoksaben, Isoprocab, Ivermektin B1a, Izokarbamid, Izoksadifen-etil, Izoksafutol, Izoksafutol dikenotril, Izoksation, Izometiozin, Izopirazam, Izoprotiolan, Izoproturon, Jodosulfuron-metil, Kadusafos, Karanjin, Karbaril (NAC), Karbendazim, Karbetamid, Karbofuran, Karboksinsulfoksid, Karbosulfan, Klethodim, Klimbazol, Klodinafop (slobodna kiselina), Klodinafop-propargil, Klofentezin, Klokintocet-meksil, Klomazon, Klorantraniliprol, Klorbromuron, Klorfenvinfos (E, Z), Klorfluazuron, Kloridazon, Kloroksuron, Klorosulfuron, Klotianidin, Krotoksifos, Krufomat, Kumafos, Kumatetralil, Kvinalfos, Kvinoksifen, Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, Mefenacet, Mefenpir-dietil, Mefentriklukonazol, Mefosfolan, Mekarbam, Meksakarb, Mepanipirim, Meptildinokap, Metabenzthiazuron, Metaflumizon, Metakrifos, Metalaksil i Metalaksil-M, Metamidofos, Metazaklor, Metazaklor OA, Metidation, Metiokarb, Metiokarb-sulfoksid, Metiokarb-sulfon, Metkonazol, Metobromuron, Metoksifenoizid, Metoksuron, Metolaklor, Metolkarb (MTMC), Metomil, Metominostrobilin (E, Z), Metoprotrin, Metrafenon, Metribuzin, Metsulfuron-metil, Miklobutanil, Molinat, Monokrotofos, Monolinuron, Monuron, N-(2,4-dimetilfenil)formamid, N,N-dietil-m-toluamid, Napropamid, N-desmetil-propamokarb, Neburon, Nikosulfuron, Nitenpiram, Norflurazon, Novaluron, Nuarimol, Ofurak, Oksadiazon,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Oksikarboksin, Ometoat, Paklobutrazol, Paraokson-etil, Pencikuron, Pendimetalin, Penflufen, Penkonazol, Pentiopirad, Petoksamid, Pikaridin (Ikaridin), Pikoksistrobin, Pikolinafen, Pimetrozin, Piperofos, Piperonil-butoksid, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Pirazokson, Piridaben, Piridafention, Piridalil, Piridat, Pirimetanil, Pirimicarb, Pirimifos-metil, Pirimikarb-desmetil, Piriofenon, Piriproksifen, Potasan, Praletrin, Pretilaklor, Profenofos, Profoksidim, Prokloraz, Prokloraz met. BTS 44595, Prokloraz met. BTS 44596, Prokvinazid, Promecarb, Prometon, Prometrin, Propaklor, Propakvizafop, Propamocarb, Propamokarb-N-oksidi, Propargit, Propazin, Propetamfos, Propizamid, Propoksur, Prosulfokarb, Prosulfuron, Protiokonazol-destio, Rabenzazol, Rimsulfuron, Rotenon, Saflufenacil, Sebutilazin, Sekbumeton, Setoksidim, Siduron, Silthiofam, Simazin, Simetrin, Spinetoram, Spinosad (zbroj spinosina A i spinosina D), Spirodiklofen, Spiroksamin, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glukozid, Spirotetramat-keto-hidroksi, Spirotetramat-mono-hidroksi, Sulfoksaflo, Sulfosulfuron, Sulfotep, Sulprofos, TCMTB, Tebufenozid, Tebufenpirad, Tebukonazol, Tebupirimfos, Tebutam, Tebutiuron, Teflubenzuron, Temephos, Tepraloksidim, Terbumeton, Terbutilazin, Terbutrin, Tetraklorinfos (CVMP), Tetrakonazol, Tetrametrin, Thionazin, Tiaklopid, Tiametoksam, Tiazuron, Tifensulfuron-metil, Tiobenkarb, Tiodicarb, Tiofanat-metil, Tiofanoks-sulfoksid, Tolfenpirad, Triadimefon, Triadimenol, Triasulfuron, Triazamat, Triazofos, Tribenuron-metil, Triciklazol, Trietazin, Triflumizol, Triflumizol metabolit FM-6-1, Triflururon, Triflururon-metil, Triflorin, Tritikonazol, Tritosulfuron, Valifenalat, Vamidotin, XMC, Zoksamid / <i>1-naphthaleneacetamide, 2,3,5-Trimethacarb, 2,4,5-T, 3,4,5-Trimethacarb, 3-Decen-2-one, Acephate, Acibenzolar-S-methyl, Alachlor, Alanycarb, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Amidosulfuron, Amisulbrom, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a, Azaconazole, Azamethiphos, Azinphos-methyl, Aziprotryne, Azoxystrobin, Benalaxyl, Bendiocarb, Benodanil, Bensulfuron-methyl, Bensulide, Bentazone, Benthialdicarb, Benthialdicarb-isopropyl, Benzovindflupyr, Benzoximate, Benzoylprop-ethyl, Benzyladenine, Bifenazate, Bifenazate-diazene, Bifenthrin, Bitertanol, Bixafen, Boscalid, Bromacil, Bromuconazole,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Bupirimate, Butafenacil, Butoxycarboxim, Butralin, Buturon, Butylate, Cadusafos, Carbaryl (NAC), Carbensulfam, Carbetamide, Carbofuran, Carbosulfan, Carboxin sulfoxide, Centrimonium chloride, Chlorantraniliprole, Chlorbromuron, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chloroxuron, Chlorsulfuron, Cinosulfuron, Clethodim, Climbazole, Clodinafop (free acid), Clodinafop-propargyl, Clofentezine, Clomazone, Cloquintocet-mexyl, Clothianidin, Coumaphos, Coumatetralyl, Crotoxyphos, Crufomate, Cyanazine, Cyanofenphos, Cyantraniliprole, Cyazofamid, Cyclanilide, Cycloxydim, Cycluron, Cyflufenamid, Cyflumetofen, Cymiazole, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Diafenthiuron, Diazinon, Dichlofenac, Diclobutrazol, Dicrotophos, Diethofencarb, Difenoconazole, Difenoxyuron, Diflubenzuron, Diflufenican, Dimefuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimetilan, Dimoxystrobin, Diniconazole, Dinoseb, Dinotefuran, Dinoterb, Dioxacarb, Diphenamid, Dipropetryn, Disulfoton-sulfone, Disulfoton-sulfoxide, Diuron (DCMU), DMST, Dodemorph, Dodine, Edifenphos, Emetectin-benzoate, Epoxiconazole, Etaconazole, Ethiofencarb, Ethion, Ethirimol, Ethoprophos, Etofenprox, Etoazole, Etrifos, Famoxadone, Famphur, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenothiocarb, Fenoxaprop-ethyl, Fenoxycarb, Fenpicoxamid, Fenpropimorph, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fensulfothion-sulfone, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfide, Fipronil-sulfone, Flamprop-isopropyl, Flamprop-methyl, Flamprop-M-isopropyl, Flazasulfuron, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluacrypyrim, Fluazifop-P-butyl, Fluazinam, Fluazuron, Fludioxonil, Flufenacet, Flufenoxuron, Fluopyram, Fluoxastrobilin, Fluquinconazole, Fluridone, Flusilazole, Fluthiacet-methyl, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Furmecyclox, Haloxyfop-2-ethoxyethyl, Heptenophos, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazaquin, Imibenconazole, Imidacloprid, Indoxacarb, Iodosulfuron-methyl, Ipconazole, Iprobenfos, Iprovalicarb,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Isazofos, Isocarbamid, Isofetamid, Isomethiozin, Isoprocab, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaben, Isoxadifen-ethyl, Isoxaflutole, Isoxaflutole dikenotriple, Isoxathion, Ivermectin B1a, Karanjin, Lenacil, Linuron, Lufenuron, Malaaxon, Malathion, Mandipropamid, Matrine, Mecarbam, Mefenacet, Mefenpyr-diethyl, Mefentrifluconazole, Mepanipyrim, Mephosfolan, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metazachlor, Metazachlor OA, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotrene, Methoxyfenozide, Metobromuron, Metolachlor, Metolcarb (MTMC), Metominostrobin (E, Z), Metoxuron, Metrafenone, Metribuzin, Metsulfuron-methyl, Mexacarbate, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N-(2,4-Dimethylphenyl)formamide, N,N-diethyl-m-toluamide, Napropamide, N-desmethyl-propamocarb, Neburon, Nicosulfuron, Nitenpyram, Norflurazon, Novaluron, Nuairimol, Ofurace, Omethoate, Oxadiazon, Oxadixyl, Oxamyl, Oxsulfuron, Oxathiapiprolin, Oxycarboxin, Oxydemeton-methyl, Paclobutrazol, Paraoxon-ethyl, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthioopyrad, Pethoxamid, Phenmedipham, Phenothrin, Phorate, Phorate-oxon, Phorate-oxon-sulfone, Phorate-sulfone, Phorate-sulfoxide, Phosalone, Phosfolan, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picaridin (Icaridin), Picolinafen, Picoxystrobin, Piperonyl-butoxide, Piperophos, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Potasan, Prallethrin, Pretilachlor, Prochloraz, Prochloraz met. BTS 44595, Prochloraz met. BTS 44596, Profenofos, Profoxydim, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propamocarb-N--oxide, Propaquizafop, Propargite, Propazine, Propetamphos, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prosulfuron, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrazoxon, Pyridaben, Pyridalyl, Pyridaphenthion, Pyridate, Pyrimethanil, Pyriofenone, Pyriproxifen, Quinalphos, Quinoxifen, Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rabenzazole, Rimsulfuron, Rotenone, Saflufenacil, Sebuthylazine, Secbumeton, Sethoxydim, Siduron, Silthiofam, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glucoside, Spirotetramat-keto-hydroxy, Spirotetramat-mono-hydroxy, Spiroxamine, Sulfosulfuron, Sulfotep,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebupirimfos, Tebutam, Tebuthiuron, Teflubenzuron, Temephos, Tepraloxymid, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiacloprid, Thiamethoxam, Thidiazuron, Thifensulfuron-methyl, Thiobencarb, Thiodicarb, Thiofanox-sulfoxide, Thionazin, Thiophanate-methyl, Tolfenpyrad, Triadimefon, Triadimenol, Triasulfuron, Triazamate, Triazophos, Tribenuron-methyl, Tricyclazole, Trietazine, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triflurosulfuron-methyl, Triforine, Triticonazole, Tritosulfuron, Valifenalate, Vamidothion, XMC, Zoxamide Ukupno/ In Total: 435				
A10	A10-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 1-Naftiloctena kiselina, 2,3,5,6-Tetrakloroanilin, 2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Akrinatri (suma izomera), Alaklor, Aldrin, Aletrin -3,4 (Bioaletrin), Alidoklor, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Benfluralin, Benziladenin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Binapacryl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Kaptan, Karbofenotion, Karbofenotion metil, Kinometionat, Klorbufam, Klordekon, Klorotoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorobenzilat, Kloroneb, Klorothalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klorzinat, Cinidon etil, cis-Kloran, cis-Permetrin, Klodinafop propargil, Klomazon, Klorvintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Ci-	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		permetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Demeton-S-metil, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklofluamid, Diklorvos, Diklobutrazol, Dikloran, Dikofol, Dieldrin, Dimefoks, Dimetaklor, Difenamid, Dimetipin, Dimetilan, Dinitramin, Difenilamid, Disulfoton, Ditalimfos, Edifenfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Etrimfos, Famfur, Fenarimol, Fenklorfos, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flupiradifuron, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutriafol, Fluvalinat-tau (suma izomera), Folpet, Formotion, Furilazol, gama-BHC (lindan), Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izodrin, Izoksaben, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lencil, Leptofos, Linuron, Mandestrobilin, Metalaksil (Mefenoksam), Metazaklor, Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitalin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Pebulat, Penkonazol, Pendimetalin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Protio-konazol, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirime-tanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebu-konazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tio-nazin, THPI (Tetrahidroftalimid), Tiokarbamil, Tolklofos-metil, Tolilfluamid, Tolilfluamid metabolit, trans-Klordan, Translutrin, trans-Permetrin, Triadime-fon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Triciklazol, Tride-morf, Trietazin, Trifluralin, Unikonazol, Valifenalat, Vernolat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 1-Naphthylacetic acid,, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Dichlorobenzophenone, 4,4'-Methoxychlor olefin,, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allethrin-3,4 (Bioallethrin), Allidochlor, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Benfluralin, Benzyladenine, beta- BHC, beta-Endosulfan, Bifenthrin, Binapacryl, Biphenyl, Binapacryl, Bitrex (Denatonium benzoat), Bromacil, Bromfeninfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Captan, Carbophenothion, Carbophenothion methyl, Chinomethionat, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozoline, Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Clomazone, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Demeton-S-methyl, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlofluanid, Dichlorvos, Diclobutrazol, Dicloran, Dicofol, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Dinitramine, Diphenamid, Diphenylamine, Disulfoton, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfuralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Etrinfos, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flupyradifurone, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutriafol, Fluvalinate (sum of isomers), Folpet, Formothion, Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isodrin, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitalin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychloridan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Ofurace, Paclobutrazol, Parathion, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phorate, Phosfolan,, Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prothioconazole, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tricyclazole, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vernolate, Vinclozolin Ukupno / In total 264				
	A10-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroanilin, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Aktrin (suma izomera), Alaklor, Aldrin, Aletrin -3,4 (Bioaletrin), alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Azinfos-methyl, Benfluralin, Benziladenin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Bupirimat, Karbofenotion, Karbofenotion metil, Kinometionat, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorbenzid, Klorobenzilat, Kloroneb, Klorothalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Cinidon etil, cis-Klordan, cis-	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Permetrin, Klodinafop propargil, Klokintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklorvos, Diklobutrazol, Dikloran, Dikofol, Dieldrin, Dimefoks, Dimetaklor, Difenamid, Dimetipin, Dimetilan, Dinitramin, Difenilamid, Disulfoton, Ditalimfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Famfur, Fenarimol, Fenklorfos, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutriafol, Fluvalinat-tau (suma izomera), Fonofos, Formotion, Fostiazat (suma izomera), Furilazol, gama-BHC (lindan), Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrion, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izofenfos metil, Izoksaben, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, Nitrilin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Orizalin, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Pebulat, Penkonazol, Pendimetalin, Pentakloroanilin, Pentakloroanisol, Pentaklorobenzen, Pentanoklor (Solan), Fosfolan, Fosmet, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Protiokonazol, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Spiromezifen, Sulkotrin, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbamil, Tolklofos-metil, Tolilfluamid, Tolilfluamid metabolit, trans-Klordan, Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Tridemorf, Trietazin, Trifluralin, Unikonazol, Valifenalat, Vernolat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Tetrachloroaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allethrin-3,4 (Bioallethrin), alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Azinphos methyl, Benfluralin, Benzyladenine, beta-BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Binapacryl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Carbophenothion, Carbophenothion methyl, Chinomethionat, Chlorbufam, Chlordacone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlorvos, Diclobutrazol, Dicloran, Dicofol, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Dinitramine, Diphenamid, Diphenylamine, Disulfoton, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutriafol, Fluvalinate (sum of isomers), Folpet, Formothion, Fosthiazate (sum of isomers), Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos-methyl, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, Nitralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Ofurace, Oryzalin, Paclobutrazol, Parathion, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phosfolan, Phosmet, Piperonyl				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prothioconazole, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesifen, Sulcotrione, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vernolate, Vinclozolin Ukupno / In total 257				
	A10-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroaniline, 2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Akrinatrin (suma izomera), Alaklor, Aldrin, Aletrin - 3,4 (Bioaletrin), Alidoklor, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Benfluralin, Benziladenin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Kaptan, Karbofenotion, Karbofenotion metil, Karfentrazon etil, Kinometionat, Klorbenzid, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorobenzilat, Kloroneb, Klorothalonil, Klorprofam,	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Klordan, cis-Permetrin, Klodinafop propargil, Klokintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklofluanid, Diklorvos, Diklobutrazol, Dikloran, Dikofol, Dieldrin, Dimefoks, Dimetaklor, Difenamid, Dimetipin, Dimetilan, Dinitramin, Difenilamid, Disulfoton, Ditalimfos, Edifenfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Famfur, Fenarimol, Fenklorfos, Fenflutrin, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flupiradifuron, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutolanil, Flutriafol, Fluvalinat-tau (suma izomera), Formotion, Furilazol, gama-BHC (lindan), Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrion, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izofenfos metil, Izoksaben, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metidation, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitrilin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Pebulat, Penkonazol, Pendimetilin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Procimidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Protiokonazol, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Spiromezifen, Sulkotrin, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbamil, Tolklofos-metil, Tolilfluanid, Tolilfluanid metabolit, trans-Klordan, Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Triklornat, Tridemorf, Trietazin, Trifluralin, Unikonazol, Valifenalat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allethrin-3,4 (Bioallethrin), Allidochlor, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Atraton, Azinphos-ethyl, Benfluralin, Benzyladenine, beta-BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Binapacryl, Bitrex (Denatonium benzoate), Bromacil, Bromfeninfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Captan, Carbophenothion, Carbophenothion methyl, Carfentrazone ethyl, Chinomethionat, Chlorbenside, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapir, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolate, Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlofluanid, Dichlorvos, Diclobutrazol, Dicloran, Dicofol, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Dinitramine, Diphenamid, Diphenylamine, Disulfoton, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Famphur, Fenarimol, Fenchlorphos, Fenfluthrin, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flupyradifurone, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutolanil, Flutriafol, Fluvalinate (sum of isomers), Formothion, Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos methyl, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitalin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		DDD, p,p'-DDE, p,p'-DDT, Ofurace, Paclobutrazol, Parathion, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phorate, Phosfolan, Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prosulfocarb, Prothioconazole, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesife, Sulcotrione, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vinclozolin Ukupno / In total 259				
	A10-4 Žitarice i proizvodi od žitarica – visoki udio škroba i/ili proteina te niski udio vode i masti <i>Cereals and cereals products – high starch and/or protein content and low water and fat content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (Z)-Pirifenox, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroaniline, 2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Akrinatri (suma izomera), Alaklor, Aldrin, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Azinfos-metil, Benfluralin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitreks (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Kadusafos, Karbofenotion, Karbofenotion metil, Kinometionat, Klorbenzid, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klor-	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		fenson, Klormefos, Klorobenzilat, Kloroneb, Kloropropilat, Klorothalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Klordan, cis-Nonaklor, cis-Permetrin, Klodinafop propargil, Klokventocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklorvos, Diklobutrazol, Dikloran, Dieldrin, Dimetaklor, Dimetipin, Dimetilan, Difenilamid, Disulfoton, Ditalimfos, Edifenfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Famfur, Fenarimol, Fenklorfos, Fenitrotion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludoksionil, Fluensulfon, Flumetrin, Fluotrimazol, Flurokloridon, Flurprimidol, Flutiacet-metil, Fluvalinat-tau (suma izomera), Formotion, Furilazol, gama-BHC (lindan), Genit, Halfenproks, Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izofenfos-metil, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metakrifos, Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitalin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Pebulat, Penkonazol, Pendimetalin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosalon, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Prodiamin, Profenofos, Profluralin, Prometrin, Propaklor, Propanil, Propazin, Propetamfos, Propizoklor, Propizamid, Prosulfokarb, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Spiromezifen, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbazil, Tolklofos-metil, Tolilfluamid, Tolilfluamid metabolit, trans-Klordan, Translutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Triciklazol, Tridemorf, Trietazin, Trifluralin, Unikonazol, Vinklozolin				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		/ (E)-Metominstrobin, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin., Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Azinphos methyl, Benfluralin, beta-BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Cadusafos, Carbophenothion, Carbophenothion methyl, Chinomethionat, Chlorbenside, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chloropropylate, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolate, Cinidon ethyl, cis-Chlordan, cis-Nonachlor, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Coumaphos, Crimidine, Crotoxyphos, Cyanazine, Cyanofenphos, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlorvos, Diclobutrazol, Dicloran, Dieldrin, Dimethachlor, Dimethipin, Dimetilan, Diphenamid, Diphenylamine, Disulfoton, Ditalimfos, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flurochloridon, Flurprimidol, Fluthiacet-methyl, Fluvalinate (sum of isomers), Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos-methyl, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methacrifos, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitalin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paclobutrazol,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Parathion, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phorate, Phosfolan,, Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Prodiamine, Profenofos, Profluralin, Prometryn, Propachlor, Propanil, Propazin, Propetamphos, Propisochlor, Propyzamide, Prosulfocarb, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesifen, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transluthrin, trans-Nonachlor, trans-Permethrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Vinclozolin Ukupno / In total 256				
	A10-5 Teški matriksi <i>Difficult commodities</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroanilin, 2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Aktrinatin (suma izomera), Alaklor, Aldrin, , Alidoklor, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Azinfos metil, Benfluralin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitrex (Decanionium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromoci-	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		klen, Bromofos, Bromofos-etil, Bromopropilat, Kadusafos, Karbofenotion, Karbofenotion metil, Karfentrazon etil, Kinometionat, Klorbenzid, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klor-mefos, Klorobenzilat, Kloroneb, Klorothalonil, Klorprofam, Klorpirifos, Klor-pirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Klordan, cis-Permetrin, Klodinafop propargil, Klokvincet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklorvos, Diklobutrazol, Dikloran, Dieldrin, Dimefoks, Dimetaklor, Difena-mid, Dimetipin, Dimetilan, Difenamid, Difenilamid, Disulfoton, Ditalimfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Famfur, Fenarimol, Fenklorfos, Fenitrotion, Fenok-saprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flam-prop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutriafol, Fluvalinat-tau (suma izomera), Formotion, Fostia-zat (suma izomera), Furilazol, gama-BHC (lindan), Genit, Heptaklor, Hepta-klor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izo-fenfos metil, Izoksaben, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lena-cil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitratin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Pebulat, Penkonazol, Pendimeta-lin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Profluralin, Prometrin, Pro-paklor, Propanil, Propetamfos, Propizoklor, Propizamid, Prosulfokarb, Protio-fos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piripoksi-fen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Silti-ofam, Simazin, Spiromesifen, Sulkotripon, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflut-rin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbazil, Tolklofos-metil, Tolilfluamid, Tolilfluamid metabolit, trans-Klordan, Transflutrin, trans-				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Tridemorf, Trietazin, Trifluralin, Unikonazol, Valifenalat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allidochlor, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Azinphos methyl, Benfluralin, beta-BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Bitrex (Denatonium benzoate), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Cadusafos, Carbophenothion, Carbophenothion methyl, Carfentrazone-ethyl, Chinomethionat, Chlorbenside, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolinat, Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanofenphos, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlorvos, Diclobutrazol, Dicloran, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Diphenamid, Diphenylamine, Disulfoton, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutriafol, Fluvalinate (sum of isomers), Formothion, Fosthiazate (sum of isomers), Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos-methyl, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Methoxychlor,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paclobutrazol, Parathion, Pebulate, Penconazole, Pendimethalin, Pentachloroaniline, Pentachloroanisole, Pentachlorobenzene, Pentanochlor (Solan), Phorate, Phosfolan, Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prosulfocarb, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesifen, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, trans-Permethrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vinclozolin Ukupno / In total 254				
	A10-6 Meso i mesni proizvodi <i>Meat and meat products</i>	Metoda za određivanje ostataka pesticida u hrani životinjskog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of animal origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,005 – 0,02 mg/kg Kloroneb, Pentaklorobenzen, Teknazen, alfa-BHC, Heksaklorobenzen, Pentakloroanisol, beta-BHC, gama-BHC (Lindane), Kvintozen, Diazinon, delta-BHC, Endosulfan eter, Klorpirifos-metil, Vinklozolin, Paration metil, Heptaklor, Pirimifos metil, Pentaklorotioanisol, Malation, Aldrin, Fention, Klorpirifos, 4,4'-Diklorobenzofenon, Paration, Fenson, Izodrin, Klorfenvinfos (suma izomera), Pendimetalin, Heptaklor-egzo-epoksd (izomer b),	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-1- 2,16/03 5. izdanje / <i>edition</i> (2022-05-17)	24/01/2005	17/05/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Oksiklordan, Fipronil, Klorbenzid, Metidation, trans-Klordan, o,p'-DDE, alfa-Endosulfan (Endosulfan I), cis-Klordan, trans-Nonaklor, Profenofos, Klorfenson, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, Fention sulfone, Klorobenzilat, beta-Endosulfan (Endosulfan II), Fention sulfoksid, p,p'-DDD, o,p'-DDT, cis-Nonaklor, Endrin aldehyd, Triazofos, 4,4'-metoksiklor olefin, Endosulfan sulfat, p,p'-DDT, 2,4'-Metoksiklor, Resmetrin (suma izomera), Fosmet, Endrin keton, Bifentrin, Tetradifon, Mireks, lambda-Cihalotrin, Pirazofos, Azinfos etil, cis-Permetrin, trans-Permetrin, Ciflutrin-1, Cifluthrin-2, Cifluthrin-3, Ciflutrin-4, Cipermetrin-1, Cipermetrin-2, Fenvalerat-1, Fenvalerat-2, Deltametrin-2, Indoksakarb / Chloroneb, Pentachlorobenzene, Tecnazene, alpha-BHC, Hexachlorobenzene, Pentachloroanisole, beta-BHC, gamma-BHC (Lindane), Quintozene, Diazinon, delta-BHC, Endosulfan ether, Chlorpyrifos-methyl, Vinclozolin, Parathion methyl, Heptachlor, Pirimiphos methyl, Pentachlorothioanisole, Malathion, Aldrin, Fenthion, Chlorpyrifos, 4,4'-Dichlorobenzophenone, Parathion, Fenson, Isodrin, Chlorfenvinphos (sum of isomers), Pendimethalin, Heptachlor-exo-epoxide, Oxychlordan, Fipronil, Chlorbenside, Methidation, trans-Chlordane, o,p'-DDE, alpha-Endosulfan, cis-Chlordane, trans-Nonachlor, Profenofos, Klorfenson, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, Fenthion sulfone, Chlorobenzilate, beta-Endosulfan, Fenthion sulfoxide, p,p'-DDD, o,p'-DDT, cis-Nonachlor, Endrin aldehyde, Triazophos, 4,4'-methoxychlor olefin, Endosulfan sulfate, p,p'-DDT, 2,4'-Methoxychlor, Resmethrin (sum of isomers), Phosmet, Endrin ketone, Bifenthrin, Tetradifon, Mirex, lambda-Cyhalothrin, Pyrazophos, Azinphos ethyl, cis-Permethrin, trans-Permethrin, Cyfluthrin 1, Cyfluthrin 2, Cyfluthrin 3, Cyfluthrin 4, Cypermethrin 1, Cypermethrin 2, Fenvalerate 1, Fenvalerate 2, Deltamethrin 2, Indoxacarb Ukupno / In total 78				
	A10-7 Mlijeko i mliječni proizvodi, sladoled i puding	Metoda za određivanje ostataka pesticida u hrani životinjskog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of animal origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i>	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-1-2,16/03 5. izdanje /	01/12/2004	17/05/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	Milk and dairy products, icecream and pudding	0,0008 – 0,008 mg/kg Kloroneb, Pentaklorobenzen, Teknazen, alfa-BHC, Heksaklorobenzen, Pentakloroanisol, beta-BHC, gama-BHC (Lindane), Kvintozen, Diazinon, delta-BHC, Endosulfan eter, Klorpirifos-metil, Vinklozolin, Paration metil, Heptaklor, Pirimifos metil, Pentaklorotioanisol, Malation, Aldrin, Fention, Klorpirifos, 4,4'-Diklorobenzofenon, Paration, Fenson, Izodrin, Klorfenvinfos (suma izomera), Pendimetalin, Heptaklor-egzo-epoksd (izomer b), Oksiklordan, Fipronil, Klorbenzid, Metidation, trans-Klordan, o,p'-DDE, alfa-Endosulfan (Endosulfan I), cis-Klordan, trans-Nonaklor, Profenofos, Klorfenson, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, Fention sulfone, Klorobenzilat, beta-Endosulfan (Endosulfan II), Fention sulfoksid, p,p'-DDD, o,p'-DDT, cis-Nonaklor, Endrin aldehid, Triazofos, 4,4'-metoksiklor olefin, Endosulfan sulfat, p,p'-DDT, 2,4'-Metoksiklor, Resmetrin (suma izomera), Fosmet, Endrin keton, Bifentrin, Tetradifon, Mireks, lambda-Cihalotrin, Pirazofos, Azinfos etil, cis-Permetrin, trans-Permetrin, Ciflutrin-1, Cifluthrin-2, Cifluthrin-3, Ciflutrin-4, Cipermetrin-1, Cipermetrin-2, Fenvalerat-1, Fenvalerat-2, Deltametrin-2, Indoksakarb / Chloroneb, Pentachlorobenzene, Tecnazene, alpha-BHC, Hexachlorobenzene, Pentachloroanisole, beta-BHC, gamma-BHC (Lindane), Quintozene, Diazinon, delta-BHC, Endosulfan ether, Chlorpyrifos-methyl, Vinclozolin, Parathion methyl, Heptachlor, Pirimiphos methyl, Pentachlorothioanisole, Malathion, Aldrin, Fenthion, Chlorpyrifos, 4,4'-Dichlorobenzophenone, Parathion, Fenson, Isodrin, Chlorfenvinphos (sum of isomers), Pendimethalin, Heptachlor-exo-epoxide, Oxychlordan, Fipronil, Chlorbenside, Methidation, trans-Chlordane, o,p'-DDE, alpha-Endosulfan, cis-Chlordane, trans-Nonachlor, Profenofos, Chlorfenson, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, Fenthion sulfone, Chlorobenzilate, beta-Endosulfan, Fenthion sulfoxide, p,p'-DDD, o,p'-DDT, cis-Nonachlor, Endrin aldehyde, Triazophos, 4,4'-methoxychlor olefin, Endosulfan sulfate, p,p'-DDT, 2,4'-Methoxychlor, Resmethrin (sum of isomers), Phosmet, Endrin ketone, Bifenthrin, Tetradifon, Mirex, lambda-Cyhalothrin, Pyrazophos, Azinphos ethyl, cis-Permethrin, trans-Permethrin, Cyfluthrin 1, Cyfluthrin 2, Cyfluthrin 3, Cyfluthrin 4, Cypermethrin 1, Cypermethrin 2, Fenvalerate 1, Fenvalerate 2, Deltamethrin 2, Indoxacarb		edition (2022-05-17)		

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Ukupno / In total 78				
	A10-8 Jaja Eggs	Metoda za određivanje ostataka pesticida u hrani životinjskog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of animal origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,005 – 0,012 mg/kg Kloroneb, Pentaklorobenzen, Teknazen, alfa-BHC, Heksaklorobenzen, Pentakloroanisol, beta-BHC, gama-BHC (Lindane), Kvintozen, Diazinon, delta-BHC, Endosulfan eter, Klorpirifos-metil, Vinklozolin, Paration metil, Heptaklor, Pirimifos metil, Pentaklorotioanisol, Malation, Aldrin, Fention, Klorpirifos, 4,4'-Diklorobenzofenon, Paration, Fenson, Izodrin, Klorfenvinfos (suma izomera), Pendimetalin, Heptaklor-egzo-epoksd (izomer b), Oksiklordan, Fipronil, Klorbenzid, Metidation, trans-Klordan, o,p'-DDE, alfa-Endosulfan (Endosulfan I), cis-Klordan, trans-Nonaklor, Profenofos, Klorfenson, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, Fention sulfone, Klorobenzilat, beta-Endosulfan (Endosulfan II), Fention sulfoksid, p,p'-DDD, o,p'-DDT, cis-Nonaklor, Endrin aldehid, Triazofos, 4,4'-metoksiklor olefin, Endosulfan sulfat, p,p'-DDT, 2,4'-Metoksiklor, Resmetrin (suma izomera), Fosmet, Endrin keton, Bifentrin, Tetradifon, Mireks, lambda-Cihalotrin, Pirazofos, Azinfos etil, cis-Permetrin, trans-Permetrin, Ciflutrin-1, Cifluthrin-2, Cifluthrin-3, Ciflutrin-4, Cipermetrin-1, Cipermetrin-2, Fenvalerat-1, Fenvalerat-2, Deltametrin-2, Indoksakarb / <i>Chloroneb, Pentachlorobenzene, Tecnazene, alpha-BHC, Hexachlorobenzene, Pentachloroanisole, beta-BHC, gamma-BHC (Lindane), Quintozene, Diazinon, delta-BHC, Endosulfan ether, Chlorpyrifos-methyl, Vinclozolin, Parathion methyl, Heptachlor, Pirimiphos methyl, Pentachlorothioanisole, Malathion, Aldrin, Fenthion, Chlorpyrifos, 4,4'-Dichlorobenzophenone, Parathion, Fenson, Isodrin, Chlorfenvinphos (sum of isomers), Pendimethalin, Heptachlor-exo-epoxide, Oxychlordan, Fipronil, Chlorbenside, Methidation, trans-Chlordane, o,p'-DDE, alpha-Endosulfan, cis-Chlordane, trans-</i>	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-1- 2,16/03 5. izdanje / <i>edition</i> (2022-05-17)	15/06/2021	17/05/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Nonachlor, Profenofos, Chlorfenson, p,p'-DDE, Dieldrin, o,p'-DDD, Endrin, Fenthion sulfone, Chlorobenzilate, beta-Endosulfan, Fenthion sulfoxide, p,p'-DDD, o,p'-DDT, cis-Nonachlor, Endrin aldehyde, Triazophos, 4,4'-methoxychlor olefin, Endosulfan sulfate, p,p'-DDT, 2,4'-Methoxychlor, Resmethrin (sum of isomers), Phosmet, Endrin ketone, Bifenthrin, Tetradifon, Mirex, lambda-Cyhalothrin, Pyrazophos, Azinphos ethyl, cis-Permethrin, trans-Permethrin, Cyfluthrin 1, Cyfluthrin 2, Cyfluthrin 3, Cyfluthrin 4, Cypermethrin 1, Cypermethrin 2, Fenvalerate 1, Fenvalerate 2, Deltamethrin 2, Indoxacarb</i> Ukupno / In total 78				
A11	A11-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, podzemne vode, vode za kupanje, otpadne, procjedne i površinske vode, eluati otpada <i>Water for human consumption, natural spring water, natural mineral water, table water, ground water, bathing water, waste water, surface water and leachate water, waste eluates</i>	Određivanje žive <i>Determination of mercury</i> 0,007 µg/L – vode, eluati otpada / water, leachate water 0,001 mg/kg – mulj, sediment, biota, tlo, otpad / Sludge, sediment, biota, soil, waste 0,010 mg/kg – kruta oporabljena goriva / solid recovered fuels	AAS (FIMS)	HRN EN ISO 12846:2012 (ISO 12846:2012; EN ISO 12846:2012)	13/09/2013	23/02/2022
	A11-2 Mulj, sediment, biota, tlo, otpad,					

Oznaka/ Identification	Materijali/ Proizvodi <i>Materials/ Products</i>	Vrsta ispitivanja/Svojstvo <i>Type of test / Property</i> Granica kvantifikacije / <i>Limit of quantification</i>	Tehnika / <i>Technique</i>	Metoda ispitivanja <i>Test method</i>	U primjeni od/ <i>In use from</i>	Datum posljednje izmjene/ <i>Date of last change</i>																																												
	kruta oporabljena goriva <i>Sludge, sediment, biota, soil, waste, solid recovered fuels</i>																																																	
A12	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, površinske i podzemne vode <i>Water for human consumption, natural spring water, natural mineral water, table water, surface water, and ground water</i>	Određivanje odabranih elemenata metodom ICP-MS <i>Determination of selected elements by ICP-MS method</i> µg/l <table><tr><th>Element/ <i>Element</i></th><th>Granica kvantifikacija / <i>Limit of quantification</i></th><th>Element/ <i>Element</i></th><th>Granica kvantifikacija / <i>Limit of quantification</i></th></tr><tr><td>Al</td><td>0,360</td><td>Ni</td><td>0,065</td></tr><tr><td>As</td><td>0,042</td><td>Pb</td><td>0,005</td></tr><tr><td>Ba</td><td>0,039</td><td>Sb</td><td>0,002</td></tr><tr><td>Be</td><td>0,042</td><td>Se</td><td>0,034</td></tr><tr><td>Cd</td><td>0,005</td><td>Sn</td><td>0,037</td></tr><tr><td>Co</td><td>0,004</td><td>Te</td><td>0,005</td></tr><tr><td>Cr</td><td>0,030</td><td>Tl</td><td>0,001</td></tr><tr><td>Cu</td><td>0,046</td><td>V</td><td>0,007</td></tr><tr><td>Mn</td><td>0,046</td><td>Zn</td><td>0,300</td></tr><tr><td>Mo</td><td>0,002</td><td></td><td></td></tr></table>	Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>	Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>	Al	0,360	Ni	0,065	As	0,042	Pb	0,005	Ba	0,039	Sb	0,002	Be	0,042	Se	0,034	Cd	0,005	Sn	0,037	Co	0,004	Te	0,005	Cr	0,030	Tl	0,001	Cu	0,046	V	0,007	Mn	0,046	Zn	0,300	Mo	0,002			ICP-MS	HRN EN ISO 17294-2:2016 <i>(ISO 17294-2:2016; EN ISO 17294-2:2016)</i>	20/09/2019	11/05/2021
Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>	Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>																																															
Al	0,360	Ni	0,065																																															
As	0,042	Pb	0,005																																															
Ba	0,039	Sb	0,002																																															
Be	0,042	Se	0,034																																															
Cd	0,005	Sn	0,037																																															
Co	0,004	Te	0,005																																															
Cr	0,030	Tl	0,001																																															
Cu	0,046	V	0,007																																															
Mn	0,046	Zn	0,300																																															
Mo	0,002																																																	
A13	Hrana <i>Food</i>	Određivanje odabranih elemenata metodom ICP-MS nakon digestije nitratnom kiselinom <i>Determination of selected elements by ICP-MS method after nitric acid digestion</i> 0,01 mg/kg As, Cd, Pb	ICP-MS	HRN EN 15763:2010 <i>(EN 15763:2009)</i>	27/10/2020	25/02/2022																																												
A14	Hrana i hrana za životinje <i>Food and animal</i>	Određivanje odabranih elemenata u hrani i stočnoj hrani optičkom emisijskom spektrometrijom induktivno vezane	ICP-OES	HRN EN 16943:2017 <i>(EN 6943:2017)</i>	18/04/2016	04/05/2022																																												

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	feeding stuffs	plazme (ICP-OES) <i>Determination of sodium and calcium content in food and animal feeding stuffs by inductively coupled plasma optical emission (ICP-OES)</i> 83 mg/kg Ca, Na				
A15	A15-1 Hrana biljnog podrijetla (uključujući mlinske i pekarske proizvode, tjesteninu i brzo smrznuta tijesta) <i>Food of plant origin (including mill and bakery products, pasta and quick- frozen dough)</i>	Određivanje mikotoksina metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of mycotoxins - liquid chromatography with mass spectrometry (LC-MS/MS)</i> 1 µg/kg Aflatoksin B1, B2, G1, G2, HT2- toksin, T2-toksin, ohratoksin A/aflatoxin B1, B2, G1, G2, HT2-toxin, T2-toxin, ochratoxin 10 µg/kg Zearalenon/zearalenone 50 µg/kg Deoksinivalenol/deoxynivalenol	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 02,41/203 IV. izdanje / edition (2022-02-21)	07/09/2018	21/02/2022
	A15-2 Mlijeko Milk	Određivanje mikotoksina metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of mycotoxins - liquid chromatography with mass spectrometry (LC-MS/MS)</i> 0,02 µg/kg Aflatoksin M1/aflatoxin M1	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 02,41/203 IV. izdanje / edition (2022-02-21)	21/02/2022	21/02/2022
A16	Tlo, sediment / Soil, sediment	Određivanje pesticida metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of pesticide - liquid chromatography with mass spectrometry (LC-MS/MS)</i>	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 38,39/200 IV. izdanje /	10/09/2018	14/03/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		0,002 mg/kg Atrazin, simazin, kinoksifen / Atrazine, simazine, quinoxifen		edition (2022-03-14)		
A17	Tlo, sediment / Soil, sediment	Određivanje pesticida metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Determination of pesticide - gas chromatography with mass spectrometry (GC-MS/MS)</i> 0,005 mg/kg Kloroneb, pentaklorbenzen (PeCB), alfa-BHC (alfa-HCH), heksaklorbenzen (HCB), pentakloroanisol, beta-BHC (beta-HCH), gama-BHC (gama-HCH), delta-BHC (delta-HCH), endosulfan eter, heptaklor, pentaklorotioanisol, aldrin, 4,4-diklorobenzofenon, fenson, izodrin, heptaklor epoksid (izomer B), klorbenzid, trans-klordan, 2,4-DDE (o,p-DDE), endosulfan I (endosulfan alfa), cis-klordane, trans-nonaklor, klorfenson (Ovex), dieldrin, 4,4-DDE (p,p-DDE), 2,4-DDD (o,p-DDD), endrin, endosulfan II (endosulfan beta), 4,4-DDD (p,p-DDD), 2,4-DDT (o,p-DDT), cis-nonaklor, endrin aldehid, 4,4-metoksiklor olefin, endosulfan sulfat, 4,4-DDT (p,p-DDT), 2,4-metoksiklor, endrin keton, tetradifon, mirex, dikofol, cipermetrin / <i>Chloroneb, Pentachlorobenzene (PeCB), alpha-BHC (alpha-HCH), Hexachlorobenzene (HCB), Pentachloroanisole, beta-BHC (beta-HCH), gamma-BHC (gamma-HCH), delta-BHC (delta-HCH), Endosulfan ether, Heptachlor, Pentachlorothioanisole, Aldrin, 4,4-dichlorobenzophenone, Fenson, Isodrin, Heptachlor epoxide (isomer B), Chlorbenside, trans-chlordane, 2,4-DDE (o,p-DDE), Endosulfan I (Endosulfan alpha), cis-Chlordane, trans-Nonachlor, Chlorfenson (Ovex), Dieldrin, 4,4-DDE (p,p-DDE), 2,4-DDD (o,p-DDD), Endrin, Endosulfan II (Endosulfan beta), 4,4-DDD (p,p-DDD), 2,4-DDT (o,p-DDT), cis-Nonachlor, Endrin aldehyde, 4,4-methoxychlor olefin, Endosulfan sulfate, 4,4-DDT (p,p-DDT), 2,4-methoxychlor, Endrin ketone, Tetradifon, Mirex, Dicofol, Cipermetryn</i> Ukupno / In total 41	GC-MS/MS	Vlastita metoda In-house method SOP-LEK- 38,39/03 III. izdanje / edition (2022-04-13)	15/01/2013	13/04/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
A18	A18-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje matrina u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC- MS/MS) <i>Determination of Matrine in food of plant origin by liquid chromatography- mass spectrometry (LC-MS/MS)</i> 0,005 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 41/229 II. izdanje / <i>edition</i> (2021-06-10)	15/09/2020	10/06/2021
A19	A19-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> A19-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water</i>	Određivanje maleinskog hidrazida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Maleic-hydrazide in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,05 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 41/234 I. izdanje / <i>edition</i> (2021-04-20)	20/04/2021	20/04/2021
A20	A20-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje ditianona u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC- MS/MS) <i>Determination of Dithianon by liquid chromatography-mass spectrometry (LC- MS/MS)</i>	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 41/220 II. izdanje / <i>edition</i> (2021-11-03)	11/09/2019	03/11/2021

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	A20-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> A20-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i> A20-4 Žitarice i proizvodi od žitarica – visoki udio škroba i/ili proteina te niski udio vode i masti <i>Cereals and cereals products – high starch and/or protein content and low water and fat content</i>	0,01 mg/kg				
A-21	A21-1 Voće i povrće s	Određivanje klormekvata i mepikvata u hrani biljnog	LC-MS/MS	Vlastita metoda	21/11/2017	21/11/2017

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	visokim udjelom vode <i>Fruits and vegetables – high water content</i>	podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Chlormequat and Mepiquat by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg		<i>In-house method</i> SOP-LKH-44/189 I. izdanje / edition (2017-11-21) QuPPE-PO, V12		
A-22	A22-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> A22-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Određivanje Etefona u hrani bilnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Ethephon in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/210 I. izdanje / edition (2018-11-22) QuPPE-PO, V12	22/11/2018	22/11/2018
A-23	A23-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje glifosata u hrani bilnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Glyphosate by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/221 I. izdanje / edition (2019-07-16)	16/07/2019	16/07/2019

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
				QuPPe-PO, V12		
A-24	A24-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> A24-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> A24-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	Određivanje fosetila i fosfonske kiseline QuPPe metodom u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Fosetyl and Phosphonic acid following QuPPe method in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fosetil/Fosetyl 0,02 mg/kg Fosfonska kiselina/Phosphonic acid	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/236 I. izdanje / <i>edition</i> (2021-12-02) QuPPe-PO, V12	02/12/2021	02/12/2021
A-25	A25-1 Voće i povrće s visokim udjelom vode	Određivanje morfolina, dietanolamina i trietanolamina metodom tekućinske kromatografije s masenom spektrometrijom LC-MS/MS <i>Determination of morpholine, diethanolamine and triethanolamine by liquid</i>	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	28/02/2022	28/02/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	<i>Fruits and vegetables – high water content</i> A25-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> A25-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	<i>chromatography with LC-MS / MS mass spectrometry</i> 0,01 mg/kg				
A-26	A26-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> A26-2 Voće i povrće s visokim udjelom	Određivanje klorata i perklorata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Chlorate and Perchlorate in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH- 41/247 I. izdanje / <i>edition</i> (2022-12-02) <i>QuPPe-PO, V12</i>	02/12/2022	02/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>					
A-27	A27-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> A27-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Određivanje parakvata i dikvata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Paraquat and Diquat in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/245 I. izdanje / <i>edition</i> (2022-12-02) <i>QuPPE-PO, V12</i>	02/12/2022	02/12/2022

Pripreme analitičkih uzorka / Preparation of analytical samples:

- Tlo / Soil: modificirana / modified HRN EN ISO 54321:2021 (ISO 54321:2020; EN ISO 54321:2021); HRN ISO 14869-2:2004 (ISO 14869-2:2002)
- Otpad / Waste: modificirana / modified HRN EN 13657:2008 (EN 13657:2002); modificirana / modified HRN EN 12457-4:2005 (EN 12457-4:2002)
- Muljevi, sedimenti, biota / Sludges, sediments, biota: modificirana / modified HRN EN ISO 54321:2021 (ISO 54321:2020; EN ISO 54321:2021)
- Kruta oporabljena goriva / Solid recovered fuels: HRN EN 15411:2011 (EN 15411:2011)
- Mineralne sirovine / Mineral raw materials: HRN ISO 14869-2:2004 (ISO 14869-2:2002)
- Hrana / Food: HRN EN 13805:2014 (EN 13805:2014)
- Hrana biljnog podrijetla / Food of plant origin: HRN EN 15662 (EN 15662:2018)
- Hrana biljnog podrijetla / Food of plant origin: HRN EN 15662 modificirana sa alkalnom hidrolizom / modified (EN 15662:2018) with alkaline hydrolisis
- Hrana biljnog podrijetla / Food of plant origin: QuPPE-PO-Method



Napomene / Notes:

ICP-OES - optička emisijska spektrometrija induktivno vezane plazme / ICP OES - *inductively coupled plasma atomic emission spectroscopy*

AAS (FIMS) – “Flow Injecton Mercury System” atomska apsorpcijska spektrometrija / AAS (FIMS) – *flow injecton mercury system atomic absorption spectroscopy*

GC-MS – plinska kromatografija-spektrometrija masa / GC-MS - *gas chromatography-mass spectrometry*

GC-MS/MS – plinska kromatografija / tandemna spektrometrija masa / GC-MS/MS - *Gas chromatography-tandem mass spectrometry*

LC-MS/MS – tekućinska kromatografija / tandemna spektrometrija masa / LC-MS/MS – *Liquid chromatography-tandem mass spectrometry*

ICP-MS - masena spektrometrija induktivno vezane plazme / ICP - MS - *inductively coupled plasma / mass spectrometry*

Početak fleksibilnog područja akreditacije 2012.godina. / The beginning of flexible scope of accreditation 2012.

Praćenje izmjena

Datum izmjene	Opis izmjene
02/08/2019	Izmjena metoda A1 i A2 Određivanje odabranih elemenata . Promjena granica kvantifikacije za Be, Mn, P, Sr, Te, Ti, Tl i U.
03/03/2020	Izmjena metode A4 Određivanje odabranih elemenata nakon digestije kiselinama Proširenje metode po matrici: A4-2 Određivanje odabranih elemenata nakon digestije obrnutom zlatotopkom (<i>kruta uporabljena goriva</i>)
16/06/2022	Izmjena metode A6 Određivanje hlapljivih organskih spojeva metodom analize para iznad otopine (headspace) plinskom kromatografijom s masenom spektrometrijom (GC-MS/MS) Proširenje metode po matrici: A6-2 Određivanje BTEX-a i drugih lakohlapivih spojeva metodom Headspace GC-MS/MS (<i>tlo, otpad, mulj i sediment</i>)
06/06/2022	Izmjena metode A8 Određivanje odabranih sredstava za zaštitu bilja u vodi metodom tekućinske kromatografije s masenom spektrometrijom Proširenje metode po parametrima: A8-2 Određivanje kiselih pesticida (Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode i eluati otpada)
23/02/2022	Proširenje metode A-11 Određivanje sadržaja žive po matricama: - A11-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, podzemne vode, vode za kupanje, otpadne, procjedne i površinske vode, eluati otpada - A11-2 Biot, kruta uporabljena goriva
26/07/2022	Metode prebačene iz fiksnog u fleksibilno područje: - A-12 Određivanje odabranih elemenata metodom ICP-MS (Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, površinske i podzemne vode) - A-13 Određivanje odabranih elemenata metodom ICP-MS nakon digestije nitratnom kiselinom (Hrana) - A-14 Određivanje odabranih elemenata u hrani i stočnoj hrani optičkom emisijskom spektrometrijom induktivno vezane plazme (ICP-OES) (Hrana i hrana za životinje) - A-15 Određivanje mikotoksina LC-MS/MS tehnikom: Aflatoksin B1, B2, G1, G2, M1, HT2- toksin, T2-toksin, zearalenon, ohratoksin A, deoksinivalenol (Hrana) - A-16 Određivanje pesticida u tlu i sedimentu LC-MS/MS tehnikom - A-17 Određivanje ostataka pesticida u krutim matricama u okolišu metodom plinske kromatografije s masenom spektrometrijom GC-MS/MS

	• A-18 Određivanje matrina u hrani biljnog podrijetla metodom tekućinske kromatografije s masenim spektrometrom (LC-MS/MS) • A-19 Određivanje maleinskog hidrazida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenim spektrometrom (LC-MS/MS) • A-20 Određivanje ditianona u hrani biljnog podrijetla metodom tekućinske kromatografije s masenim spektrometrom (LC-MS/MS) • A-21 Određivanje klormekvata i mepikvata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenim spektrometrom (LC-MS/MS) • A-22 Određivanje Etefona u hrani biljnog podrijetla metodom tekućinske kromatografije s masenim spektrometrom (LC-MS/MS) • A-23 Određivanje glifosata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenim spektrometrom (LC-MS/MS) • A-24 Metoda za određivanje fosetila i fosfonske kiseline QuPPE metodom u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) • A-25 Određivanje morfolina, dietanolamina i trietanolamina metodom tekućinske kromatografije s masenom spektrometrijom LC-MS/MS
22/03/2023	Proširenje metoda: • A6 Određivanje hlapljivih organskih spojeva metodom analize para iznad otopine (headspace) plinskom kromatografijom s masenom spektrometrijom (GC-MS/MS) Proširenje metode po pokazateljima i matrici: <i>A6-2 Sediment</i> • A9 Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A Proširenje metode po pokazateljima i matrici: A9-5 Teški matriks • A10 Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) - proširenje metode po pokazateljima Metode dodane u fleksibilno područje: • A-26 Određivanje klorata i perklorata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) • A-27 Određivanje parakvata i dikvata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS)