

Popis ispitnih metoda u fleksibilnom području akreditacije, status 2026-08-24
The valid list of the test methods in the flexible scope of accreditation, status 2026-08-24
Ovaj popis odnosi se na Prilog potvrdi o akreditaciji broj: 1073 koji je dostupan na www.akreditacija.hr
This list is related to Annex to Accreditation Certificate Number: 1073 which is available at www.akreditacija.hr
Izmjene su vidljive na kraju dokumenta / Changes visible at the end of the document.
FLEKSIBILNO PODRUČJE AKREDITACIJE / FLEXIBLE SCOPE OF ACCREDITATION

Oznaka/ Identification	Materijali/ Proizvodi <i>Materials/ Products</i>	Vrsta ispitivanja/Svojstvo <i>Type of test / Property</i> Granica kvantifikacije / <i>Limit of quantification</i>				Tehnika / <i>Technique</i>	Metoda ispitivanja <i>Test method</i>	U primjeni od/ <i>In use from</i>	Datum posljednje izmjene/ <i>Date of last change</i>
A1	A1-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, površinske i podzemne vode <i>Water for human consumption, natural spring water, natural mineral water, table water, bathing water, surface water and ground water</i>	Određivanje odabranih elemenata <i>Determination of selected elements</i>				ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	03/01/2012	14/08/2026
		mg/l							
		Element/ <i>Element</i>	Granica kvantifikacije / <i>Limit of quantification</i>	Element/ <i>Element</i>	Granica kvantifikacije / <i>Limit of quantification</i>				
		Ag	0,001	Mg	0,075				
		Al	0,009	Mn	0,019				
		B	0,012	Na	0,084				
		Ba	0,020	Si	0,046				
		Ca	0,015	Sr	0,017				
		Fe	0,003	Ti	0,001				
	K	0,061	V	0,012					
	A1-2 Otpadne i procjedne vode <i>Waste water and, leachate water</i>	Određivanje odabranih elemenata <i>Determination of selected elements</i>				ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	03/01/2012	17/08/2026
		mg/l							
		Element/ <i>Element</i>	Granica kvantifikacije / <i>Limit of quantification</i>	Element/ <i>Element</i>	Granica kvantifikacije / <i>Limit of quantification</i>				
		Ag	0,001	Ni	0,006				
		Al	0,010	P	0,028				

Oznaka/ Identification	Materijali/ Proizvodi <i>Materials/ Products</i>	Vrsta ispitivanja/Svojstvo <i>Type of test / Property</i> Granica kvantifikacije / <i>Limit of quantification</i>				Tehnika / Technique	Metoda ispitivanja <i>Test method</i>	U primjeni od/ <i>In use from</i>	Datum posljednje izmjene/ <i>Date of last change</i>																																				
		As	0,014	Sb	0,010																																								
		B	0,006	Se	0,079																																								
		Ba	0,006	Si	0,046																																								
		Be	0,001	Sn	0,031																																								
		Cd	0,004	Sr	0,029																																								
		Co	0,004	Te	0,019																																								
		Cr	0,002	Ti	0,003																																								
		Cu	0,004	Tl	0,031																																								
		Fe	0,008	U	0,035																																								
		Mn	0,001	V	0,015																																								
		Mo	0,007	Zn	0,006																																								
		Pb	0,022																																										
A2	Otpadne i procjedne vode <i>Waste water and leachate water</i>	Određivanje žive <i>Determination of mercury</i> 0,007 µg/L				AAS (FIMS)	HRN EN ISO 12846:2012 <i>(ISO 12846:2012; EN ISO 12846:2012)</i>	01/06/2015	08/10/2024																																				
A3	A3-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, površinske i podzemne vode <i>Water for human consumption, natural spring water, natural mineral water, table water, surface water, and ground water</i>	Određivanje odabranih elemenata metodom ICP-MS <i>Determination of selected elements by ICP-MS method</i> µg/l <table><tr><th>Element/ <i>Element</i></th><th>Granica kvantifikacija / <i>Limit of quantification</i></th><th>Element/ <i>Element</i></th><th>Granica kvantifikacija / <i>Limit of quantification</i></th></tr><tr><td>Al</td><td>0,687</td><td>Pb</td><td>0,058</td></tr><tr><td>As</td><td>0,096</td><td>Sb</td><td>0,080</td></tr><tr><td>Ba</td><td>0,166</td><td>Se</td><td>0,292</td></tr><tr><td>Be</td><td>0,021</td><td>Sn</td><td>0,055</td></tr><tr><td>Cd</td><td>0,015</td><td>Te</td><td>0,053</td></tr><tr><td>Co</td><td>0,014</td><td>Tl</td><td>0,013</td></tr><tr><td>Cr</td><td>0,058</td><td>U</td><td>0,147</td></tr><tr><td>Cu</td><td>0,096</td><td>V</td><td>0,138</td></tr></table>				Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>	Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>	Al	0,687	Pb	0,058	As	0,096	Sb	0,080	Ba	0,166	Se	0,292	Be	0,021	Sn	0,055	Cd	0,015	Te	0,053	Co	0,014	Tl	0,013	Cr	0,058	U	0,147	Cu	0,096	V	0,138	ICP-MS	HRN EN ISO 17294-2:2023 <i>(ISO 17294- 2:2023; EN ISO 17294-2:2023)</i>	20/09/2019	21/08/2026
Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>	Element/ <i>Element</i>	Granica kvantifikacija / <i>Limit of quantification</i>																																										
Al	0,687	Pb	0,058																																										
As	0,096	Sb	0,080																																										
Ba	0,166	Se	0,292																																										
Be	0,021	Sn	0,055																																										
Cd	0,015	Te	0,053																																										
Co	0,014	Tl	0,013																																										
Cr	0,058	U	0,147																																										
Cu	0,096	V	0,138																																										

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Mn	0,092	Zn	0,609				
		Mo	0,446	Hg	0,013				
		Ni	0,054						
	A3-2 Otpadne i procjedne vode Waste water and, leachate water	Određivanje odabranih elemenata metodom ICP-MS Determination of selected elements by ICP-MS method				ICP-MS	HRN EN ISO 17294-2:2023 (ISO 17294- 2:2023; EN ISO 17294-2:2023)	06/08/2026	21/08/2026
		mg/l							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		Al	0,0009	Ni	0,0002				
		As	0,0004	Pb	0,0001				
		Ba	0,0175	Sb	0,0001				
		Be	0,0001	Se	0,0006				
		Cd	0,0001	Sn	0,0001				
		Co	0,0001	Te	0,0001				
		Cr	0,0002	Tl	0,0001				
		Cu	0,0005	V	0,0004				
		Mn	0,0001	Zn	0,0027				
		Mo	0,0001	Hg	0,0001				
	A3-3 Eluat otpada Waste leachate	Određivanje odabranih elemenata metodom ICP-MS Determination of selected elements by ICP-MS method				ICP-MS	HRN EN ISO 17294-2:2023 (ISO 17294- 2:2023; EN ISO 17294-2:2023)	06/08/2026	21/08/2026
		mg/kg suhe tvari							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		Al	0,009	Ni	0,002				
		As	0,004	Pb	0,001				
		Ba	0,018	Sb	0,001				
		Be	0,001	Se	0,006				
		Cd	0,001	Sn	0,001				
		Co	0,001	Te	0,001				
		Cr	0,002	Tl	0,001				
		Cu	0,005	V	0,004				
		Mn	0,001	Zn	0,003				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Mo	0,001	Hg	0,001				
A4	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode, eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more) Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)	Određivanje hlapljivih organskih spojeva metodom analize para iznad otopine (headspace) plinskom kromatografijom s masenom spektrometrijom (GC-MS/MS) <i>Determination of volatile organic compounds in water by headspace method gas chromatography / mass spectrometry (GC-MS/MS)</i>				GC-MS/MS	HRN EN ISO 10301:2002 (ISO 10301:1997; EN ISO 10301:1997) HRN ISO 11423-1:2002 (ISO 11423-1:1997)	23/07/2014	10/01/2023
		ug/L							
		Spoj /Compound		Granica kvantifikacija / Limit of quantification					
		kloroform / chloroform		0,09					
		tetraklorugljik / carbon tetrachloride		0,06					
		1,2-dikloreten / 1,2-dichloroethane		0,05					
		1,1,1-trikloreten / 1,1,1-trichloroethane		0,03					
		trikloreten (1,1,2-trikloreten) / trichloroethene (1,1,2- trichloroethene)		0,02					
		tetrakloreten (1,1,2,2-tetrakloreten) / tetrachloroethene (1,1,2,2- tetrachloroethene)		0,04					
		benzen / benzene		0,04					
		toluen / toluene		0,02					
		ksileni / xylene		0,04					
		etilbenzen / ethylbenzene		0,02					
		1,2,4- triklorbenzen / 1,2,4-trichlorobenzene		0,02					
		1,2,3-triklorbenzen / 1,2,3-trichlorobenzene		0,02					
		1,3,5-triklorbenzen / 1,3,5-trichlorobenzene		0,02					
		bromdiklormetan / bromodichloromethane		0,04					
		dibromklormetan / dibromochloromethane		0,05					
		bromoform / bromoform		0,04					
		diklormetan / dichloromethane		0,07					
		1,1-dikloreten / 1,1-dichloroethene		0,04					
		stiren / styrene		0,05					
		heksaklorbutadien / hexachlorobutadiene		0,04					
		vinil-klorid / vinyl-chloride		0,03					
		1,2,4-trimetilbenzen /1,2,4-trimethylbenzene		0,04					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification		Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		1,3,5-trimetilbenzen / 1,3,5-trimethylbenzene	0,02				
		tetrahidrofuran / tetrahydrofuran	0,02				
		tetrahidrotiofen / tetrahydrothiophen	0,04				
		1,1-dikloretnan / 1,1-dichloroethane	0,03				
		cis-1,2-dikloreten / cis-1,2-dichloroethene	0,03				
		trans-1,2-dikloreten / trans-1,2-dichloroethene	0,05				
		1,1,2-trikloretnan / 1,1,2-trichloroethane	0,04				
		1,1,2,2-tetrakloretnan / 1,1,2,2-tetrachloroethane	0,04				
		heksakloretnan / hexachloroethane	0,03				
		epiklorhidrin / epichlorohydrin	0,02				
		n-heksan / n-hexane	0,02				
		1,2-dibrometan / 1,2-dibromethane	0,05				
	Voda Water	Određivanje prioriternih tvari u vodama metodom plinske kromatografije s masenom spektrometrijom <i>Determination of priority substances in water – gas chromatographic method with mass spectrometry</i>		GC-MS/MS			
A5	A5-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode, eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more)	Određivanje pesticida u vodama metodom plinske kromatografije s masenom spektrometrijom <i>Determination of pesticide residues in water – gas chromatographic method with mass spectrometry</i> 0,0001 µg/l Cibutrin (irgarol), diklorvos / <i>cibutrin (irgarol), dichlorvos</i> 0,0006 µg/l Heptaklor i heptaklorepoksid / <i>heptachlor and heptachlor epoxide</i> 0,0007 µg/l Dikofol / <i>Dicofol</i> 0,001 µg/l		GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31-33,37/181 V. izdanje / <i>edition</i> (2020-05-28)	15/11/2017	31/05/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)	Kloroneb, alfa-BHC (alfa-HCH), heksaklorbenzen (HCB), pentakloroanisol, beta-BHC (beta-HCH), gama-BHC (gama-HCH ili Lindan), delta-BHC (delta-HCH), klorbenzid, endosulfan eter, pentaklorotioanisol, aldrin, 4,4'-diklorobenzofenon, fenson, izodrin, cis-klordan, trans-klordan, o,p'-DDE (2,4-DDE), alfa-endosulfan (endosulfan I), klorfenson, p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), cis-nonaklor, trans-nonaklor, endrin aldehyd, 4,4'-metoksiklor olefin, endosulfan sulfat, p,p'-DDT (4,4-DDT), 2,4'-metoksiklor (metoksiklor), tetradifon, mireks, pentaklorbenzen, endrin keton, etilan (pertan) / <i>Chloroneb, alph-BHC (alpha-HCH), hexachlorobenzene (HCB), pentachloroanisole, beta-BHC (beta-HCH), gamma-BHC (gamma-HCH or Lindane), chlorbenside, endosulfan ether, pentachlorothioanisole, aldrin, 4,4'-dichlorobenzophenone, fenson, isodrin, cis-chlordane, trans-chlordane, o,p'-DDE (2,4-DDE), alpha-endosulfan (endosulfan I), chlorfenson (Ovex), p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), cis-nonachlor, trans-nonachlor, endrin aldehyde, 4,4'-metoxychlor olefin, endosulfan sulfate, p,p'-DDT (4,4-DDT), 2,4'-metoxychlor (metoxychlor), tetradifon, mirex, pentachlorobenzene, endrin ketone, ethylan (perthane)</i> 0,002 µg/l Bifenoks / Bifenox 0,008 µg/l Aklonifen / Aclonifen 0,01 µg/l Klorpirifos metil, pirimifos metil, klorpirifos etil (klorpirifos), kvinalfos, piridafention, EPN, fosalon, azinfos metil, azinfos etil, isazofos, pirimifos etil, difenilamin, kvintozen, diklofluanid, tolilfluanid, pendimetalin, nitrofen, nitalin, terbufos, fenklorfos, fention, bromofos etil, protiofos, etion, klortiofos, sulprofos, karbofenotion, leptofos, forat, fonofos, paration, triazofos, piperonil butoksid, alidoklor, etridiazol, pebulat, propaklor, cikloat, di-alat I i II (cis i trans), klomazon, terbutilazin, propizamid, pirimetanil, teflutrin, tri-alat, dimetaklor, acetoklor, transflutrin, metolaklor, linuron, triadimefon, MGK-264 I i II, difenamid, ciprodinil, penkonazol, triadimenol, procimidon, triflumizol,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		paklobutrazol, flutriafol, flutolanil, fludioksonil, pretilaklor, oksadiazon, miklobutanil, bupirimat, klorfenapir, heksazinon, tebukonazol, resmetrin I i II, fenpropatrin, tebufenpirad, fenotrin I i II (cis i trans), piriproksifen, akrinatrin I i II, fenarimol, trans-permetrin, flukvinkonazol, ciflutrin I – IV, piridaben, cipermetrin I i II, flucitrat I i II, fluvalinat I i II, deltametrin, tributilfosfat (TBP), trikloroetilfosfat (TCEP), trikloropropilfosfat (TCPP) / <i>Chlorpyrifos methyl, pirimiphos methyl, chlorpyrifos ethyl (chlorpyrifos), quinalphos, pyridaphenthion, EPN, phosalone, azinphos methyl, azinphos ethyl, isazophos, pirimiphos ethyl, diphenylamine, quintozone, dichlofluanid, tolylfluanid, pendimethalin, nitrofen, nitalin, terbufos, fenchlorphos, fenthion, bromophos ethyl, prothiofos, ethion, chlorthiophos, sulprofos, carbophenothion, leptophos, phorate, fonofos, parathion, triazophos, piperonyl butoxide, allidochlor, etridiazole, pebulate, propachlor, cycloate, diallate I i II (cis i trans), clomazone, terbutylazine, propyzamide, pyrimethanil, tefluthrin, triallate, dimethachlor, acetochlor, transfluthrin, metolachlor, linuron, triadimefon, MGK-264 I i II, diphenamid, cyprodinil, penconazole, triadimenol, procymidone, triflumizole, paclobutrazol, flutriafol, flutolanil, fludioxonil, pretilachlor, oxadiazon, myclobutanil, bupirimate, chlorfenapyr, hexazinone, tebuconazole, resmethrin I i II, fenpropathrin, tebufenpyrad, phenothrin I i II (cis i trans), pyriproxyfen, acrinathrin I i II, fenarimol, trans-permethrin, fluquinconazole, cyfluthrin I – IV, pyridaben, cypermethrin I i II, flucythrinate I i II, fluvalinate I i II, deltamethrin, tributyl phosphate (TBP), trichloroethyl phosphate (TCEP), tris chloropropyl phosphate (TCPP)</i> 0,02 µg/l Prometrin / Prometryn 0,021 µg/l Folpet / Folpet 0,025 µg/l Mankozebe / Mancozeb Ukupno / In total 147				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
A6	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, površinske i podzemne vode, morska voda (prijelazne i priobalne vode, teritorijalno more) Water for human consumption, natural spring water, natural mineral water, table water, bathing water, surface water and ground water, sea water (transitional and coastal water, territorial sea)	Određivanje određenih alkilfenola i njihovih etoksilata, bisfenola A i pentaklorfenola metodom plinske kromatografije s masenom spektrometrijom i μDrop pripremom <i>Determination of selected alkylphenols, their ethoxylates, bisphenol A and pentachlorophenol using gas chromatographic method with mass spectrometry and μDrop sample preparation</i> 0,0002 μg/L Bisfenol A / <i>Bisphenol A</i> 4-(1,1,3,3-Tetrametilbutil)fenol (4-tert-oktilfenol) / 4-(1,1,3,3-Tetramethylbutyl)phenol (4-tert-Octylphenol) 4-(1,1,3,3-Tetrametilbutil)fenol-monoetoksilat, (4-tert-oktilfenol-monoetoksilat) / 4-(1,1,3,3-Tetramethylbutyl)phenol-monoethoxylate (4-tert-Octylphenol-monoethoxylate) 4-(1,1,3,3-Tetrametilbutil)fenol-dietoksilat (4-tert-oktilfenol-dietoksilat) / 4-(1,1,3,3-Tetramethylbutyl)phenol-dietoxylate (4-tert-Octylphenol-dietoxylate) 0,001 μg/L 4-Nonilfenol (iso-nonilfenol) / 4-Nonylphenol (iso-Nonylphenol) 4-Nonilfenol-monoetoksilat (iso-nonilfenol-monoetoksilat) / 4-Nonylphenol-monoethoxylate (iso-Nonylphenol-monoethoxylate) 4-Nonilfenol-dietoksilat (iso-nonilfenol-dietoksilat) / 4-Nonylphenol-dietoxylate (iso-Nonylphenol-dietoxylate) Pentaklorfenol / <i>pentachlorophenol</i> 0,1 μg/L Fenol / <i>Phenol</i>	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31-33,37/249 II. izdanje / <i>edition</i> (2025-01-30)	25/06/2023	25/06/2023
A7	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i	Određivanje organokositrenih spojeva (tributilkositreni, dibutilkositreni, monobutilkositreni, difenilkositreni, monofenilkositreni, monooktilkositreni, tetrabutilkositreni, dioktilkositreni, tricikloheksilkositreni, trifenilkositreni spojevi) metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) i μDrop pripremom	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31-33,37/204a	10/04/2023	10/04/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	podzemne vode, eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more) <i>Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)</i>	<i>Determination of organotin compounds (tributyltin, dibutyltin, monobutyltin, diphenyltin, monophenyltin, mono-octyltin, tetrabutyltin, dioctyltin, tricyclohexyltin, triphenyltin compounds) using gas chromatographic method with mass spectrometry (GC-MS/MS) and μDrop sample preparation</i> 0,00005 μg/L		I. izdanje / edition (2023-04-10)		
A8	A8-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, površinske i podzemne vode <i>Water for human consumption, natural spring water, natural mineral water, table water, bathing water, surface water and ground water</i>	Određivanje poluhlapljivih organskih spojeva (SVOC) u vodama metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) i μDrop pripremom <i>Determination of semi-volatile organic compounds (SVOCs) in water samples using gas chromatographic method with mass spectrometry (GC-MS/MS) and μDrop sample preparation</i> 0,00005 μg/l Cibutrin (irgarol), heptaklor, heptaklorepoksid, dikofol, kloroneb, alfa-BHC (alfa-HCH), heksaklorbenzen, pentakloroanisol, gama-BHC (gama-HCH ili Lindan), delta-BHC (delta-HCH), pentaklorotioanisol, aldrin, 4,4'-diklorobenzofenon, fenon, izodrin, o,p'-DDE (2,4-DDE), alfa-endosulfan (endosulfan I), klorfenon, p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), 4,4'-metoksiklor olefin, endosulfan sulfat, p,p'-DDT (4,4-DDT), 2,4'-metoksiklor (metoksiklor), pentaklorbenzen, etilan (perten), bifenoks, akonifen, klorpirifos metil, piridafention, EPN, pirimifos etil,	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31-33,37/248 II. izdanje / edition (2024-10-28)	10/04/2023	09/04/2025

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		difenilamin, pendimetalin, nitrofen, nitralin, terbufos, fenklorfos, fention, protiofos, sulprofos, fonofos, triazofos, etridiazol, tri-alat, dimetaklor, acetoklor, transflutrin, triadimefon, ciprodinil, penkonazol, oksadiazon, tebukonazol, fenpropatrin, piripoksifen, akrinatri I i II, flukvinkonazol, ciflutrin I-IV (suma izomera), cipermetrin I-IV (suma izomera), alaklor, butilat, diklobenil, EPTC, etalfluralin, fenitrothion, fenvalerat I, fluridon, flusilazol, jodofenfos, malation, metazaklor, metakrifos, metoksiklor, napropamid, norflurazon, pentakloroanilin, pentaklorobenzonitril, profenofos, piraklofos, pirazofos, tolklofos-metil / <i>cibutrin (irgarol), heptachlor and heptachlor epoxide, dicofol, chloroneb, alfa-BHC (alfa-HCH), hexachlorbenzene, pentachloroanisole, gamma-BHC (gamma-HCH or Lindane), delta-BHC (delta-HCH), pentachlorothioanisole, aldrin, 4,4'-dichlorobenzophenone, fenson, isodrin, o,p'-DDE (2,4-DDE), alpha-endosulfan (endosulfan I), chlofenson (Ovex), p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), 4,4'-metoxychlor olefin, endosulfan sulfate, p,p'-DDT (4,4-DDT), 2,4'-metoxychlor (metoxychlor), pentachlorobenzene, ethylan (perthane), bifenox, aclonifen, chlorpyrifos-methyl, pyridaphenthion, EPN, pirimiphos ethyl, iphenylamine, pendimethalin, nitrofen, nitralin, terbufos, fenchlorphos, fenthion, prothiofos, sulprofos, fonofos, triazophos, etridiazole, triallate, dimethachlor, acetochlor, transfluthrin, triadimefon, cyprodinil, penconazole, oxadiazon, tebuconazole, fenpropathrin, pyriproxyfen, acrinathrin I i II, fluquinconazole, cyfluthrin I-IV (sum of isomers), cypermethrin I-IV (sum of isomers), alachlor, butylate, dichlobenil, EPTC, ethalfluralin, fenitrothion, fenvalerate I, fluridone, flusilazole, iodofenphos, malathion, metazachlor, methacrifos, methoxychlor, napropamide, norflurazon, pentachloroaniline, pentachlorobenzonitril, profenofos, pyraclofos, pyrazophos, tolclofos-methyl</i> 0,0001 µg/l Diklorvos, beta-BHC (beta-HCH), klorbenzid, endosulfan eter, cis-klordan, trans-klordan, cis-nonaklor, trans-nonaklor, tetradifon, mireks, pirimifos metil, klorpirifos etil (klorpirifos), fosalon, azinfos etil, isazofos, bromofos etil, klortiofos, karbofenotion, leptofos, forat, paration, piperonil butoksid, propaklor, di-alat I i II (cis i trans), klomazon, terbutilazin, propizamid, pirimetanil, teflutrin, metolaklor, linuron, difenamid, procimidon, paklobutrazol, flutolanil, fludioksonil, miklobutanil, klorfenapir, heksazinon,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		tebufenpirad, fenotrin I i II (cis i trans), fenarimol, piridaben, flucitrinat I i II, fluvalinat I i II, deltametrin, prometrin, 2-fenilfenol, 2,3,5,6-tetrakloroanilin, 3,4-dikloroanilin, antrakvinon, benfluralin, bifentrin, bifenil, bromfenvinfos-metil, bromfenvinfos, bromofos-metil, bromopropilat, karfentrazon-etil, klorfenvinfos, klorobenzilat, klorprofam, klortal-dimetil, cis-permetrin, kumafos, diazinon, dikloran, disulfoton, edifenfos, etoprofos, etofenproks, etrimfos, fensulfotion, fipronil, fluazifop-P-butil, flukloralin, iprodion, izopralin (izopropalin), lambda-ihalotrin, lenacil, metalaksil, metidation, molinat, oksiklordan, oksifluorfen, paration-metil, fosmet, prodiamin, propargit, propazin, propizoklor, sulfotep, teknazen, tetraklorvinfos, trifluralin, vinklozolin / <i>dichlorvos, beta-BHC (beta-HCH), chlordane, endosulfan ether, cis-chlordan, trans-chlordan, cis-nonachlor, trans-nonachlor, tetradifon, mirex, pirimiphos-methyl, chlorpyrifos-ethyl (chlorpyrifos), phosalone, azinphos-ethyl, isazofos, bromophos-ethyl, chlorthiophos, carbophenothion, leptophos, phorate, parathion, piperonyl butoxide, propachlor, diallate I i II, (cis i trans), clomazone, terbuthylazine, propyzamide, pyrimethanil, tefluthrin, metolachlor, linuron, diphenamid, procymidone, paclobutrazol, flutolanil, fludioxonil, myclobutanil, chlorfenapyr, hexazinon, tebufenpyrad, phenothrin I i II (cis i trans), fenarimol, pyridaben, flucythrinate I i II, fluvalinate I i II, deltamethrin, prometryn, 2-phenylphenol, 2,3,5,6-tetrachloroaniline, 3,4-dichloroaniline, anthraquinone, benfluralin, bifenthrin, biphenyl, bromfenvinfos-methyl, bromfenvinfos, bromophos-methyl, bromopropylate, carfentrazone-ethyl, chlorfenvinfos, chlorobenzilate, chlorpropham, chlorthal-dimethyl, cis-permethrin, coumaphos, diazinon, dicloran, disulfoton, edifenfos, ethoprophos, etofenprox, etrimfos, fensulfotion, fipronil, fluazifop-P-butyl, fluchloralin, iprodione, isopralin (isopropalin), lambda-cyhalothrin, lenacil, metalaxyl, methidathion, molinate, oxychlordan, oxyfluorfen, parathion-methyl, phosmet, prodiamine, propargite, propazine, propisochlor, sulfotep, tecnazene, tetrachlorvinphos, trifluralin, vinclozolin</i> 0,0002 µg/l Etion, fenamifos / <i>Ethion, fenamiphos</i> 0,001 µg/l Klorotalonil / <i>Chlorothalonil</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		0,00005 µg/l Polibromirani difenil-eteri (PBDE): BDE-028, BDE-047, BDE-099, BDE-100, BDE-153, BDE-154, BDE-183 / <i>polybrominated diphenyl ethers (PBDE): BDE-028, BDE-047, BDE-099, BDE-100, BDE-153, BDE-154, BDE-183</i> 0,00005 µg/l Poliklorirani bifenili (PCB): PCB-28, PCB-31, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153, PCB-180, PCB-194 / <i>polychlorinated biphenyls (PCB): PCB-28, PCB-31, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153, PCB-180, PCB-194</i> 0,00005 µg/l Policiklički aromatski ugljikovodici (PAH): naftalen, acenaftilen, acenaften, fluoren, fenantren, antracen, fluoranten, piren, benzo[a]antracen, krizen, benzo[b]fluoranten, benzo[k]fluoranten, benzo[a]piren, indeno[1,2,3-c,d]piren, dibenzo[a,h]antracen, benzo[g,h,i]perilen / <i>Polycyclic aromatic hydrocarbons (PAH): Naphthalene, Acenaphthylen, Acenaphthene, Fluorene, Phenanthrene, Anthracene, Fluoranthene, Pyrene, Benzo[a]anthracene, Chrysene, Benzo[b]fluoranthene, Benzo[k]fluoranthene, Benzo[a]pyrene, Indeno[1,2,3-cd]pyrene, Dibenzo[a,h]anthracene, Benzo[g,h,i]perylene</i> 0,005 µg/l Ftalat esteri: dimetil ftalat (DMP), dibutil ftalat (DBP), butil benzil ftalat (BBzP), di (2-etilheksil) ftalat (DEHP), dietil ftalat (DEP), di-n-oktil ftalat (DNOP) / <i>Phthalate esters: dimethyl phthalate (DMP), dibutyl phthalate (DBP), butyl benzyl phthalate (BBzP), bis (2-ethylhexyl) phthalate (DEHP), diethyl phthalate (DEP), di-n-octyl phthalate (DNOP)</i> 0,005 µg/l Fosfat esteri: tributilfosfat (TBP), trikloroetilfosfat (TCEP), trikloropropilfosfat (TCPP) / <i>Phosphate esters: tributyl phosphate (TBP), trichloroethyl phosphate (TCEP), tris chloropropyl phosphate (TCPP)</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Ukupno / In total 226 0,05 µg/l Kratkolančani poliklorirani alkani C10-C13 / <i>short-chain polychlorinated alkanes C10-C13</i>				
	A8-2 Morska voda (prijelazne i priobalne vode, teritorijalno more) <i>Sea water</i> <i>(transitional and</i> <i>coastal water;</i> <i>territorial sea)</i>	Određivanje poluhlapljivih organskih spojeva (SVOC) u vodama metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) i µDrop pripremom <i>Determination of semi-volatile organic compounds (SVOCs) in water samples using gas chromatographic method with mass spectrometry (GC-MS/MS) and µDrop sample preparation</i> 0,00005 µg/l Cibutrin (irgarol), heptaklor, heptaklorepoksid, dikofol, kloroneb, alfa-BHC (alfa-HCH), heksaklorbenzen, pentakloroanisol, gama-BHC (gama-HCH ili Lindan), delta-BHC (delta-HCH), pentaklorotioanisol, aldrin, 4,4'-diklorobenzofenon, fenon, izodrin, o,p'-DDE (2,4-DDE), alfa-endosulfan (endosulfan I), klorfenon, p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), 4,4'-metoksiklor olefin, endosulfan sulfat, p,p'-DDT (4,4-DDT), 2,4'-metoksiklor (metoksiklor), pentaklorbenzen, etilan (perten), bifenoks, aktonifen, klorpirifos metil, piridafention, EPN, pirimifos etil, difenilamin, pendimetalin, nitrofen, nitratin, terbufos, fenklorfos, fention, protiofos, sulprofos, fonofos, triazofos, etridiazol, tri-alat, dimetaklor, acetoklor, transflutrin, triadimefon, ciprodinil, penkonazol, oksadiazon, tebukonazol, fenpropatrin, piriproksifen, akrinatri I i II, flukvinkonazol, ciflutrin I-IV, cipermetrin I i II, alaklor, butilat, diklobenil, EPTC, etalfluralin, fenitroton, fenvalerat I, fluridon, flusilazol, jodofenfos, malation, metazaklor, metakrifos, metoksiklor, napropamid, norflurazon, pentakloroanilin, pentaklorobenzonitril, profenofos, piraklofos, pirazofos, tolklofos-metil, diklorvos, beta-BHC (beta-HCH), klorbenzid, endosulfan eter, cis-klordan, trans-klordan, cis-nonaklor, trans-nonaklor, tetradifon, mireks, pirimifos metil, klorpirifos etil (klorpirifos), fosalon, azinfos etil, isazofos, bromofos etil, klortiofos, karbofenotion, leptofos, forat, paration, piperonil butoksid, propaklor, di-alat I i II (cis i trans), klomazon, terbutilazin, propizamid,	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31- 33,37/248 II. izdanje / <i>edition</i> (2024-10-28)	10/04/2023	28/10/2024

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		pirimetanil, teflutrin, metolaklor, linuron, difenamid, procimidon, paklobutrazol, flutolanil, fludioksonil, miklobutanil, klorfenapir, heksazinon, tebufenpirad, fenotrin I i II (cis i trans), fenarimol, piridaben, flucitrinat I i II, fluvalinat I i II, deltametrin, prometrin, 2-fenilfenol, 2,3,5,6-tetrakloroanilin, 3,4-dikloroanilin, antrakvinon, benfluralin, bifentrin, bifenil, bromfenvinfos-metil, bromfenvinfos, bromofos-metil, bromopropilat, karfentrazon-etil, klorfenvinfos, klorobenzilat, klorprofam, klortal-dimetil, cis-permetrin, kumafos, diazinon, dikloran, disulfoton, edifenfos, etoprofos, etofenproks, etrimfos, fensulfoton, fipronil, fluazifop-P-butyl, flukloralin, iprodion, izopralin (izopropalin), lambda-cihalotrin, lenacil, metalaksil, metidation, molinat, oksiklordan, oksifluorfen, paration-metil, fosmet, prodiamin, propargit, propazin, propizoklor, sulfotep, teknazen, tetraklorvinfos, trifluralin, vinklozolin, etion, fenamifos / <i>cibutrin (irgarol), heptachlor and heptachlor epoxide, dicofof, chloroneb, alfa-BHC (alfa-HCH), hexachlorbenzene, pentachloroanisole, gamma-BHC (gamma-HCH or Lindane), delta-BHC (delta-HCH), pentachlorothioanisole, aldrin, 4,4'-dichlorobenzophenone, fenson, isodrin, o,p'-DDE (2,4-DDE), alpha-endosulfan (endosulfan I), chlofenson (Ovex), p,p'-DDE (4,4-DDE), dieldrin, o,p'-DDD (2,4-DDD), endrin, beta-endosulfan (endosulfan II), p,p'-DDD (4,4-DDD), o,p'-DDT (2,4-DDT), 4,4'-metoxychlor olefin, endosulfan sulfate, p,p'-DDT (4,4-DDT), 2,4'-metoxychlor (metoxychlor), pentachlorobenzene, ethylan (perthane), bifenox, aclonifen, chlorpyrifos-methyl, pyridaphenthion, EPN, pirimiphos ethyl, iphenylamine, pendimethalin, nitrofen, nitalin, terbufos, fenclorophos, fenthion, prothiofos, sulprofos, fonofos, triazophos, etridiazole, triallate, dimethachlor, acetochlor, transfluthrin, triadimefon, cyprodinil, penconazole, oxadiazon, tebuconazole, fenpropathrin, pyriproxyfen, acrinathrin I i II, fluquinconazole, cyfluthrin I-IV, cypermethrin I i II, alachlor, butylate, dichlobenil, EPTC, ethalfluralin, fenitrothion, fenvalerate I, fluridone, flusilazole, iodofenphos, malathion, metazachlor, methacrifos, methoxychlor, napropamide, norflurazon, pentachloroaniline, pentachlorobenzonitril, profenofos, pyraclofos, pyrazophos, tolclorfos-methyl, dichlorvos, beta-BHC (beta-HCH), chlorbenside, endosulfan ether, cis-chlordan, trans-chlordan, cis-nonachlor, trans-nonachlor, tetradifon, mirex, pirimiphos-methyl, chlorpyrifos-ethyl (chlorpyrifos), phosalone, azinphos-ethyl, isazofos, bromophos-ethyl, chlorthiophos, carbophenothion, leptophos, phorate, parathion, piperonyl butoxide, propachlor, diallate I i II, (cis i trans), clomazone, terbutylazine,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		propyzamide, pyrimethanil, tefluthrin, metolachlor, linuron, diphenamid, procymidone, paclobutrazol, flutolanil, fludioxonil, myclobutanil, chlorfenapyr, hexazinon, tebufenpyrad, phenothrin I i II (cis i trans), fenarimol, pyridaben, flucythrinate I i II, fluvalinate I i II, deltamethrin, prometryn, 2-phenylphenol, 2,3,5,6-tetrachloroaniline, 3,4-dichloroaniline, anthraquinone, benfluralin, bifenthrin, biphenyl, bromfenvinfos-methyl, bromfenvinfos, bromophos-methyl, bromopropylate, carfentrazone-ethyl, chlorfenvinphos, chlorobenzilate, chlorpropham, chlorthal-dimethyl, cis-permethrin, coumaphos, diazinon, dicloran, disulfoton, edifenphos, ethoprophos, etofenprox, etrimfos, fensulfothion, fipronil, fluazifop-P-butyl, fluchloralin, iprodione, isopralin (isopropalin), lambda-cyhalothrin, lenacil, metalaxyl, methidathion, molinate, oxychlordan, oxyfluorfen, parathion-methyl, phosmet, prodiamine, propargite, propazine, propisochlor, sulfotep, tecnazene, tetrachlorvinphos, trifluralin, vinclozolin, ethion, fenamiphos 0,00005 µg/l Polibromirani difenil-eteri (PBDE): BDE-028, BDE-047, BDE-099, BDE-100, BDE-153, BDE-154, BDE-183 / <i>polybrominated diphenyl ethers (PBDE): BDE-028, BDE-047, BDE-099, BDE-100, BDE-153, BDE-154, BDE-183</i> 0,00005 µg/l Poliklorirani bifenili (PCB): PCB-28, PCB-31, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153, PCB-180, PCB-194 / <i>polychlorinated biphenyls (PCB): PCB-28, PCB-31, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153, PCB-180, PCB-194</i> 0,00005 µg/l Policiklički aromatski ugljikovodici (PAH): naftalen, acenaftilen, acenaften, fluoren, fenantren, antracen, fluoranten, piren, benzo[a]antracen, krizen, benzo[b]fluoranten, benzo[k]fluoranten, benzo[a]piren, indeno[1,2,3-c,d]piren, dibenzo[a,h]antracen, benzo[g,h,i]perilen / <i>Polycyclic aromatic hydrocarbons (PAH): Naphthalene, Acenaphthylen, Acenaphthene, Fluorene, Phenanthrene, Anthracene, Fluoranthene, Pyrene, Benzo[a]anthracene, Chrysene, Benzo[b]fluoranthene, Benzo[k]fluoranthene, Benzo[a]pyrene, Indeno[1,2,3-cd]pyrene, Dibenzo[a,h]anthracene, Benzo[g,h,i]perylene</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		0,005 µg/l Ftalat esteri: dimetil ftalat (DMP), dibutil ftalat (DBP), butil benzil ftalat (BBzP), di (2-etilheksil) ftalat (DEHP), dietil ftalat (DEP), di-n-oktil ftalat (DNOP) / <i>Phthalate esters: dimethyl phthalate (DMP), dibutyl phthalate (DBP), butyl benzyl phthalate (BBzP), bis (2-ethylhexyl) phthalate (DEHP), diethyl phthalate (DEP), di-n-octyl phthalate (DNOP)</i> 0,005 µg/l Fosfat esteri: tributilfosfat (TBP), trikloroetilfosfat (TCEP), trikloropropilfosfat (TCPP) / <i>Phosphate esters: tributyl phosphate (TBP), trichloroethyl phosphate (TCEP), tris chloropropyl phosphate (TCPP)</i> Ukupno / In total 225 0,05 µg/l Kratkolančani poliklorirani alkani C10-C13 / <i>short-chain polychlorinated alkanes C10-C13</i>				
	A8-3 Otpadna voda <i>Waste water</i>	Određivanje poluhlapljivih organskih spojeva (SVOC) u vodama metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) i µDrop pripremom <i>Determination of semi-volatile organic compounds (SVOCs) in water samples using gas chromatographic method with mass spectrometry (GC-MS/MS) and µDrop sample preparation</i> 0,00001 mg/l 2,4-DDD (o,p'-DDD), 2,4-DDE (o,p'-DDE), 2,4-DDT (o,p'-DDT), 2,4-metoksiklor, 2,4-DDD (o,p'-DDD), 4,4-DDE (p,p'-DDE), 4,4-DDT (p,p'-DDT), 4,4-metoksiklor olefin, acetoklor, aklonifen, aldrin, alfa-HCH (alfa-BHC), atrazin, beta-HCH (beta-BHC), bifenoks, bromfenvinfos-metil, bromofos etil, bromofos metil, kloroneb, klorpirifos, cis-klordan, cis-permetrin, cipermetrin, delta-HCH (delta-BHC), diklobenil, diklorvos, dikofol, dieldrin, endosulfan I, endosulfan II, endosulfan sulfat, endrin, gama-HCH (gama-BHC, lindan), heptaklor, heptaklor-endo-epoksid, heptaklor-egzo-	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-31-33,37/248 II. izdanje / <i>edition</i> (2024-10-28)	10/04/2023	28/10/2024

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		epoksid, heksaklorbenzen, irgarol (cibutrin), izodrin, metazaklor, metoksiklor, metolaklor, paration etil, paration metil, pendimetalin, pentaklorbenzen, pirimifos etil, pirimifos metil, prometrin, propazin, TBP (tributilfosfat),), trikloroetilfosfat (TCEP), trikloropropilfosfat (TCPP), TEPP (tetraetilpirofosfat), tolklofos metil, trans-klordan, trans-permetrin / 2,4-DDD (o,p'-DDD), 2,4-DDE (o,p'-DDE), 2,4-DDT (o,p'-DDT), 2,4-methoxychlor, 2,4-DDD (o,p'-DDD), 4,4-DDE (p,p'-DDE), 4,4-DDT (p,p'-DDT), 4,4-methoxychlor olefin, acetochlor, aclonifen, aldrin, alpha-HCH (alpha-BHC), atrazin, beta-HCH (beta-BHC), bifenox, bromfenvinfos-methyl, bromophos ethyl, bromophos methyl, chloroneb, chlorpyrifos, cis-chlordan, cis-permethrin, cypermethrin, delta-HCH (delta-BHC), dichobenil, dichlorvos, dicofol, dieldrin, endosulfan I, endosulfan II, endosulfan sulfate, endrin, gamma-HCH (gamma-BHC, lindane), heptachlor, heptachlor-endo-epoxide, heptachlor-egzo-epoxide, hexachlorobenzene, irgarol (cibuthrin), isodrin, metazachlor, methoxychlor, metolachlor, parathion ethyl, parathion methyl, pendimethalin, pentachlorobenzene, pirimiphos ethyl, pirimiphos methyl, prometryn, propazine, tributyl phosphate (TBP), trichloroethyl phosphate (TCEP), tris chloropropyl phosphate (TCPP), TEPP (Tetraethyl pyrophosphate), tolclorfos methyl, trans-chlordan, trans-permethrin 0,00001 mg/l Polibromirani difenil-eteri (PBDE): BDE-028, BDE-047, BDE-099, BDE-100, BDE-153, BDE-154, BDE-183 / polybrominated diphenyl ethers (PBDE): BDE-028, BDE-047, BDE-099, BDE-100, BDE-153, BDE-154, BDE-183 0,00001 mg/l Poliklorirani bifenili (PCB): PCB-28, PCB-31, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153, PCB-180, PCB-194 / polychlorinated biphenyls (PCB): PCB-28, PCB-31, PCB-52, PCB-101, PCB-118, PCB-138, PCB-153, PCB-180, PCB-194 0,00001 mg/l Policiklički aromatski ugljikovodici (PAH): naftalen, acenaftilen, acenaften, fluoren, fenantren, antracen, fluoranten, piren, benzo[a]antracen, krizen,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		benzo[b]fluoranten, benzo[k]fluoranten, benzo[a]piren, indeno[1,2,3-c,d]piren, dibenzo[a,h]antracen, benzo[g,h,i]perilen / <i>Polycyclic aromatic hydrocarbons (PAH): Naphthalene, Acenaphthylen, Acenaphthene, Fluorene, Phenanthrene, Anthracene, Fluoranthene, Pyrene, Benzo[a]anthracene, Chrysene, Benzo[b]fluoranthene, Benzo[k]fluoranthene, Benzo[a]pyrene, Indeno[1,2,3-cd]pyrene, Dibenzo[a,h]anthracene, Benzo[g,h,i]perylene</i> 0,00005 mg/l Ftalat esteri: dimetil ftalat (DMP), dibutil ftalat (DBP), butil benzil ftalat (BBP), di (2-etilheksil) ftalat (DEHP), dietil ftalat (DEP), di-n-oktil ftalat (DNOP) / <i>Phthalate esters: dimethyl phthalate (DMP), dibutyl phthalate (DBP), butyl benzyl phthalate (BBzP), bis (2-ethylhexyl) phthalate (DEHP), diethyl phthalate (DEP), di-n-octil phthalate (DNOP)</i> Ukupno / In total 95 0,00005 mg/l Kratkolančani poliklorirani alkani C10-C13 / <i>short-chain polychlorinated alkanes C10-C13</i>				
A9	A9-1 Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode, eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more)	Određivanje odabranih sredstava za zaštitu bilja u vodi metodom tekućinske kromatografije s masenom spektrometrijom <i>Determination of selected plant treatment agents in water – liquid chromatographic method with mass spectrometry</i> 0,002 µg/l Bromacil, trifluralin, kinoksifen, cianazin, alaklor, simazin, diuron, heksazinon, dimetenamid, metazaklor, metolaklor, sebumeton, N,N-dietil-m-toluamid, sebutilazin, mevinfos, atrazin desetil, diazinon, terbutilazin, terbutrin, malation, klorfenvinfos, diklofenak, propazin, klorpirifos, atrazin, izoproturon, dimetoat, klorbromuron, metamitron, metoksuron, terbumeton, metalaksil, klorotoluron, fluometuron I, fluometuron II, imidakloprid, linuron, monuron, neburon, triadimefon, buturon, desmetil klorotoluron, metribuzin, azoksistrobin, ciprokonazol, propikonazol, acetamiprid /	LC-MS/MS	Vlastita metoda In-house method SOP-LEK-31-33,37/183 VI. izdanje / edition (2022-06-01)	14/12/2017	01/06/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)	Bromacil, trifluralin, quinoxifen, cyanazine, alachlor, simazine, diuron, hexazinon, dimethenamide, metazachlor, metolachlor, sebumeton, N,N-diethyl-m-toluamide, sebutylazine, mevinphos, atrazine desethyl, diazinon, terbutylazine, terbutryn, malathion, chlorfenvinphos, dichlofenac, propazine, chlorpyrifos (Dursban), atrazine, isoproturon, dimethoate, chlorbromuron, metamidron, metoxuron, terbuteton, metalaxyl, chlorotoluron, fluometuron I, fluometuron II, imidacloprid, linuron, monuron, neburon, triadimefon, buturon, chlorotoluron-desmethyl, metribuzin, azoxystrobin, cyproconazole, propiconazole, acetamiprid 0,006 µg/l Metobromuron / metobromuron 0,008 µg/l Atrazin-desizopropil /Atrazine-desisopropyl 0,012 µg/l Monolinuron / monolinuron 0,02 µg/l Piretrin, ometoat, flzasulfuron / pyrethrin, omethoat, flzasulfuron Ukupno / In total 52				
	A9-2 Voda za ljudsku potrošnju, prirodne mineralne, prirodne izvorske, stolne vode, podzemne, površinske vode, procjedne i otpadne vode i eluati otpada Water for human consumption, natural spring water, natural	Određivanje kiselih pesticida u vodi metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) Determination of acidic pesticides in water – liquid chromatographic method with mass spectrometry (LC-MS/MS) 0,002 µg/L 2,4-DP, bentazon, dikamba, MCPB, 2,4-DB, MCPA, fenoprop (Silveks), 2,4-D, ioksinil, MCPP (mekoprop), dikloroprop, bromoksinil, 2,6-diklorobenzamid, 2,4,5-T /	LC-MS/MS	Vlastita metoda In-house method SOP-LEK-31-33,37/182 II. izdanje / edition (2020-02-03) modificirana /modified DIN 38407-35:2010	20/04/2016	10/07/2024

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	mineral water; table water; ground water; surface water; leachate water; waste water and waste eluate	2,4-DP, bentazon, dicamba, MCPB, 2,4-DB, MCPA, fenoprop (Silvex), 2,4-D, ioxynil, MCPP (mecoprop), dichloroprop, bromoxynil, 2,6-dichlorobenzamid, 2,4,5-T				
	A9-3 Površinske vode Surface water	Određivanje glifozata u vodi metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) Determination of Glyphosate in water – liquid chromatographic method with mass spectrometry (LC-MS/MS) 0,1 µg/L	LC-MS/MS	Vlastita metoda In-house method SOP-LEK-31,37/243 I. izdanje / edition (2022-06-03)	03/06/2022	25/7/2023
	A9-4 Voda za ljudsku potrošnju, prirodne mineralne, prirodne izvorske, stolne vode, podzemne i površinske vode Water for human consumption, natural spring water; natural mineral water; table water; ground water and surface water	Određivanje metabolita metolaklora i acetoklora u vodi metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) Determination of metabolites of Metolachlor and Acetochlor in water – liquid chromatographic method with mass spectrometry (LC-MS/MS) 0,020 µg/L Metolaklor OXA, metolaklor ESA, acetoklor OXA, acetoklor ESA / Metolachlor OXA, Metolachlor ESA, Acetochlor OXA, Acetochlor ESA	LC-MS/MS	Vlastita metoda In-house method SOP-LEK-31,37/265 I. izdanje / edition (2024-10-16)	16/10/2024	17/10/2024
A10	Voda za ljudsku potrošnju, prirodne izvorske, prirodne mineralne, stolne vode, vode za kupanje, otpadne vode, površinske, procjedne i podzemne vode,	Određivanje PFAS-a u vodi metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) Determination of PFAS in water by liquid chromatography-mass spectrometry (LC-MS/MS) 0,000039 µg/L	LC-MS/MS	Vlastita metoda In-house method SOP-LEK-31-33 i 37/214	24/04/2019	17/02/2025

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change			
	eluat otpada, morska voda (prijelazne i priobalne vode, teritorijalno more) Water for human consumption, natural spring water, natural mineral water, table water, bathing water, waste water, surface water, leachate water and ground water, waste eluat, sea water (transitional and coastal water, territorial sea)	PFOS 0,030 µg/L <i>PFBS, PFHxA, PFOA, PFNA, PFDA, PFNS, PFBA, 4-2 FTS, 6-2 FTS, 8-2 FTS, PFPeS, PFHpA, PFHxS, PFHpS</i>		III. izdanje / edition (2024- 11-12)					
B1	Eluat otpada Waste leachate	Određivanje odabranih elemenata <i>Determination of selected elements</i>		ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	03/01/2012	17/08/2026		
		mg/kg suhe tvari							
		Element/ Element	Granica kvantifikacije / Limit of quantification					Element/ Element	Granica kvantifikacije / Limit of quantification
		Ag	0,01					P	0,28
		Al	0,10					Pb	0,22
		As	0,14					Sb	0,10
		B	0,06					Se	0,79
		Ba	0,06					Si	0,46
		Be	0,01					Sn	0,31
		Cd	0,04					Sr	0,29
		Co	0,04					Te	0,19
		Cr	0,02					Ti	0,03
		Cu	0,04					Tl	0,31

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Fe	0,08	U	0,35				
		Mn	0,01	V	0,15				
		Mo	0,07	Zn	0,06				
		Ni	0,06	P	0,28				
B2	B2–1 Otpad Waste	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN ISO 22036:2024 (ISO 22036:2024; EN ISO 22036:2024)	13/09/2013	10/03/2026
		mg/kg							
		Element/ Element	Granica kvantifikacije / Limit of quantification	Element/ Element	Granica kvantifikacije / Limit of quantification				
		As	1,45	Ni	0,73				
		Cd	0,21	Pb	2,40				
		Co	0,39	Sb	3,15				
		Cr	0,04	Tl	2,23				
		Cu	0,36	V	1,54				
		Mn	0,18	Zn	0,73				
	Mo	0,85							
	B2–2 Otpad Waste	Određivanje ukupnih sadržaja elemenata – raščinjavanje alkalnom fuzijom Determination of total contents of elements - dissolution by alkaline fusion				ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	06/11/2019	15/01/2026
		mg/kg							
		Element/ Element	Granica kvantifikacije / Limit of quantification	Element/ Element	Granica kvantifikacije / Limit of quantification				
		Al	12,50	Mn	1,25				
		Ca	18,75	Na	105,00				
		Cu	5,00	Pb	43,75				
		Fe	10,00	Si	57,50				
K		76,25	Sr	36,25					
Mg		93,75	Zn	7,50					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		oksidi / oxides - %							
		Oksid/ Oxide	Granica kvantifikacije / Limit of quantification	Oksid/ Oxide	Granica kvantifikacije/ Limit of quantification				
		Al ₂ O ₃	0,0024	MnO ₂	0,0002				
		CaO	0,0026	Na ₂ O	0,0142				
		Fe ₂ O ₃	0,0014	PbO ₂	0,0051				
		K ₂ O	0,0092	SiO ₂	0,0123				
		MgO	0,0155	SrO	0,0043				
		B3	B3–1 Eluat otpada Waste leachate	Određivanje žive Determination of mercury 0,00007 mg/kg suhe tvari					
B3–2 Otpad Waste	Određivanje žive Determination of mercury 0,001 mg/kg				AAS (FIMS)	HRN EN ISO 12846:2012 (ISO 12846:2012; EN ISO 12846:2012)	01/06/2015	08/10/2024	
B4	Otpad Waste	Određivanje BTEXa i heksaklorbutadiena metodom analize para iznad otopine (Headspace) plinskom kromatografijom sa masenom spektrometrijom (GC-MS/MS) Determination of BTEX and hexachlorobutadiene method Headspace gas chromatography mass spectrometry (GC-MS/MS)				GC-MS/MS	Vlastita metoda In-house method SOP-LEK- 38,39/142 VI. izdanje / edition (2025-03-14)	16/6/2014	14/3/2025
		mg/kg							
		Spoj /Compound		Granica određivanja/ Limit of detection					
		benzen / benzene		0,0002					
		toluen / toluene		0,0001					
		ksileni / xylene		0,0002					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification		Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change																				
		etilbenzen / ethylbenzene	0,0001																								
		heksaklorbutadien / hexachlorobutadiene	0,0002																								
B5	Otpad / Waste	Određivanje pesticida metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Determination of pesticide - gas chromatography with mass spectrometry (GC-MS/MS)</i> Otpad / <i>waste</i> 0,005 mg/kg Alfa-BHC (alfa-HCH), heksaklorbenzen (HCB), beta-BHC (beta-HCH), gama-BHC (gama-HCH; Lindan), heptaklor, aldrin, heptaklor epoksid (izomer B), trans-klordan, endosulfan I (endosulfan alfa), cis-klordane, dieldrin, 4,4-DDE (p,p-DDE), endrin, endosulfan II (endosulfan beta), 4,4-DDD (p,p-DDD), 2,4-DDT (o,p-DDT), endosulfan sulfat, 4,4-DDT (p,p-DDT) / <i>Alpha-BHC (alpha-HCH), Hexachlorobenzene (HCB), beta-BHC (beta-HCH), gamma-BHC (gamma-HCH, Lindane), Heptachlor, Aldrin, Heptachlor epoxide (isomer B), trans-chlordane, Endosulfan I (Endosulfan alpha), cis-Chlordane, Dieldrin, 4,4-DDE (p,p-DDE), Endrin, Endosulfan II (Endosulfan beta), 4,4-DDD (p,p-DDD), 2,4-DDT (o,p-DDT), Endosulfan sulfate, 4,4-DDT (p,p-DDT)</i> Ukupno / In total 18		GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LEK-38,39/03 III. izdanje / <i>edition</i> (2022-04-13)	15/01/2013	13/08/2024																				
C1	Mulj, sediment <i>Sludge, sediment</i>	Određivanje odabranih elemenata <i>Determination of selected elements</i> mg/kg <table><tr><th>Element/ Element</th><th>Granica kvantifikacija / Limit of quantification</th><th>Element/ Element</th><th>Granica kvantifikacija / Limit of quantification</th></tr><tr><td>Al</td><td>0,83</td><td>Ni</td><td>0,44</td></tr><tr><td>As</td><td>0,83</td><td>P</td><td>8,33</td></tr><tr><td>Cd</td><td>0,16</td><td>Pb</td><td>2,00</td></tr><tr><td>Cr</td><td>0,17</td><td>Zn</td><td>0,47</td></tr></table>		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification	Al	0,83	Ni	0,44	As	0,83	P	8,33	Cd	0,16	Pb	2,00	Cr	0,17	Zn	0,47	ICP - OES	HRN EN ISO 22036:2024 (ISO 22036:2024; EN ISO 22036:2024)	13/09/2013	10/02/2025
Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification																								
Al	0,83	Ni	0,44																								
As	0,83	P	8,33																								
Cd	0,16	Pb	2,00																								
Cr	0,17	Zn	0,47																								

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification				Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Cu	0,87						
C2	Tlo Soil	Određivanje odabranih elemenata Determination of selected elements				ICP - OES	HRN EN ISO 22036:2024 (ISO 22036:2024; EN ISO 22036:2024)	23/07/2014	05/02/2025
		mg/kg							
		Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification				
		Al	0,40	Ni	0,44				
		Cd	0,16	Pb	2,00				
		Cr	0,17	Zn	0,47				
		Cu	0,87	As	2,62				
		Fe	0,93	Mo	1,97				
		Mn	0,03	Co	0,74				
C3	Mulj, sediment, biota, tlo Sludge, sediment, biota, soil	Određivanje žive Determination of mercury 0,001 mg/kg – mulj, sediment, biota, tlo / sludge, sediment, biota, soil				AAS (FIMS)	HRN EN ISO 12846:2012 (ISO 12846:2012; EN ISO 12846:2012)	01/06/2015	08/10/2024
C4	C4-1 Tlo, mulj, sediment Soil, sludge, sediment	Određivanje BTEXa i heksaklorbutadiena metodom analize para iznad otopine (Headspace) plinskom kromatografijom sa masenom spektrometrijom (GC-MS/MS) Determination of BTEX and hexachlorobutadiene method Headspace gas chromatography mass spectrometry (GC-MS/MS)				GC-MS/MS	Vlastita metoda In-house method SOP-LEK- 38,39/142 VI. izdanje / edition (2025-03-14)	16/6/2014	14/3/2025
		mg/kg							
		Spoj /Compound		Granica određivanja/ Limit of detection					
		benzen / benzene		0,0002					
		toluen / toluene		0,0001					
		ksileni / xylene		0,0002					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification		Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		etilbenzen / ethylbenzene	0,0001				
		heksaklorbutadien / hexachlorobutadiene	0,0002				
C4	C4-2 Biota Biota	Određivanje heksaklorbutadiena metodom analize para iznad otopine (Headspace) plinskom kromatografijom sa masenom spektrometrijom (GC-MS/MS) Determination of hexachlorobutadiene method Headspace gas chromatography mass spectrometry (GC-MS/MS)		GC-MS/MS	Vlastita metoda In-house method SOP-LEK-38,39/142 VI. izdanje / edition (2025-03-14)	14/3/2025	14/3/2025
		µg/kg m.t.					
		Spoj /Compound	Granica određivanja/ Limit of detection				
		heksaklorbutadien / hexachlorobutadiene	10,0				
C5	Tlo, sediment / Soil, sediment	Određivanje pesticida metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) Determination of pesticide - liquid chromatography with mass spectrometry (LC-MS/MS)		LC-MS/MS	Vlastita metoda In-house method SOP-LKH-38,39/200 V. izdanje / edition (2026-04-29)	10/09/2018	29/4/2026
		0,002 mg/kg Atrazin, simazin, kinoksifen / Atrazine, simazine, quinoxifen					
C6	C6-1 Tlo, sediment / Soil, sediment	Određivanje pesticida metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) Determination of pesticide - gas chromatography with mass spectrometry (GC-MS/MS)		GC-MS/MS	Vlastita metoda In-house method SOP-LEK-38,39/03 III. izdanje / edition (2022-04-13)	15/01/2013	13/08/2024
		Tlo i sediment / soil, sediment 0,005 mg/kg Kloroneb, pentaklorbenzen (PeCB), alfa-BHC (alfa-HCH), heksaklorbenzen (HCB), pentakloroanisol, beta-BHC (beta-HCH), gama-BHC (gama-HCH), delta-BHC (delta-HCH), endosulfan eter, heptaklor, pentaklorotioanisol, aldrin, 4,4-diklorobenzofenon, fenson, izodrin, heptaklor epoksid (izomer B), klorbenzid, trans-klordan, 2,4-DDE (o,p-DDE), endosulfan I (endosulfan alfa), cis-klordane, trans-nonaklor, klorfenson (Ovex), dieldrin, 4,4-DDE (p,p-DDE), 2,4-DDD (o,p-DDD), endrin, endosulfan II (endosulfan beta), 4,4-DDD (p,p-					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		DDD), 2,4-DDT (o,p-DDT), cis-nonaklor, endrin aldehid, 4,4-metoksiklor olefin, endosulfan sulfat, 4,4-DDT (p,p-DDT), 2,4-metoksiklor, endrin keton, tetradifon, mirex, dikofol, cipermetrin / Chloroneb, Pentachlorobenzene (PeCB), alpha-BHC (alpha-HCH), Hexachlorobenzene (HCB), Pentachloroanisole, beta-BHC (beta-HCH), gamma-BHC (gamma-HCH), delta-BHC (delta-HCH), Endosulfan ether, Heptachlor, Pentachlorothioanisole, Aldrin, 4,4-dichlorobenzophenone, Fenson, Isodrin, Heptachlor epoxide (isomer B), Chlorbenside, trans-chlordane, 2,4-DDE (o,p-DDE), Endosulfan I (Endosulfan alpha), cis-Chlordane, trans-Nonachlor, Chlorfenson (Ovex), Dieldrin, 4,4-DDE (p,p-DDE), 2,4-DDD (o,p-DDD), Endrin, Endosulfan II (Endosulfan beta), 4,4-DDD (p,p-DDD), 2,4-DDT (o,p-DDT), cis-Nonachlor, Endrin aldehyde, 4,4-methoxychlor olefin, Endosulfan sulfate, 4,4-DDT (p,p-DDT), 2,4-methoxychlor, Endrin ketone, Tetradifon, Mirex, Dicofol, Cipermetryn Ukupno / In total 41				
	C6-2 Biota / Biota	Određivanje pesticida metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Determination of pesticide - gas chromatography with mass spectrometry (GC-MS/MS)</i> 0,02 µg/kg Heksaklorobenzen (HCB), heptaklor, heptaklor epoksid / Hexachlorobenzene (HCB), heptachlor, heptachlor epoxide 1,0 µg/kg Dikofol / Dicofol Ukupno / In total 4	GC-MS/MS	Vlastita metoda In-house method SOP-LEK-44/03 Izdavanje / I.edition (2024-10-16)	16/10/2024	16/10/2024
C7	Sediment, biota / Sediment, biota	Određivanje PFAS-a metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS)	LC-MS/MS	Vlastita metoda In-house method	05/11/2024	27/01/2025

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change																																				
		Determination of PFAS by liquid chromatography-mass spectrometry (LC-MS/MS) 2,33 µg/kg mokre težine / wet weight PFOS, PFBS, PFHxA, PFOA, PFNA, PFDA, PFNS, PFBA, 4-2 FTS, 6-2 FTS, 8-2 FTS, PFPeS, PFHpA, PFHxS, PFHpS, PFPeA, N-MeFOSAA, FOSA, PFDS, PFUnA, N-EtFOSAA, PFDoA, PFTrDA, PFTeDA		SOP-LEK-38,44/214 II. izdanje / edition (2025-01-27)																																						
C8	Mulj, sediment, tlo / Sludge, sediment, soil	Određivanje polikloriranih bifenila (PCB 28, PCB 52, PCB 101, PCB 118, PCB 138, PCB 141, PCB 153, PCB 180) metodom plinske kromatografije s masenom spektrometrijom (GC-MS) Determination of polychlorinated biphenyls (PCB 28, PCB 52, PCB 101, PCB 118, PCB 141, PCB 138, PCB 153, PCB 180) by using capillary gas chromatography with mass spectrometric detection (GC-MS) 0,010 mg/kg	GC-MS	HRN EN ISO 18475:2025 (EN ISO 18475:2025)	14/3/2014	15/1/2024																																				
D1	Kruta oporabljena goriva Solid recovered fuels	Određivanje odabranih elemenata Determination of selected elements mg/kg <table><thead><tr><th>Element/ Element</th><th>Granica kvantifikacija / Limit of quantification</th><th>Element/ Element</th><th>Granica kvantifikacija / Limit of quantification</th></tr></thead><tbody><tr><td>As</td><td>7,09</td><td>Mo</td><td>0,68</td></tr><tr><td>Ba</td><td>0,60</td><td>Ni</td><td>0,78</td></tr><tr><td>Be</td><td>0,07</td><td>Pb</td><td>4,19</td></tr><tr><td>Cd</td><td>3,92</td><td>Sb</td><td>4,74</td></tr><tr><td>Co</td><td>0,24</td><td>Tl</td><td>3,39</td></tr><tr><td>Cr</td><td>0,26</td><td>V</td><td>1,13</td></tr><tr><td>Cu</td><td>5,34</td><td>Zn</td><td>5,27</td></tr><tr><td>Mn</td><td>0,05</td><td></td><td></td></tr></tbody></table>	Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification	As	7,09	Mo	0,68	Ba	0,60	Ni	0,78	Be	0,07	Pb	4,19	Cd	3,92	Sb	4,74	Co	0,24	Tl	3,39	Cr	0,26	V	1,13	Cu	5,34	Zn	5,27	Mn	0,05			ICP - OES	HRN EN ISO 11885:2010 (ISO 11885:2007, EN ISO 11885:2009)	05/05/2014	21/03/2023
Element/ Element	Granica kvantifikacija / Limit of quantification	Element/ Element	Granica kvantifikacija / Limit of quantification																																							
As	7,09	Mo	0,68																																							
Ba	0,60	Ni	0,78																																							
Be	0,07	Pb	4,19																																							
Cd	3,92	Sb	4,74																																							
Co	0,24	Tl	3,39																																							
Cr	0,26	V	1,13																																							
Cu	5,34	Zn	5,27																																							
Mn	0,05																																									

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
D2	Kruta oporabljena goriva Solid recovered fuels	Određivanje žive Determination of mercury 0,010 mg/kg – kruta oporabljena goriva / solid recovered fuels	AAS (FIMS)	HRN EN ISO 12846:2012 (ISO 12846:2012; EN ISO 12846:2012)	01/06/2015	08/10/2024
E1	E1-1 Hrana i hrana za životinje Food and animal feeding stuffs	Određivanje odabranih elemenata optičkom emisijskom spektrometrijom induktivno vezane plazme (ICP-OES) Determination of selected elements by inductively coupled plasma optical emission spectrometry (ICP-OES) 2,4 mg/kg - Na 1,5 mg/kg - Ca 4,8 mg/kg - P 0,62 mg/kg - Cu	ICP-OES	HRN EN 16943:2017 (EN 6943:2017)	18/04/2016	26/05/2026
	E1-2 Meso i mesni proizvodi	Određivanje ukupne količine fosfora u mesu i mesnim proizvodima Determination of total phosphorus content in meat and meat products by inductively coupled plasma optical emission (ICP-OES) and spectrometric method 26 mg/kg P (ICP-OES) 541 mg/kg P₂O₅ (UV-vis spektrofotometrija / spectrophotometry)	ICP-OES UV-Vis spektrofotometrija / spectrophotometry	HRN ISO 23776:2022 (ISO 23776:2021)	25/10/2023	15/06/2026
E2	Hrana Food	Određivanje žive Determination of mercury 0,02 mg/kg – food / food	AAS (FIMS)	Vlastita metoda In-house method SOP-LKH-01- 28/07a-1 4. izdanje / edition (2026- 04-13)	16/03/2021	13/04/2026

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
E3	Hrana Food	Određivanje odabranih elemenata metodom ICP-MS nakon digestije nitratnom kiselinom <i>Determination of selected elements by ICP-MS method after nitric acid digestion</i> 0,01 mg/kg - As, Pb 0,005 mg/kg - Cd	ICP-MS	HRN EN 15763:2010 (EN 15763:2009)	27/10/2020	16/06/2023
E4	E4-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,005 mg/kg 1-naftalenacetamid, 2,3,5-Trimetakarb, 2,4,5-T, 2,4-D, 2,4-dimetilanilin, 3,4,5-Trimetakarb, 3-Decen-2-on, 3-hidroksikarbofuran, 4-bromofenilurea, 4-klorofenoksiotena kiselina, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Aktrin, Alaklor, Alanikarb, Aldikarb, Aldikarb-sulfoksid, Aldikarb-sulfon (Aldoksikarb), Alidoklor, Ametotradin, Ametrin, Amidosulfuron, Aminokarb, Amisulfuron, Ancimidol, Anilazin, Atraton, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a (Abamektin), Azadirachtin, Azakonazol, Azametifos, Azinfos-etil, Azinfos-metil, Azinprotrifos, Azoksistrobilin, Benalaksil, Bendiokarb, Benfurakarb, Benodanil, Benomil, Bensulfuron-metil, Bensulid, Bentazon, Bentiavalikarb-izopropil, Bentivalikarb, Benziladenin, Benzoilprop-etil, Benzoksimat, Benzovindflupir, Bifenazat-diazen, Bifentrin, Biksafen, Bispiribak-natrij, Bitertanol, Boskalid, Bromacil, Bromofos-etil, Bromoksinil, Bromukonazol, Bupirimat, Buprofezin, Butafenacil, Butaklor, Butilat, Butoksikarboksime, Butralin, Buturon, Centrimonium klorid, Cianazin, Cianfenfos, Cianofos, Ciantraniliprol, Ciazofamid, Ciflufenamid, Ciflumetofen, Cihalofof-butil, Ciklanilid, Cikloat, Cikloksidim, Cikluron, Cimiazol, Cimoksanil, Cinidon-etil, Ciprodinil, Ciprokonazol, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfoksid,	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	17/09/2024

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Demeton-S-metil-sulfon, Desmedifam, Diafentiuron, Diallylate, Diazinon, Dietofenkarb, Difenamid, Difenokonazole, Difenoksuron, Difetialon, Diflubenzuron, Diflufenikan, Diklobutrazol, Diklofenak, Diklofention, Diklofop, Diklofop-metil, Diklorprop, Diklorvos, Dikrotofos, Dimefuron, Dimetaklor, Dimetenamid, Dimetilan, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinoseb, Dinotefuran, Dinoterb, Dioksakarb, Dioksation, Dipropetrin, Disulfoton-sulfoksid, Disulfoton-sulfon, Diuron (DCMU), DMST, DNOC, Dodemorf, Dodin, Edifenfos, Enamektin-benzoat, Endosulfan sulfat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etiofenkarb-sulfon, Etion, Etrimol, Etofenproks, Etofumezat, Etoksakonazol, Etoprofos, Etrimfos, Famoksadon, Fampur, Fenamidon, Fenamifos, Fenamifos-sulfoksid, Fenamifos-sulfon, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenmedifam, Fenobukarb, Fenoksaprop-etil, Fenoksikarb, Fenoprop, Fenotiokarb, Fenotrin, Fenpikoksamid, Fenpirazamin, Fenpiroksimat, Fenpropatrin, Fenpropidin, Fenpropimorf, Fensulfotion, Fensulfotion-sulfon, Fention, Fention-okson, Fention-okson-sulfoksid, Fention-okson-sulfon, Fention-sulfoksid, Fention-sulfon, Fentoat, Fenuron, Fipronil, Fipronil-sulfid, Fipronil-sulfon, Flamprop-izopropil, Flamprop-metil, Flamprop-M-izopropil, Flazasulfuron, Flazuron, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazifop-P-butil, Fluazinam, Flubendiamid, Flucipirim, Flucitrinat, Fludioksonil, Flufenacet, Flufenacet OA, Flufenacet sulfonska kiselina-ESA, Flufenacet tioglikolat sulfoksid, Flufenoksuron, Fluksapiroksad, Flukvinkonazol, Flumetralin, Fluoksastrobin, Fluometuron, Fluopikolid, Flupiradifuron, Fluopiram, Fluridon, Flurokloridon, Fluroksipir, Flusilazol, Flutiacet-metil, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Foksim, Fonofos, Forat, Forat-okson, Forat-okson-sulfon, Forat-sulfoksid, Forat-sulfon, Forklorfenuron, Formetanat hidroklorid, Fosalon, Fosfamidon, Fosfolan, Fosmet, Fosmet-okson, Fostiazat, Fuberidazol, Furalaksil, Furatiokarb, Furmecikloks, Haloksifop, Haloksifop-2-etoksietil, Haloksifop-metil, Haloksifop-P-metil, Heksaflumuron, Heksakonazol, Heksazinon, Heksitiazoks, Heptenofos, Hidrametilnon, Imazakvin, Imazalil, Imazamoks, Imazapir, Imazetapir, Imibenkonazol, Imidakloprid, Indoksakarb, Ioksonil, Ipkonazol, Iprobenfos, Iprovalikarb, Ivermektin B1a, Izazofos, Izofenfos-metil, Izofetamid, Izokarbamid, Izokarbofos, Izoksaben, Izoksadifen-etil, Izoksafutol, Izoksafutol-dikenotril, Izoksation, Izometiozin, Izopirazam, Izoprokarb, Izopropalin, Izoprotiolan, Izoproturon, Jodosulfuron-metil,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Kadusafos, Karanjin, Karbaril (NAC), Karbendazim, Karbetamid, Karbofenotion, Karbofuran, Karboksini, Karboksini sulfoksidi, Karbosulfan, Karfentrazon-etil, Kinosulfuron, Kletodim, Klimbazol, Klodinafop (slobodna kiselina), Klodinafop-propargil, Klofentezin, Klokintocet-meksil, Klomazon, Klopuralid, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Kloridazon-desfenil, Klorantraniliprol, Kloroksuron, Klorotoluron, Klorpirifos, Klorpirifos-metil, Klorprofam, Klorosulfuron, Klortiamid, Klotianidin, Krezoksime-metil, Krofumat, Krotoksifos, Kumafos, Kumatetrilil, Kvinalfos, Kvinklorak, Kvinmerak, Kvinoklamin, Kinoksifen, Kvizalofop (slobodna kiselina), Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), MCPB, Mefenacet, Mefenpir-dietil, Mefentriklonazol, Mefosfolan, Mekarbam, Meksakarbat, Mepanipirim, Mepronil, Meptildinokap, Metabenzotiazuron, Metaflumizon, Metakrifos, Metalaksil i Metalaksil-M, Metamidofos, Metamitron, Metazaklor, Metazaklor ESA, Metazaklor OA, Metidation, Metiokarb, Metiokarb-sulfoksidi, Metiokarb-sulfon, Metkonazol, Metobromuron, Metoksifenozi, Metoksuron, Metolaklor, Metolcarb (MTMC), Metomil, Metominostrobin (E, Z), Metoprotrilin, Metrafenon, Metribuzin, Metsulfuron-metil, Mevinfos, Miklobutanil, Milbemektin A3, Milbemektin A4, Molinat, Monokrotofos, Monolinuron, Monuron, N-(2,4-dimetilfenil)-N-metilformamidin, N-(2,4-Dimetilfenil) formamid, N,N-dietil-m-toluamid, Napropamid, N-desmetil-propamokarb, Neburon, Nereistoksin oksalat, Nikosulfuron, Nitenpiram, Norflurazon, Novaluron, Nuairimol, Ofurak, Oksadiargil, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Oksikarboksini, Oksimatrini, Ometoat, Orizalin, Paklobutrazol, Paraokson-etil, Paraokson-metil, Pebulat, Pencikuron, Pendimetalin, Penflufen, Penkonazol, Pentiopirad, Petoksamid, Picaridin (Icaridin), Pikloram, Pikoksistrobin, Pikolinafen, Pimetrozin, Piperofos, Piperonil-butoksidi, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Pirazokson, Piretrini, Piridaben, Piridafention, Piridafol, Piridalil, Piridat, Pirifeonks, Pirimetanil, Pirimifos-metil, Pirimikarb, Pirimikarb-desmetil, Piriofenon, Piriprosifen, Potasan, Praletrin, Pretilaklor, Profam, Profenofos, Profoksidi, Prokloraz, Prokloraz met. BTS 44595, Prokloraz met. BTS 44596, Prokinazid, Promekarb, Prometon, Prometrin, Propaklor, Propakvizafop, Propamokarb, Propamokarb-N-oksidi, Propanil, Propargit, Propazin, Propetamfos,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Propikonazol, Propizamid, Propoksur, Prosulfokarb, Prosulfuron, Protiokonazol-destio, Rabenzazol, Resmetrin, Rimsulfuron, Rotenon, Saflufenacil, Sebutilazin, Sekbumeton, Setoksidim, Siduron, Silafluofen, Siltiofam, Simazin, Simetrin, Spinetoram, Spinosad (suma Spinosin A i Spinsin D), Spirodiklofen, Spiroksamin, Spiromezifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glukozid, Spirotetramat-keto-hidroksi, Spirotetramat-mono-hidroksi, Sulfoksaflor, Sulfosulfuron, Sulfotep, Sulkotrion, Sulprofos, TCMTB, Tebufenozid, Tebufenpirad, Tebukonazol, Tebupirimifos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrion, Temefos, Tepraloksidim, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetrakonazol, Tetraklorvinfos (CVMP), Tetrametrin, Tiabendazol, Tiaklopid, Tiametoksam, Tidiazuron, Tifensulfuron-metil, Tiobenkarb, Tiociklam hidrogen oksalat, Tiodikarb, Tiofanat-metil, Tiofanoks, Tiofanoks-sulfoksid, Tiofanoks-sulfon, Tiometon, Tionazin, Tiram, Tolfenpirad, Tolilfluanid, Tolklofos-metil, Tralometrin, Triadimefon, Triadimenol, Trialat, Triazamat, Triazofos, Tribenuron-metil, Triciklazol, Trietazin, Trifloksistrobin, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Triflusulfuron-metil, Triforin, Triklopir, Triklorfon, Trisulfuron, Tritikonazol, Tritosulfuron, Valifenalat, Vamidotion, Vernolat, Vinklozolin, XMC, Zoksamid / 1-naphthaleneacetamide, 2,3,5-Trimethacarb, 2,4,5-T, 2,4-D, 2,4-dimethylaniline, 3,4,5-Trimethacarb, 3-Decen-2-one, 3-Hydroxycarbofuran, 4-Bromophenylurea, 4-Chlorophenoxyacetic acid, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Acrinathrin, Alachlor, Alanycarb, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Amidosulfuron, Aminocarb, Amisulbrom, Ancymidol, Anilazine, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a (Abamectin), Azaconazole, Azadirachtin, Azamethiphos, Azinphos-ethyl, Azinphos-methyl, Aziprotryne, Azoxystrobin, Benalaxyl, Bendiocarb, Benfuracarb, Benodanil, Benomyl, Bensulfuron-methyl, Bensulide, Bentazone, Bentiavalicarb, Bentiavalicarb-isopropyl, Benzovindflupyr, Benzoximate, Benzoylprop-ethyl, Benzyladenine, Bifenazate-diazene, Bifenthrin, Bispyribac-sodium, Bitertanol, Bixafen, Boscalid, Bromacil, Bromophos-ethyl, Bromoxynil, Bromuconazole, Bupirimate, Buprofezin, Butachlor, Butafenacil, Butoxycarboxim, Butralin, Buturon, Butylate, Cadusafos, Carbaryl (NAC), Carbendazim,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Carbetamide, Carbofuran, Carbophenothion, Carbosulfan, Carboxin, Carboxin sulfoxide, Carfentrazone-ethyl, Centrimonium chloride, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chloridazon-desphenyl, Chlorotoluron, Chloroxuron, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorsulfuron, Chlorthiamid, Cinidon-ethyl, Cinosulfuron, Clethodim, Climbazole, Clodinafop (free acid), Clodinafop-propargyl, Clofentezine, Clomazone, Clopyralid, Cloquintocet-mexyl, Clothianidin, Coumaphos, Coumatetralyl, Crotoxyphos, Crufomate, Cyanazine, Cyanofenphos, Cyanophos, Cyantraniliprole, Cyazofamid, Cyclanilide, Cycloate, Cycloxydim, Cycluron, Cyflufenamid, Cyflumetofen, Cyhalofop-butyl, Cymiazole, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Di-allate, Diazinon, Dichlofenac, Dichlofenthion, Dichlorprop, Dichlorvos, Diclobutrazol, Diclofop, Diclofop-methyl, Dicrotophos, Diethofencarb, Difenconazole, Difenoxuron, Difethialone, Diflubenzuron, Diflufenican, Dimefuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimetilan, Dimoxystrobin, Diniconazole, Dinoseb, Dinotefuran, Dinoterb, Dioxacarb, Dioxathion, Diphenamid, Dipropetryn, Disulfoton-sulfone, Disulfoton-sulfoxide, Diuron (DCMU), DMST, DNOC, Dodemorph, Dodine, Edifenphos, Eamectin-benzoate, Endosulfan-sulfate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethiofencarb-sulfone, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoazole, Etrimfos, Famoxadone, Famphur, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenothiocarb, Fenoxaprop-ethyl, Fenoxycarb, Fepicoxamid, Fenpropathrin, Fenpropidin, Fenpropimorph, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fensulfothion-sulfone, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfide, Fipronil-sulfone, Flamprop-isopropyl, Flamprop-methyl, Flamprop-M-isopropyl, Flazasulfuron, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluacrypyrim, Fluazifop, Fluazifop-P-butyl, Fluazinam, Fluazuron, Flubendiamide, Flucythrinate, Fludioxonil, Flufenacet, Flufenacet OA, Flufenacet sulfonic acid-ESA, Flufenacet tioglycolate sulfoxide, Flufenoxuron, Flumetralin, Fluometuron,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Fluopicolide, Fluopyram, Flupyradifurone, Fluoxastrobin, Fluquinconazole, Fluridone, Flurochloridone, Fluroxypyr, Flusilazole, Fluthiacet-methyl, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Furmecyclox, Haloxypop, Haloxypop-2-ethoxyethyl, Haloxypop-methyl, Haloxypop-P-methyl, Heptenophos, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazamox, Imazapyr, Imazaquin, Imazethapyr, Imibenconazole, Imidacloprid, Indoxacarb, Iodosulfuron-methyl, Ioxynil, Ipconazole, Iprobenfos, Iprovalicarb, Isazofos, Isocarbamid, Isocarbophos, Isafenphos-methyl, Isofetamid, Isomethiozin, Isoprocab, Isopropalin, Isoprothiolane, Isoproturon, Isoprazam, Isoxaben, Isoxadifen-ethyl, Isoxaflutole, Isoxaflutole dikenotril, Isoxathion, Ivermectin B1a, Karanjin, Kresoxim-methyl, Lenacil, Linuron, Lufenuron, Malaoxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), MCPB, Mecarbam, Mefenacet, Mefenpyr-diethyl, Mefentrifluconazole, Mepanipyrim, Mephosfolan, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metazachlor ESA, Metazachlor OA, Metconazole, Methabenzthiazuron, Methacryfos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotiryne, Methoxyfenozide, Metobromuron, Metolachlor, Metolcarb (MTMC), Metominostrobin, Metoxuron, Metrafenone, Metribuzin, Metsulfuron-methyl, Mevinphos, Mexacarbate, Milbemectin A3, Milbemectin A4, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N-(2,4-Dimethylphenyl)formamide, N-(2,4-dimethylphenyl)-N-methylformamidine, N,N-diethyl-m-toluamide, Napropamide, N-desmethyl-propamocarb, Neburon, Nereistoxin oxalate, Nicosulfuron, Nitenpyram, Norflurazon, Novaluron, Nuarimol, Ofurace, Omethoate, Oryzalin, Oxadiargyl, Oxadiazon, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapirolin, Oxyacboxin, Oxydemeton-methyl, Oxymatrine, Paclobutrazol, Paraoxon-ethyl, Paraoxon-methyl, Pebulate, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthioopyrad, Pethoxamid, Phenmedipham, Phenothrin, Phenthoate, Phorate, Phorate-oxon, Phorate-oxon-sulfone, Phorate-sulfone, Phorate-sulfoxide, Phosalone, Phosfolan, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picaridin (Icaridin), Picloram, Picolinafen, Picoxystrobin, Piperonyl-butoxide, Piperophos, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Potasan, Prallethrin,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Pretilachlor, Prochloraz, Prochloraz met. BTS 44595, Prochloraz met. BTS 44596, Profenofos, Profoxydim, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propamocarb-N--oxide, Propanil, Propaquizafof, Propargite, Propazine, Propetamphos, Propham, Propiconazole, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prosulfuron, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrazoxon, Pyrethrins, Pyridaben, Pyridafol, Pyridalyl, Pyridaphenthion, Pyridate, Pyrifenox, Pyrimethanil, Pyriofenone, Pyriproxyfen, Quinalphos, Quinclorac, Quinmerac, Quinoclamine, Quinoxifen, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rabenzazole, Resmethrin, Rimsulfuron, Rotenone, Saflufenacil, Sebuthylazine, Sebumeton, Sethoxydim, Siduron, Silafluofen, Silthiofam, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spiromesifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glucoside, Spirotetramat-keto-hydroxy, Spirotetramat-mono-hydroxy, Spiroxamine, Sulcotrione, Sulfosulfuron, Sulfotep, Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebupirimfos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrione, Temephos, Tepraloxymid, Terbacil, Terbumeton, Terbuthylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiamethoxam, Thidiazuron, Thifensulfuron-methyl, Thiobencarb, Thiocyclam hydrogene oxalate, Thiodicarb, Thiofanox, Thiofanox-sulfone, Thiofanox-sulfoxide, Thiometon, Thionazin, Thiophanate-methyl, Thiram, Tolclofos-methyl, Tolfenpyrad, Tolyfluanid, Tralomethrin, Triadimefon, Triadimenol, Tri-allate, Triasulfuron, Triazamate, Triazophos, Tribenuron-methyl, Trichlorfon, Triclopyr, Tricyclazole, Trietazine, Trifloxystrobin, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triflusulfuron-methyl, Triforine, Triticonazole, Tritosulfuron, Valifenalate, Vamidothion, Vernolate, Vinclozolin, XMC, Zoxamide</i> Ukupno / In total 562				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-DB; 2,4,5-T; Bentazone				
	E4-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,005 mg/kg 1-naftalenacetamid, 2,3,5-Trimetakarb, 2,4,5-T, 2,4-D, 2,4-DB, 2,4-dimetilanilin, 3,4,5-Trimetakarb, 3-Decen-2-on, 3-hidroksikarbofuran, 4-bromofenilurea, 4-klorofenoksiocetna kiselina, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Alaklor, Alanikarb, Aldikarb, Aldikarb-sulfoksid, Aldikarb-sulfon (Aldoksikarb), Alidoklor, Ametoktradin, Ametrin, Amidosulfuron, Aminokarb, Amisulbrom, Ancimidol, Anilazin, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a (Abamektin), Azadirachtin, Azakonazol, Azametifos, Azinfos-etil, Azinfos-metil, Aziprotrin, Azoksistrobin, Benalaksil, Benfurakarb, Benodanil, Benomil, Bensulfuron-metil, Bensulid, Bentazon, Bentiavalikarb-izopropil, Bentivalikarb, Benziladenin, Benzoilprop-etil, Benzoksimat, Benzovindflupir, Bifenazat, Bifenazat-diazen, Bifentrin, Biksafen, Bispiribak-natrij, Bitertanol, Boskalid, Bromacil, Bromofos-etil, Bromoksinil, Bromukonazol, Bupirimat, Buprofezin, Butafenacil, Butaklor, Butilat, Butoksikarboksim, Butralin, Buturon, Centrimonium klorid, Cianazin, Cianfenfos, Cianofos, Ciantraniliprol, Ciazofamid, Ciflufenamid,	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	17/09/2024

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Ciflumetofen, Cihalofof-butil, Ciklanilid, Cikloat, Cikloksidim, Cikluron, Cimiazol, Cimoksanil, Cinidon-etil, Ciprodinil, Ciprokonazol, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfoksid, Demeton-S-metil-sulfon, Desmedifam, Diafentiuron, Diallylate, Diazinon, Dietofenkarb, Difenokonazole, Difenoksuron, Difetialon, Diflubenzuron, Diflufenikan, Diklobutrazol, Diklofenak, Diklofention, Diklofop, Diklofop-metil, Diklorprop, Diklorvos, Dikrotofos, Dimefuron, Dimetaklor, Dimetenamid, Dimetilan, Dimetoat, Dimetomorf (E,Z), Dimoksisstrobin, Dinikonazol, Dinoseb, Dinotefuran, Dinoterb, Dioksation, Dipropetrin, Disulfoton-sulfoksid, Disulfoton-sulfon, Diuron (DCMU), DMST, DNOC, Dodemorf, Dodin, Edifenfos, Enamektin-benzoat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etiofenkarb-sulfon, Etion, Etirimol, Etofenproks, Etofumezat, Etoksakonazol, Etoprofos, Etrimfos, Famoksadon, Fampur, Fenamidon, Fenamifos, Fenamifos-sulfoksid, Fenamifos-sulfon, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenobukarb, Fenoksaprop-etil, Fenoksikarb, Fenoprop, Fenotiokarb, Fenotrin, Fempikoksamid, Fempirazamin, Fempiroksimat, Fenpropatrin, Fenpropidin, Fensulfotion, Fensulfotion-sulfon, Fention, Fention-okson, Fention-okson-sulfoksid, Fention-okson-sulfon, Fention-sulfoksid, Fention-sulfon, Fentoat, Fenuron, Fipronil, Fipronil-sulfid, Fipronil-sulfon, Flamprop-izopropil, Flamprop-metil, Flamprop-M-izopropil, Flazasulfuron, Flazuron, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazifop-P-butyl, Fluazinam, Flubendiamid, Flucipirim, Flucitrinat, Fludioksonil, Flufenacet, Flufenacet OA, Flufenacet sulfonska kiselina-ESA, Flufenacet tioglikolat sulfoksid, Flufenoksuron, Fluksapiroksad, Flukvinkonazol, Flumetralin, Fluoksastrobin, Fluometuron, Fluopikolid, Flupiradifuron, Fluopiram, Fluridon, Fluorokloridon, Fluoroksipir, Flusilazol, Flutiacet-metil, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Foksim, Fonofos, Forat, Forat-okson, Forat-okson-sulfon, Forat-sulfoksid, Forat-sulfon, Forklorfenuron, Formetanat hidroklorid, Fosalon, Fosfamidon, Fosfolan, Fosmet, Fosmet-okson, Fostiazat, Fuberidazol, Furalaksil, Furatiokarb, Furmecikloks, Haloksifop, Haloksifop-2-etoksietil, Haloksifop-metil, Haloksifop-P-metil, Heksaflumuron, Heksakonazol, Heksazinon, Heksitiazoks, Heptenofos, Hidrametilnon, Imazakvin, Imazalil, Imazamoks, Imazapir, Imazetapir, Imibenkonazol, Imidakloprid, Indoksakarb, Ioksonil, Ipkonazol, Iprobenfos, Iprovalikarb, Ivermektin B1a, Izazofos, Izofenfos-metil, Izofetamid, Izokarbamid, Izokarbofos, Izoksaben, Izoksadifen-etil,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Izoksaf lutol, Izoksaf lutol dikenotril, Izoksation, Izometiozin, Izopirazam, Izoprokarb, Izopropalin, Izoprotiolan, Izoproturon, Jodosulfuron-metil, Kadusafos, Karanjin, Karbendazim, Karbetamid, Karbofenotion, Karbofuran, Karboksini, Karboksini sulfoksid, Karbosulfan, Karfentrazon-etil, Kinosulfuron, Kletodim, Klimbazol, Klodinafop (slobodna kiselina), Klodinafop-propargil, Klofentezin, Klokintocet-meksil, Klomazon, Klopuralid, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Kloridazon-desfenil, Klorantraniliprol, Kloroksuron, Klorotoluron, Klorpirifos, Klorpirifos-metil, Klorprofam, Klorsulfuron, Klortiamid, Klotianidin, Krezoksim-metil, Krofumat, Krotoksifos, Kumafos, Kumatetrilil, Kvinalfos, Kvinklorak, Kvinmerak, Kvinoklamin, Kinoksifen, Kvizalofop (slobodna kiselina), Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), MCPB, Mefenacet, Mefenpir-dietil, Mefentriklukonazol, Mefosfolan, Mekarbam, Meksakarbat, Mepanipirim, Mepronil, Meptildinokap, Metabenzthiazuron, Metaflumizon, Metakrifos, Metalaksil i Metalaksil-M, Metamidofos, Metamitron, Metazaklor, Metazaklor ESA, Metazaklor OA, Metidation, Metkonazol, Metobromuron, Metoksifenoizid, Metoksuron, Metolaklor, Metolcarb (MTMC), Metomil, Metominostrobin (E, Z), Metoprotin, Metrafenon, Metribuzin, Metsulfuron-metil, Mevinfos, Miklobutanol, Milbemektin A3, Milbemektin A4, Molinat, Monokrotofos, Monolinuron, Monuron, N-(2,4-dimetilfelil)-N-metilformamidin, N-(2,4-Dimetilfenil) formamid, N,N-dietil-m-toluamid, Napropamid, N-desmetil-propamokarb, Neburon, Nereistoksin oksalat, Nikosulfuron, Nikotin, Nitenpiram, Norflurazon, Novaluron, Nuarimol, Ofurak, Oksadiargil, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapiiprolin, Oksidemeton-metil, Oksikarboksini, Oksimatin, Ometoat, Orizalin, Paklobutrazol, Paraokson-etil, Paraokson-metil, Pebulat, Pencikuron, Pendimetalin, Penflufen, Penkonazol, Pentiopirad, Petoksamid, Picaridin (Icaridin), Pikloram, Pikoksistrobin, Pikolinafen, Pimetrozin, Piperofos, Piperonil-butoksid, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Pirazokson, Piretrini, Piridaben, Piridafention, Piridafol, Piridalil, Pirifeonks, Pirimetanol, Pirimifos-metil, Pirimikarb, Pirimikarb-desmetil, Piriifenon, Piriiprosifen, Potasan, Praletrin, Pretilaklor, Profam, Profenofos, Profoksidim, Prokloraz, Prokloraz met. BTS 44595, Prokloraz met. BTS 44596, Prokinazid, Promekarb, Prometon, Prometrin, Propaklor,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Propakvizafop, Propamokarb, Propamokarb-N-oksidi, Propanil, Propargit, Propazin, Propetamfos, Propikonazol, Propizamid, Propoksur, Prosulfokarb, Prosulfuron, Protiokonazol-destio, Rabenzazol, Resmetrin, Rimsulfuron, Rotenon, Saflufenacil, Sebutilazin, Sekbumeton, Setoksidim, Siduron, Silafluofen, Siltiofam, Simazin, Simetrin, Spinetoram, Spinosad (suma Spinosin A i Spinosin D), Spirodiklofen, Spiroksamin, Spiromezifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glukozid, Spirotetramat-keto-hidroksi, Spirotetramat-mono-hidroksi, Sulfoksaflo, Sulfosulfuron, Sulfotep, Sulkotriol, Sulprofos, TCMTB, Tebufenozid, Tebufenpirad, Tebukonazol, Tebupirimifos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrion, Temefos, Tepraloksidim, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetrakonazol, Tetraklorvinfos (CVMP), Tetrametrin, Tiabendazol, Tiaklopirid, Tiametoksam, Tiazuron, Tifensulfuron-metil, Tiobenkarb, Tiociklam hidrogen oksalat, Tiodikarb, Tiofanat-metil, Tiofanoks, Tiofanoks-sulfoksidi, Tiofanoks-sulfon, Tiometon, Tionazin, Tolfenpirad, Tolilfluorid, Tolklofos-metil, Triadimefon, Triadimenol, Trialat, Triazamat, Triazofos, Tribenuron-metil, Triciklazol, Trietazin, Trifloksistrobin, Triflumizol, Triflumizol metabolit FM-6-1, Triflururon, Triflururon-metil, Triflorin, Triklorpir, Triklorfon, Trisulfuron, Tritikonazol, Tritosulfuron, Valifenalat, Vamidotin, Vernolat, XMC, Zoksamid / 1-naphthaleneacetamide, 2,3,5-Trimethacarb, 2,4,5-T, 2,4-D, 2,4-DB, 2,4-dimethylaniline, 3,4,5-Trimethacarb, 3-Decen-2-one, 3-Hydroxycarbofuran, 4-Bromophenylurea, 4-Chlorophenoxyacetic acid, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Alachlor, Alanycarb, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Amidosulfuron, Aminocarb, Amisulbrom, Ancymidol, Anilazine, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a (Abamectin), Azaconazole, Azadirachtin, Azamethiphos, Azinphos-ethyl, Azinphos-methyl, Aziprotryne, Azoxystrobin, Benalaxyl, Benfuracarb, Benodanil, Benomyl, Bensulfuron-methyl, Bensulide, Bentazone, Benthiavalicarb, Benthiavalicarb-isopropyl, Benzovindflupyr, Benzoximate, Benzoylprop-ethyl, Benzyladenine, Bifenazate, Bifenazate-diazene, Bifenthrin, Bispyribac-sodium, Bitertanol, Bixafen, Boscalid, Bromacil, Bromophos-ethyl, Bromoxynil, Bromuconazole, Bupirimate, Buprofezin,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Butachlor, Butafenacil, Butoxycarboxim, Butralin, Buturon, Butylate, Cadusafos, Car bendazim, Carbetamide, Carbofuran, Carbophenothion, Carbosulfan, Carboxin, Carboxin sulfoxide, Carfentrazone-ethyl, Centrimonium chloride, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chloridazon-desphenyl, Chlorotoluron, Chloroxuron, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorsulfuron, Chlorthiamid, Cinidon-ethyl, Cinosulfuron, Clethodim, Climbazole, Clodinafop (free acid), Clodinafop-propargyl, Clofentezine, Clomazone, Clopyralid, Cloquintocet-mexyl, Clothianidin, Coumaphos, Coumatetralyl, Crotoxyphos, Crufomate, Cyanazine, Cyanofenphos, Cyanophos, Cyantraniliprole, Cyazofamid, Cyclanilide, Cycloate, Cycloxydim, Cycluron, Cyflufenamid, Cyflumetofen, Cyhalofop-butyl, Cymiazole, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Di-allate, Diazinon, Dichlofenac, Dichlofenthion, Dichlorprop, Dichlorvos, Diclobutrazol, Diclofop, Diclofop-methyl, Dicrotophos, Diethofencarb, Difenconazole, Difenoxuron, Difethialone, Diflubenzuron, Diflufenican, Dimefuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimetilan, Dimoxystrobin, Diniconazole, Dinoseb, Dinotefuran, Dinoterb, Dioxathion, Dipropetryn, Disulfoton-sulfone, Disulfoton-sulfoxide, Diuron (DCMU), DMST, DNOC, Dodemorph, Dodine, Edifenphos, Eamectin-benzoate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethiofencarb-sulfone, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrifos, Famoxadone, Famphur, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenothiocab, Fenoxaprop-ethyl, Fenoxycarb, Fenpicoxamid, Fenpropathrin, Fenpropidin, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fensulfothion-sulfone, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfide, Fipronil-sulfone, Flamprop-isopropyl, Flamprop-methyl, Flamprop-M-isopropyl, Flazasulfuron, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluacrypyrim, Fluazifop, Fluazifop-P-butyl, Fluazinam, Flubendiamide, Fluazuron, Flucythrinate, Fludioxonil, Flufenacet, Flufenacet OA, Flufenacet sulfonic acid-ESA, Flufenacet tioglycolate sulfoxide, Flufenoxuron,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Flumetralin, Fluometuron, Fluopicolide, Fluopyram, Flupyradifurone, Fluoxastrobin, Fluquinconazole, Fluridone, Flurochloridone, Fluroxypyr, Flusilazole, Fluthiacet-methyl, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Furmecyclox, Haloxypop, Haloxypop-2-ethoxyethyl, Haloxypop-methyl, Haloxypop-P-methyl, Heptenophos, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazamox, Imazapyr, Imazaquin, Imazethapyr, Imibenconazole, Imidacloprid, Indoxacarb, Iodosulfuron-methyl, Ioxynil, Ipconazole, Iprobenfos, Iprovalicarb, Isazofos, Isocarbamid, Isocarbophos, Isofenphos-methyl, Isofetamid, Isomethiozin, Isoprocab, Isopropalin, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaben, Isoxadifen-ethyl, Isoxaflutole, Isoxaflutole dikenotril, Isoxathion, Ivermectin B1a, Karanjin, Kresoxim-methyl, Linuron, Lufenuron, Malaoxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), MCPB, Mecarbam, Mefenacet, Mefenpyr-diethyl, Mefentrifluconazole, Mepanipyrim, Mephosfolan, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metazachlor ESA, Metazachlor OA, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methomyl, Methoprotryne, Methoxyfenozide, Metobromuron, Metolachlor, Metolcarb (MTMC), Metominostrobin (E, Z), Metoxuron, Metrafenone, Metribuzin, Metsulfuron-methyl, Mevinphos, Mexacarbate, Milbemectin A3, Milbemectin A4, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N-(2,4-Dimethylphenyl)formamide, N-(2,4-dimethylphenyl)-N-methylformamidine, N,N-diethyl-m-toluamide, Napropamide, N-desmethyl-propamocarb, Neburon, Nereistoxin oxalate, Nicosulfuron, Nicotine, Nitenpyram, Norflurazon, Novaluron, Nuarimol, Ofurace, Omethoate, Oryzalin, Oxadiargyl, Oxadiazon, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapiprolin, Oxyacboxin, Oxydemeton-methyl, Oxymatrine, Paclobutrazol, Paraoxon-ethyl, Paraoxon-methyl, Pebulate, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthiopyrad, Pethoxamid, Phenothrin, Phenthoate, Phorate, Phorate-oxon, Phorate-oxon-sulfone, Phorate-sulfone, Phorate-sulfoxide, Phosalone, Phosfolan, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picaridin (Icaridin), Picloram, Picolinafen, Picoxystrobin, Piperonyl-butoxide, Piperophos, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Potasan, Prallethrin, Pretilachlor,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Prochloraz, Prochloraz met. BTS 44595, Prochloraz met. BTS 44596, Profenofos, Profoxydim, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propamocarb-N--oxide, Propanil, Propaquizafop, Propargite, Propazine, Propetamphos, Propham, Propiconazole, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prosulfuron, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrazoxon, Pyrethrins, Pyridaben, Pyridafol, Pyridalyl, Pyridaphenthion, Pyrifenox (E, Z), Pyrimethanil, Pyriofenone, Pyriproxyfen, Quinalphos, Quinclorac, Quinmerac, Quinoclamine, Quinoxyfen, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rabenzazole, Resmethrin, Rimsulfuron, Rotenone, Saflufenacil, Sethuthylazine, Secbumeton, Sethoxydim, Siduron, Silafluofen, Silthiofam, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spiromesifen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glucoside, Spirotetramat-keto-hydroxy, Spirotetramat-mono-hydroxy, Spiroxamine, Sulcotrione, Sulfosulfuron, Sulfotep, Sulfoxaflo, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebupirimfos, Tebutam, Tebutiuron, Teflubenzuron, Tembotrione, Temephos, Tepraloxym, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiamethoxam, Thidiazuron, Thifensulfuron-methyl, Thiobencarb, Thiocyclam hydrogene oxalate, Thiodicarb, Thiofanox, Thiofanox-sulfone, Thiofanox-sulfoxide, Thiometon, Thionazin, Thiophanate-methyl, Tolclofos-methyl, Tolfenpyrad, Tolyfluanid, Triadimefon, Triadimenol, Tri-allate, Triasulfuron, Triazamate, Triazophos, Tribenuron-methyl, Trichlorfon, Triclopyr, Tricyclazole, Trietazine, Trifloxystrobin, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triflusulfuron-methyl, Triflorine, Triticonazole, Tritosulfuron, Valifenalate, Vamidothion, Vernolate, XMC, Zoxamide Ukupno / In total 548 Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS)				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-DB; 2,4,5-T; Bentazone				
	E4-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,005 mg/kg 1-naftalenacetamid, 2,4-D, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Alaklor, Aldikarb, Aldikarb-sulfon (Aldoksikarb), Aldikarb-sulfoksid, Alidoklor, Ametoktradin, Ametrin, Aminokarb, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a (Abamektin), Azinfos-etil, Azinfos-metil, Azoksistrobin, Benalaksil, Bendiokarb, Bentazon, Benzovindflupir, Bifenazat, Bifentrin, Bitertanol, Biksafen, Boskalid, Bromacil, Bromukonazol, Bupirimat, Butafenacil, Butoksikarboksini, Buturon, Kadusafos, Karbaril (NAC), Karbendazim, Karbetamid, Karbosulfan, Karboksini, Karfentrazon-etil, Klorantraniliprol, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Klorotoluron, Kloroksuron, Klorpirifos-metil, Kletodim, Klofentezin, Klotianidin, Kumafos, Kumatetralil, Ciantraniliprol, Ciazofamid, Cikloksidim, Cikluron, Ciprokonazol, Ciprodinil, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfon, Demeton-S-metil-sulfoksid, Desmedifam, Diafentiuuron, Diazinon, Diklofenak, Diklorprop, Diklobutrazol, Dikrotofos, Dietofenkarb, Difenokonazole, Diflubenzuron, Dimetenamid, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinotefuran, Dioksakarb, Diuron (DCMU), DMST, Dodin, Edifenfos, Enamektin-benzoat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etion, Etirimol, Etofumezat, Etoprofos, Etofenproks, Etoksazol, Etrifos, Famoksadon, Fenamidon, Fenamifos, Fenamifos-sulfon, Fenamifos-	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	21/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		sulfoksid, Fenarimol, Fenbukonazol, Fenheksamid, Fenobukarb, Fenoprop, Fenoksikarb, Fenpikoksamid, Fenpropidin, Fenpirazamin, Fenpiroksimat, Fensulfotion, Fention, Fention-okson, Fention-okson-sulfon, Fention-okson-sulfoksid, Fention-sulfon, Fention-sulfoksid, Fipronil, Fipronil-sulfon, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazinam, Fludioksonil, Flufenacet, Flufenoksuron, Fluometuron, Fluopikolid, Fluopiram, Fluoksastrobin, Flukvinkonazol, Fluridon, Flusilazol, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Fluksapiroksad, Fonofos, Forklorfenuron, Formetanat hidroklorid, Fostiazat, Fuberidazol, Furalaksil, Heksakonazol, Heksaflumuron, Heksazinon, Heksitiazoks, Hidrametilnon, Imazalil, Imazetapir, Imidaklopid, Indoksakarb, Ioksonil, Ipkonazol, Iprovalikarb, Isazofos, Izokarbofos, Izofenfos-metil, Izofetamid, Izoprokarb, Izopropalin, Izoproturon, Izopirazam, Izoksafutol, Izoksafutol-dikenotril, Ivermektin B1a, Krezoksim-metil, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), Mefenacet, Mefentriflukonazol, Mepronil, Meptildinokap, Metaflumizon, Metalaksil i Metalaksil-M, Metamitron, Metazaklor, Metkonazol, Metabenzthiazuron, Metakrifos, Metamidofos, Metidation, Metiokarb, Metiokarb-sulfon, Metiokarb-sulfoksid, Metomil, Metoprotrin, Metoksifenozyd, Metobromuron, Metolaklor, Metoksuron, Metrafenon, Metribuzin, Mevinfos, Meksakarbat, Molinat, Monokrotofos, Monolinuron, Monuron, Miklobutanil, N,N-dietil-m-toluamid, Napropamid, Neburon, Novaluron, Ometoat, Oksadiargil, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Paklobutrazol, Penkonazol, Pencikuron, Pendimetalin, Penflufen, Pentiopirad, Fenmedifam, Fenotrin, Fentoat, Forat, Fosalon, Fosmet, Fosmet-okson, Fosfamidon, Foksim, Pikoksistrobin, Pirimikarb, Pirimikarb-desmetil, Pirimifos-metil, Prokloraz, Profenofos, Promekarb, Prometon, Prometrin, Propamokarb, Propargit, Propazin, Profam, Propikonazol, Propoksir, Prosulfokarb, Protiokonazol-destio, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Piridaben, Piridalil, Piridafention, Piridat, Kvinalfos, Kvinoklamin, Kvizalofop (slobodna kiselina), Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Rotenon, Sebutilazin, Sekbumeton, Siduron, Simazin, Simetrin, Spinetoram L, Spinosad (suma Spinosin A i Spinsin D), Spiromezifen, Spirotetramat, Spirotetramat-enol-glukozid, Spirotetramat-mono-hidroksi, Spiroksamin, Sulfotep, Sulfoksaflor, Sulprofos, TCMTB, Tebukonazol, Tebufenozyd, Tebufenpirad, Tebutiuron, Teflubenzuron, Temefos, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetraklorinfos				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		(CVMP), Tetrakonazol, Tetrametrin, Tiabendazol, Tiaklopid, Tiobenkarb, Tiodikarb, Tolfenpirad, Tolilfluani, Triadimefon, Triadimenol, Triazofos, Triklorfon, Triklorpir, Triciklazol, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Tritikonazol, Tritosulfuron, Vamidotion, Vernolat, Vinklozolin, XMC, Zoksamid / <i>1-naphthaleneacetamide, 2,4-D, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Alachlor, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametoctradin, Ametryn, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a (Abamectin), Azinphos-ethyl, Azinphos-methyl, Azoxystrobin, Benalaxyl, Bendiocarb, Bentazone, Benzovindflupyr, Bifenazate, Bifenthrin, Bitertanol, Bixafen, Boscalid, Bromacil, Bromuconazole, Bupirimate, Butafenacil, Butoxycarboxim, Buturon, Cadusafos, Carbaryl (NAC), Carbendazim, Carbetamide, Carbosulfan, Carboxin, Carfentrazone-ethyl, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon, Chlorotoluron, Chloroxuron, Chlorpyrifos-methyl, Clethodim, Clofentezine, Clothianidin, Coumaphos, Coumatetralyl, Cyantraniliprole, Cyazofamid, Cycloxydim, Cycluron, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Diazinon, Dichlofenac, Dichlorprop, Diclobutrazol (stereo isomer), Dicrotophos, Diethofencarb, Difenoconazole (isomer), Diflubenzuron, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimoxystrobin, Diniconazole, Dinotefuran, Dioxacarb, Diuron (DCMU), DMST, Dodine, Edifenphos, Emamectin-benzoate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrimfos, Famoxadone, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenoxycarb, Fencicoxamid, Fenpropidin, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fipronil, Fipronil-sulfone, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Flupyrauxifen-benzyl, Fluazifop, Fluazinam, Fludioxonil, Flufenacet, Flufenoxuron, Fluometuron, Fluopicolide, Fluopyram, Fluoxastrobin, Fluquinconazole, Fluridone, Flusilazole, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazethapyr, Imidacloprid, Indoxacarb, Ioxynil, Ipconazole, Iprovalicarb, Isazofos, Isocarbophos, Isofenphos-methyl, Isofetamid, Isoprocarb, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaflutole, Isoxaflutole dikenotriple, Ivermectin B1a, Kresoxim-methyl, Lenacil, Linuron, Lufenuron, Malaaxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), Mefenacet, Mefentrifluconazole, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotryne, Methoxyfenozide, Metobromuron, Metolachlor, Metoxuron, Metrafenone, Metribuzin, Mevinphos, Mexacarbate, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N,N-diethyl-m-toluamide, Napropamide, Neburon, Novaluron, Omethoate, Oxadiargyl, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapiprolin, Oxydemeton-methyl, Paclobutrazol, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthiopyrad, Phenmedipham, Phenothrin, Phenthoate, Phorate, Phosalone, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picoxystrobin, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Prochloraz, Profenofos, Promecarb, Prometon, Prometryn, Propamocarb, Propargite, Propazine, Propham, Propiconazole, Propoxur, Prosulfocarb, Prothioconazole-desthio, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyridaben, Pyridalyl, Pyridaphenthion, Pyridate, Quinalphos, Quinoclamine, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rotenone, Sebuthylazine, Secbumeton, Siduron, Simazine, Simetryn, Spinetoram L, Spinosad (sum of spinosyn A and spinosyn D), Spiromesifen, Spirotetramat, Spirotetramat-enol-glucoside, Spirotetramat-mono-hydroxy, Spiroxamine, Sulfotep, Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebuthiuron, Teflubenzuron, Temephos, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiobencarb, Thiodicarb, Tolfenpyrad, Tolyfluanid, Triadimefon, Triadimenol, Triazophos, Trichlorfon, Triclopyr, Tricyclazole, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triticonazole, Tritosulfuron, Vamidothion, Vernolate, Vinclozolin, XMC, Zoxamide Ukupno / In total 327				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-DB; 2,4,5-T; Bentazone				
	E4-4 Povrće i žitarice s visokim udjelom škroba i/ili proteina te niskim udjelom vode i masti <i>Vegetables and cereals - high starch and/or protein content and low water and fat content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,01 mg/kg 1-naftalenacetamid, 2,4-D, 2,4-dimetilanilin, 3-hidroksikarbofuran, 4-klorofenoksioktenu kiselina, Acefat, Acetamiprid, Acetoklor, Acibenzolar-S-metil, Akrinatri, Alaklor, Aldikarb, Aldikarb-sulfon (Aldoksikarb), Aldikarb-sulfoksid, Alidoklor, Ametoktradin, Ametrin, Aminokarb, Atraton, Atrazin, Atrazin desetil, Atrazin-desizopropil, Avermektin B1a (Abamektin), Azinfos-etil, Azinfos-metil, Azoksistrobin, Benalaksil, Bendiokarb, Benfurakarb, Benomil, Bentazon, Benzovindflupir, Benzoksimat, Bifenazat, Bifentrin, Bitertanol, Biksafen, Boskalid, Bromacil, Bromoksinil, Bromukonazol, Bupirimat, Buprofezin, Butaklor, Butafenacil, Butokarboks, Butoksikarboks, Buturon, Butilat, Karbaril (NAC), Karbendazim, Karbetamid, Karbofuran, Karbofenotio, Karbosulfan, Karboks, Karfentrazon-etil, Klorantraniliprol, Klorbromuron, Klordimeform, Klorfenvinfos (E,Z), Klorflurazon, Kloridazon, Klorotoluron, Kloroksuron, Klorprofam, Klorpirifos, Klorpirifos-metil, Kletodim, Klotofentzin, Klotiazon, Klotianidin, Kumafos, Kumatetralil, Cianazin, Ciantraniliprol, Ciazofamid,	LC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	07/12/2018	21/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Cikloat, Cikloksidim, Cikluron, Ciflufenamid, Cimoksanil, Ciprokonazol, Ciprodinil, Ciromazin, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfon, Demeton-S-metil-sulfoksid, Desmedifam, Diafentiuron, Diallate, Diazinon, Diklofenak, Diklorprop, Diklorvos, Diklobutrazol, Dikrotofos, Dietofenkarb, Difenokonazole, Difetialon, Diflubenzuron, Dimetaklor, Dimetenamid, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin, Dinikonazol, Dinotefuran, Dioksakarb, Difenamid, Diuron (DCMU), DMST, Dodin, Edifenfos, Enamektin-benzoat, Epoksikonazol, EPTC, Etakonazol, Etiofenkarb, Etion, Etrimol, Etofumezat, Etoprofos, Etofenproks, Etoksakonazol, Etrimfos, Famoksadon, Fenamidon, Fenamifos, Fenamifos-sulfon, Fenamifos-sulfoksid, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenobukarb, Fenoprop, Fenoksikarb, Fenpikoksamid, Fenpropatrin, Fenpropidin, Fenpirazamin, Fenpiroksimat, Fensulfotion, Fention, Fention-okson, Fention-okson-sulfon, Fention-okson-sulfoksid, Fention-sulfon, Fention-sulfoksid, Fenuron, Fipronil, Fipronil-sulfon, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluazifop, Fluazifop-P-butil, Fluazinam, Flucitrinat, Fludioksonil, Flufenacet, Flufenoksuron, Fluometuron, Fluopikolid, Fluopiram, Fluoksastrobin, Flukvinkonazol, Fluridon, Flusilazol, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Fluksapiroksad, Fonofos, Forklorfenuron, Formetanat hidroklorid, Fostiazat, Fuberidazol, Furalaksil, Furatiokarb, Haloksifop, Heksakonazol, Heksaflumuron, Heksazinon, Heksitiazoks, Hidrametilnon, Imazalil, Imidaklopid, Indoksakarb, Ioksonil, Ipkonazol, Iprovalikarb, Izazofos, Izofenfos-metil, Izofetamid, Izoprokarb, Izopropalin, Izoprotiolan, Izoproturon, Izopirazam, Izoksafutol, Izoksafutol-dikenotril, Ivermektin B1a, Krezoksim-metil, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, MCPA (MCP), Mefenacet, Mefentriflukonazol, Mepanipirim, Mepronil, Meptildinokap, Metaflumizon, Metalaksil i Metalaksil-M, Metamitron, Metazaklor, Metkonazol, Metabenzthiazuron, Metakrifos, Metamidofos, Metidation, Metiokarb, Metiokarb-sulfon, Metiokarb-sulfoksid, Metomil, Metoprotrin, Metoksifenozyd, Metobromuron, Metolaklor, Metoksuron, Metrafenon, Metribuzin, Mevinfos, Meksakarbat, Molinat, Monokrotofos, Monolinuron, Monuron, Miklobutanil, N,N-dietil-m-toluamid, Naled, Napropamid, Neburon, Nereistoksin oksalat, Nikotin, Nitenpiram, Norflurazon, Novaluron, Nuarimol, Ometoat, Oksadiargil, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapiiprolin, Oksidemeton-metil, Paklobutrazol, Pebulat, Penkonazol, Pencikuron,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Pendimetalin, Penflufen, Pentiopirad, Fenmedifam, Fenotrin, Fentoat, Forat, Fosalon, Fosmet, Fosmet-okson, Fosfamidon, Foksim, Pikolinafen, Pikoksistrobin, Piperonil-butoksid, Pirimikarb, Pirimikarb-desmetil, Pirimifosmetil, Pretilaklor, Prokloraz, Profenofos, Promekarb, Prometon, Prometrin, Propaklor, Propamokarb, Propanil, Propakvizafop, Propargit, Propazin, Profam, Propikonazol, Propoksur, Propizamid, Prokinazid, Prosulfokarb, Protiokonazol-destio, Pimetrozin, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Piretrini, Piridaben, Piridalil, Piridafention, Piridat, Pirimetanil, Piriopenon, Pirioproksifen, Kvinalfos, Kvinoklamin, Kinoksifen, Kvizalofofop (slobodna kiselina), Kvizalofofop-etil, Kvizalofofop-metil, Kvizalofofop-P-etil, Kvizalofofop-P-tefuril, Resmetrin, Rotenon, Sebutilazin, Sekbumeton, Siduron, Silafluofen, Simazin, Simetrin, Spinetoram, Spinosad (suma Spinosin A i Spinsin D), Spirodiklofen, Spiromezifen, Spirotetramat, Spirotetramat-enol-glukozid, Spirotetramat-mono-hidroksi, Spiroksamin, Sulfotep, Sulfoksaflor, Sulprofos, TCMTB, Tebukonazol, Tebufenozid, Tebufenpirad, Tebutiuron, Teflubenzuron, Temefos, Terbacil, Terbumeton, Terbutilazin, Terbutrin, Tetraklorinfos (CVMP), Tetrakonazol, Tetrametrin, Tiabendazol, Tiaklopid, Tiametoksam, Tidiazuron, Tiobenkarb, Tiociklam hidrogen oksalat, Tiodikarb, Tiometon, Tiofanat-metil, Tiram, Tolklofos-metil, Tolfenpirad, Tralometrin, Triadimefon, Triadimenol, Trialat, Triazofos, Triklorfon, Triklorpir, Triciklazol, Trifloksistrobin, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Tritikonazol, Tritosulfuron, Vamidotion, Vernolat, Vinklozolin, XMC, Zoksamid / <i>1-naphthaleneacetamide, 2,4-D, 2,4-dimethylaniline, 3-Hydroxycarbofuran, 4-Chlorophenoxyacetic acid, Acephate, Acetamiprid, Acetochlor, Acibenzolar-S-methyl, Acrinathrin, Alachlor, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor, Ametocetradin, Ametryn, Aminocarb, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a (Abamectin), Azinphos-ethyl, Azinphos-methyl, Azoxystrobin, Benalaxyl, Bendiocarb, Benfuracarb, Benomyl, Bentazone, Benzovindflupyr, Benzoximate, Bifenazate, Bifenthrin, Bitertanol, Bixafen, Boscalid, Bromacil, Bromoxynil, Bromuconazole, Bupirimate, Buprofezin, Butachlor, Butafenacil, Butocarboxim, Butoxycarboxim, Buturon, Butylate, Carbaryl (NAC), Carbazim, Carbetamide, Carbofuran, Carbophenothion, Carbosulfan, Carboxin, Carfentrazon-ethyl, Chlorantraniliprole, Chlorbromuron, Chlordimeform, Chlorfenvinphos (E, Z), Chlorfluazuron, Chloridazon,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Chlorotoluron, Chloroxuron, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Clethodim, Clofentezine, Clomazone, Clothianidin, Coumaphos, Coumatetralyl, Cyanazine, Cyantraniliprole, Cyazofamid, Cycloate, Cycloxydim, Cycluron, Cyflufenamid, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Desmedipham, Diafenthiuron, Di-allate, Diazinon, Dichlofenac, Dichlorprop, Dichlorvos, Diclobutrazol, Dicrotophos, Diethofencarb, Difenconazole, Difethialone, Diflubenzuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimoxystrobin, Diniconazole, Dinotefuran, Dioxacarb, Diphenamid, Diuron (DCMU), DMST, Dodine, Edifenphos, Emamectin-benzoate, Epoxiconazole, EPTC, Etaconazole, Ethiofencarb, Ethion, Ethirimol, Ethofumesate, Ethoprophos, Etofenprox, Etoxazole, Etrifos, Famoxadone, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenoprop, Fenoxycarb, Fenpicoxamid, Fenpropathrin, Fenpropidin, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fenthion, Fenthion-oxon, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfone, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfone, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluazifop, Fluazifop-P-butyl, Fluazinam, Flucythrinate, Fludioxonil, Flufenacet, Flufenoxuron, Fluometuron, Fluopicolide, Fluopyram, Fluoxastrobin, Fluquinconazole, Fluridone, Flusilazole, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Fonofos, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Haloxyfop, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imidacloprid, Indoxacarb, Ioxynil, Ipconazole, Iprovalicarb, Isazofos, Isafenphos-methyl, Isofetamid, Isoprocarb, Isopropalin, Isoprothiolane, Isoproturon, Isopyrazam, Isoxaflutole, Isoxaflutole dikenotriple, Ivermectin B1a, Kresoxim-methyl, Lenacil, Linuron, Lufenuron, Malaoxon, Malathion, Mandipropamid, Matrine, MCPA (MCP), Mefenacet, Mefentrifluconazole, Mepanipyrim, Mepronil, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metamitron, Metazachlor, Metconazole, Methabenzthiazuron, Methacrifos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotetryne, Methoxyfenozide, Metobromuron, Metolachlor, Metoxuron, Metrafenone, Metribuzin, Mevinphos, Mexacarbate, Molinate, Monocrotophos,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Monolinuron, Monuron, Myclobutanil, N,N-diethyl-m-toluamide, Naled, Napropamide, Neburon, Nereistoxin oxalate, Nicotine, Nitenpyram, Norflurazon, Novaluron, Nuairimol, Omethoate, Oxadiargyl, Oxadiazon, Oxadixyl, Oxamyl, Oxa-sulfuron, Oxathiapiprolin, Oxydemeton-methyl, Paclobutrazol, Pebulate, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthio-pyrad, Phenmedipham, Phenothrin, Phenthoate, Phorate, Phosalone, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picolinafen, Picoxystrobin, Piperonyl-butoxide, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Pretilachlor, Prochloraz, Profenofos, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propanil, Propaquizafop, Propargite, Propazine, Propham, Propiconazole, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrethrins, Pyridaben, Pyridalyl, Pyridaphenthion, Pyridate, Pyrimethanil, Pyriofenone, Pyriproxyfen, Quinalphos, Quinoclamine, Quinoxifen, Quizalofop (free acid), Quizalofop-ethyl, Quizalofop-methyl, - Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Resmethrin, Rotenone, Se-buthylazine, Se-cbumeton, Siduron, Silafluofen, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spiro-diclofen, Spiromesifen, Spirotetramat, Spirotetramat-enol-glucoside, Spirotetramat-mono-hydroxy, Spiroxamine, Sulfotep, Sulfoxaflo-r, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebuthiuron, Teflubenzuron, Temephos, Terbacil, Terbumeton, Terbutylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiabendazole, Thiacloprid, Thiamethoxam, Thidiazuron, Thiobencarb, Thiocyclam hydrogene oxalate, Thiodicarb, Thiometon, Thiophanate-methyl, Thiram, Tolclofos-methyl, Tolfenpyrad, Tralomethrin, Triadimefon, Triadimenol, Tri-allate, Triazophos, Trichlorfon, Triclopyr, Tricyclazole, Trifloxystrobin, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triticonazole, Tritosulfuron, Vamidothion, Vernolate, Vinclozolin, XMC, Zoxamide Ukupno/ In Total: 392 Određivanje kiselih pesticida nakon alkalne hidrolize u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS)				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Determination of acidic pesticide residues following alkaline hydrolysis in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fluazifop; Haloxyfop; 2,4-D; Dichlorprop; MCPA; MCPB; Quizalofop; 2,4-DB; 2,4,5-T; Bentazone				
	E4-5 "Teški matriksi" <i>"Difficult or unique commodities"</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) verzija A <i>Method for the determination of pesticide residues in foods of plant origin – liquid chromatographic method with mass spectrometry (LC-MS/MS) version A</i> 0,01 mg/kg 1-naftalacetamid, 2,3,5-trimetakarb, 2,4,5-T, 3,4,5-trimetakarb, 3-Decen-2-on, Acefat, Acibenzolar-S-metil, Alaklor, Alanikarb, Aldikarb, Aldikarb-sulfoksid, Aldikarb-sulfon (Aldoksikarb), Alidoklor, Ametoktradin, Ametrin, Amidosulfuron, Amisulbrom, Atraton, Atrazin, Atrazin-desetil, Atrazin-desizopropil, Avermektin B1a (Abamektin), Azakonazol, Azametifos, Azinfos-metil, Aziprotrin, Azoksistrobin, Benalaksil, Bendiokarb, Benodanil, Bensulfuron-metil, Bensulid, Bentazon, Benthiavalikarb, Benthiavalikarb-izopropil, Benziladenin, Benzoilprop-etil, Benzoksimat, Benzovindflupir, Bifenazat, Bifenazat-diazen, Bifentrin, Biksafen, Bitertanol, Boskalid, Bromacil, Bromukonazol, Bupirimat, Butafenacil, Butilat, Butoksikarboksime, Butralin, Buturon, Centrimonijev klorid, Ciflufenamid, Ciflumetofen, Cijanazin, Cijanofenfos, Ciantraniliprol, Ciazofamid, Ciklanilid, Cikloksidim, Cikluron, Cimiazol, Cimoksanil, Cinosulfuron, Ciprodinil, Ciprokonazol, Ciromazin, Cyazofamid, Deltametrin, Demeton-O, Demeton-S-metil, Demeton-S-metil-sulfoksid, Demeton-S-metil-sulfon, Diazinon, Dietofenkarb, Difenamid, Difenokonazol, Difenoksuron, Diflubenzuron, Diflufenikan, Dijafentiuron, Diklobutrazol, Diklofenak, Dikrotofos, Dimefuron, Dimetaklor, Dimetenamid, Dimetilan, Dimetoat, Dimetomorf (E,Z), Dimoksistrobin,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Dinikonazol, Dinoseb, Dinotefuran, Dinoterb, Dioksakarb, Dipropetrin, Disulfoton-sulfoksid, Disulfoton-sulfon, Diuron (DCMU), DMST, Dodemorf, Dodin, Edifenfos, Eamektin-benzoat, Epoksikonazol, Etakonazol, Etiofenkarb, Etion, Etirimol, Etofenproks, Etoksazol, Etoprofos, Etrimfos, Famoksadon, Famphur, Fenamidon, Fenamifos, Fenamifos-sulfoksid, Fenamifos-sulfon, Fenarimol, Fenazakvin, Fenbukonazol, Fenheksamid, Fenmedifam, Fenobukarb, Fenoksaprop-etil, Fenoksikarb, Fenotiokarb, Fenotrin, Fepikoksamid, Fepirazamin, Fepiroksimat, Fenpropimorf, Fensulfotion, Fensulfotion-sulfon, Fention-okson-sulfoksid, Fention-okson-sulfon, Fention-sulfoksid, Fenuron, Fipronil, Fipronil-sulfid, Fipronil-sulfon, Flamprop-izopropil, Flamprop-metil, Flamprop-M-izopropil, Flazasulfuron, Flonikamid, Flonikamid TFNA, Flonikamid TFNG, Florpirauksifen-benzil, Fluakripirim, Fluazifop-P-butil, Fluazinam, Fluazuron, Fludioksonil, Flufenacet, Flufenoksuron, Fluksapiroksad, Flukvinkonazol, Fluoksastrobin, Fluopiram, Fluridon, Flusilazol, Flutiacet-metil, Flutianil, Flutolanil, Flutriafol, Fluvalinat, Foksim, Forat, Forat-okson, Forat-okson-sulfon, Forat-sulfoksid, Forat-sulfon, Forklorfenuron, Formetanat hidroklorid, Fosalon, Fosfamidon, Fosfolan, Fosmet, Fosmet-okson, Fostiazat, Fuberidazol, Furalaksil, Furathiokarb, Furmecikloks, Haloksifop-2-etoksietil, Heksaflumuron, Heksakonazol, Heksazinon, Heksitiazoks, Heptenofos, Hidrametilnon, Imazakvin, Imazalil, Imibenkonazol, Imidaklopid, Indoksakarb, Ipkonazol, Iprobenfos, Iprovalikarb, Isazofos, Isofetamid, Isoksaben, Isoprocarb, Ivermektin B1a, Izokarbamid, Izoksadifen-etil, Izoksaflutol, Izoksaflutol dikenotril, Izoksation, Izometiozin, Izopirazam, Izoprotiolan, Izoproturon, Jodosulfuron-metil, Kadusafos, Karanjin, Karbaril (NAC), Karbendazim, Karbetamid, Karbofuran, Karboksinsulfoksid, Karbosulfan, Klethodim, Klimbazol, Klodinafop (slobodna kiselina), Klodinafop-propargil, Klofentezin, Klokintocet-meksil, Klomazon, Klorantraniliprol, Klorbromuron, Klorfenvinfos (E, Z), Klorfluazuron, Kloridazon, Kloroksuron, Klorosulfuron, Klotianidin, Krotoksifos, Krufomat, Kumafos, Kumatetralil, Kvinalfos, Kinoksifen, Kvizalofop-etil, Kvizalofop-metil, Kvizalofop-P-etil, Kvizalofop-P-tefuril, Lenacil, Linuron, Lufenuron, Malaokson, Malation, Mandipropamid, Matrin, Mefenacet, Mefenpir-dietil, Mefentriklonazol, Mefosfolan, Mekarbam, Meksakarb, Mepanipirim, Meptildinokap, Metabenziazuron, Metaflumizon, Metakrifos, Metalaksil i Metalaksil-M, Metamidofos, Metazaklor, Metazaklor OA, Metidation,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Metiokarb, Metiokarb-sulfoksid, Metiokarb-sulfon, Metkonazol, Metobromuron, Metoksifenoizid, Metoksuron, Metolaklor, Metolkarb (MTMC), Metomil, Metominostrobin (E, Z), Metoprotrin, Metrafenon, Metribuzin, Metsulfuron-metil, Miklobutanil, Molinat, Monokrotofos, Monolinuron, Monuron, N-(2,4-dimetilfenil)formamid, N,N-dietil-m-toluamid, Napropamid, N-desmetil-propamokarb, Neburon, Nikosulfuron, Nitenpiram, Norflurazon, Novaluron, Nuairimol, Ofurak, Oksadiazon, Oksadiksil, Oksamil, Oksasulfuron, Oksatiapirolin, Oksidemeton-metil, Oksikarboksin, Ometoat, Paklobutrazol, Paraokson-etil, Pencikuron, Pendimetalin, Penflufen, Penkonazol, Pentiopirad, Petoksamid, Pikaridin (Ikaridin), Pikoksistrobin, Pikolinafen, Pimetrozin, Piperofos, Piperonil-butoksid, Pirakarbolid, Piraklofos, Piraklostrobin, Pirazofos, Pirazokson, Piridaben, Piridafention, Piridalil, Piridat, Pirimetanil, Pirimicarb, Pirimifos-metil, Pirimikarb-desmetil, Piriofenon, Piriproksifen, Potasan, Praletrin, Pretlaklor, Profenofos, Profoksidim, Prokloraz, Prokloraz met. BTS 44595, Prokloraz met. BTS 44596, Prokinazid, Promecarb, Prometon, Prometrin, Propaklor, Propakvizafop, Propamocarb, Propamokarb-N-oksidi, Propargit, Propazin, Propetamfos, Propizamid, Propoksur, Prosulfokarb, Prosulfuron, Protiokonazol-destio, Rabenzazol, Rimsulfuron, Rotenon, Saflufenacil, Sebutilazin, Sekbumeton, Setoksidim, Siduron, Silthiofam, Simazin, Simetrin, Spinetoram, Spinosad (zbroj spinosina A i spinosina D), Spirodiklofen, Spiroksamin, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glukozid, Spirotetramat-keto-hidroksi, Spirotetramat-mono-hidroksi, Sulfoksaflor, Sulfosulfuron, Sulfotep, Sulprofos, TCMTB, Tebufenoizid, Tebufenpirad, Tebukonazol, Tebupirimfos, Tebutam, Tebutiuron, Teflubenzuron, Temephos, Tepraloksidim, Terbumeton, Terbutilazin, Terbutrin, Tetraklorvinfos (CVMP), Tetrakonazol, Tetrametrin, Thionazin, Tiaklopid, Tiametoksam, Tidiazuron, Tifensulfuron-metil, Tiobenkarb, Tiodikarb, Tiofanat-metil, Tiofanoks-sulfoksid, Tolfenpirad, Triadimefon, Triadimenol, Triasulfuron, Triazamat, Triazofos, Tribenuron-metil, Triciklazol, Trietazin, Triflumizol, Triflumizol metabolit FM-6-1, Triflumuron, Triflusulfuron-metil, Triforin, Tritikonazol, Tritosulfuron, Valifenalat, Vamidotion, XMC, Zoksamid / <i>1-naphthaleneacetamide, 2,3,5-Trimethacarb, 2,4,5-T, 3,4,5-Trimethacarb, 3-Decen-2-one, Acephate, Acibenzolar-S-methyl, Alachlor, Alanycarb, Aldicarb, Aldicarb-sulfone (Aldoxycarb), Aldicarb-sulfoxide, Allidochlor,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Ametoctradin, Ametryn, Amidosulfuron, Amisulbrom, Atraton, Atrazine, Atrazine-desethyl, Atrazine-desisopropyl, Avermectin B1a (Abamectin), Azaconazole, Azamethiphos, Azinphos-methyl, Aziprotryne, Azoxystrobin, Benalaxyl, Bendiocarb, Benodanil, Bensulfuron-methyl, Bensulide, Bentazone, Benthiavalicarb, Benthiavalicarb-isopropyl, Benzovindflupyr, Benzoximate, Benzoylprop-ethyl, Benzyladenine, Bifenazate, Bifenazate-diazene, Bifenthrin, Bitertanol, Bixafen, Boscalid, Bromacil, Bromuconazole, Bupirimate, Butafenacil, Butoxycarboxim, Butralin, Buturon, Butylate, Cadusafos, Carbaryl (NAC), Carbendazim, Carbetamide, Carbofuran, Carbosulfan, Carboxin sulfoxide, Centrimonium chloride, Chlorantraniliprole, Chlorbromuron, Chlorfenvinphos (E, Z), Chlorfluzuron, Chloridazon, Chloroxuron, Chlorsulfuron, Cinosulfuron, Clethodim, Climbazole, Clodinafop (free acid), Clodinafop-propargyl, Clofentezine, Clomazone, Cloquintocet-mexyl, Clothianidin, Coumaphos, Coumatetralyl, Crotoxyphos, Crufomate, Cyanazine, Cyanofenphos, Cyantraniliprole, Cyazofamid, Cyclanilide, Cycloxydim, Cycluron, Cyflufenamid, Cyflumetofen, Cymiazole, Cymoxanil, Cyproconazole, Cyprodinil, Cyromazine, Deltamethrin, Demeton-O, Demeton-S-methyl, Demeton-S-methyl-sulfone, Demeton-S-methyl-sulfoxide, Diafenthiuron, Diazinon, Dichlofenac, Diclobutrazol, Dicrotophos, Diethofencarb, Difenconazole, Difenoxuron, Diflubenzuron, Diflufenican, Dimefuron, Dimethachlor, Dimethenamide, Dimethoate, Dimethomorph (E,Z), Dimetilan, Dimoxystrobin, Diniconazole, Dinoseb, Dinotefuran, Dinoterb, Dioxacarb, Diphenamid, Dipropetryn, Disulfoton-sulfone, Disulfoton-sulfoxide, Diuron (DCMU), DMST, Dodemorph, Dodine, Edifenphos, Emamectin-benzoate, Epoxiconazole, Etaconazole, Ethiofencarb, Ethion, Ethirimol, Ethoprophos, Etofenprox, Etoxazole, Etrinfos, Famoxadone, Famphur, Fenamidone, Fenamiphos, Fenamiphos-sulfone, Fenamiphos-sulfoxide, Fenarimol, Fenazaquin, Fenbuconazole, Fenhexamid, Fenobucarb, Fenothiocarb, Fenoxaprop-ethyl, Fenoxycarb, Fencicoxamid, Fenpropimorph, Fenpyrazamine, Fenpyroximate, Fensulfothion, Fensulfothion-sulfone, Fenthion-oxon-sulfone, Fenthion-oxon-sulfoxide, Fenthion-sulfoxide, Fenuron, Fipronil, Fipronil-sulfide, Fipronil-sulfone, Flamprop-isopropyl, Flamprop-methyl, Flamprop-M-isopropyl, Flazasulfuron, Flonicamid, Flonicamid TFNA, Flonicamid TFNG, Florpyrauxifen-benzyl, Fluacrypyrim, Fluazifop-P-butyl, Fluazinam, Fluazuron, Fludioxonil, Flufenacet, Flufenoxuron, Fluopyram,</i>				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Fluoxastrobin, Fluquinconazole, Fluridone, Flusilazole, Fluthiacet-methyl, Flutianil, Flutolanil, Flutriafol, Fluvalinate, Fluxapyroxad, Forchlorfenuron, Formetanate hydrochloride, Fosthiazate, Fuberidazole, Furalaxyl, Furathiocarb, Furmecyclox, Haloxypop-2-ethoxyethyl, Heptenophos, Hexaconazole, Hexaflumuron, Hexazinone, Hexythiazox, Hydramethylnon, Imazalil, Imazaquin, Imibenconazole, Imidacloprid, Indoxacarb, Iodosulfuron-methyl, Ipconazole, Iprobenfos, Iprovalicarb, Isazofos, Isocarbamid, Isofetamid, Isomethiozin, Isoprocab, Isoprothiolane, Isoproturon, Isoprazam, Isoxaben, Isoxadifen-ethyl, Isoxaflutole, Isoxaflutole dikenotriple, Isoxathion, Ivermectin B1a, Karanjin, Lenacil, Linuron, Lufenuron, Malaoxon, Malathion, Mandipropamid, Matrine, Mecarbam, Mefenacet, Mefenpyr-diethyl, Mefentrifluconazole, Mepanipyrim, Mephosfolan, Meptyldinocap, Metaflumizone, Metalaxyl and Metalaxyl-M, Metazachlor, Metazachlor OA, Metconazole, Methabenzthiazuron, Methacryfos, Methamidophos, Methidathion, Methiocarb, Methiocarb-sulfone, Methiocarb-sulfoxide, Methomyl, Methoprotiryne, Methoxyfenozide, Metobromuron, Metolachlor, Metolcarb (MTMC), Metominostrobin (E, Z), Metoxuron, Metrafenone, Metribuzin, Metsulfuron-methyl, Mexacarbate, Molinate, Monocrotophos, Monolinuron, Monuron, Myclobutanil, N-(2,4-Dimethylphenyl)formamide, N,N-diethyl-m-toluamide, Napropamide, N-desmethyl-propamocarb, Neburon, Nicosulfuron, Nitenpyram, Norflurazon, Novaluron, Nuairimol, Ofurace, Omethoate, Oxadiazon, Oxadixyl, Oxamyl, Oxasulfuron, Oxathiapiprolin, Oxycarboxin, Oxydemeton-methyl, Paclobutrazol, Paraaxon-ethyl, Penconazole, Pencycuron, Pendimethalin, Penflufen, Penthioopyrad, Pethoxamid, Phenmedipham, Phenothrin, Phorate, Phorate-oxon, Phorate-oxon-sulfone, Phorate-sulfone, Phorate-sulfoxide, Phosalone, Phosfolan, Phosmet, Phosmet-oxon, Phosphamidon, Phoxim, Picaridin (Icaridin), Picolinafen, Picoxystrobin, Piperonyl-butoxide, Piperophos, Pirimicarb, Pirimicarb-desmethyl, Pirimiphos-methyl, Potasan, Prallethrin, Pretilachlor, Prochloraz, Prochloraz met. BTS 44595, Prochloraz met. BTS 44596, Profenofos, Profoxydim, Promecarb, Prometon, Prometryn, Propachlor, Propamocarb, Propamocarb-N--oxide, Propaquizafop, Propargite, Propazine, Propetamphos, Propoxur, Propyzamide, Proquinazid, Prosulfocarb, Prosulfuron, Prothioconazole-desthio, Pymetrozine, Pyracarbolid, Pyraclofos, Pyraclostrobin, Pyrazophos, Pyrazoxon, Pyridaben, Pyridalyl, Pyridaphenthion, Pyridate, Pyrimethanil,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Pyriofenone, Pyriproxyfen, Quinalphos, Quinoxifen, Quizalofop-ethyl, Quizalofop-methyl, Quizalofop-P-ethyl, Quizalofop-P-tefuryl, Rabenzazole, Rimsulfuron, Rotenone, Saflufenacil, Sebuthylazine, Secbumeton, Sethoxydim, Siduron, Silthiofam, Simazine, Simetryn, Spinetoram, Spinosad (sum of spinosyn A and spinosyn D), Spirodiclofen, Spirotetramat, Spirotetramat-enol, Spirotetramat-enol-glucoside, Spirotetramat-keto-hydroxy, Spirotetramat-mono-hydroxy, Spiroxamine, Sulfosulfuron, Sulfotep, Sulfoxaflor, Sulprofos, TCMTB, Tebuconazole, Tebufenozide, Tebufenpyrad, Tebupirimfos, Tebutam, Tebuthiuron, Teflubenzuron, Temephos, Tepraloxym, Terbumeton, Terbuthylazine, Terbutryn, Tetrachlorvinphos (CVMP), Tetraconazole, Tetramethrin, Thiacloprid, Thiamethoxam, Thidiazuron, Thifensulfuron-methyl, Thiobencarb, Thiodicarb, Thiofanox-sulfoxide, Thionazin, Thiophanate-methyl, Tolfenpyrad, Triadimefon, Triadimenol, Triasulfuron, Triazamate, Triazophos, Tribenuron-methyl, Tricyclazole, Trietazine, Triflumizole, Triflumizole metabolite FM-6-1, Triflumuron, Triflusulfuron-methyl, Triforine, Triticonazole, Tritosulfuron, Valifenalate, Vamidothion, XMC, Zoxamide Ukupno/ In Total: 435				
E5	E5-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 1-Naftiloctena kiselina, 2,3,5,6-Tetrakloroanilin, 2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Akrinatrín (suma izomera), Alaklor, Aldrin, Aletrin -3,4 (Bioaletrin), Alidoklor, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Benfluralin, Benziladenin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Binapacryl, Bitrex (Denatonium benzoat), Bromacil, Bromfeninfos-metil, Bromfeninfos, Bromociklen, Bromofos, Bromofos-etil,	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	17/09/2024

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Bromopropilat, Kaptan, Karbofenotion, Karbofenotion metil, Kinometionat, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorobenzilat, Kloroneb, Klortalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Klordan, cis-Permetrin, Klodinafop propargil, Klomazon, Klokintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cinerin I, Cinerin II, Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Demeton-S-metil, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklofluanid, Diklorvos, Diklobutrazol, Dikloran, Dikofol, Dieltrin, Dimefoks, Dimetaklor, Difenamid, Dimetipin, Dimetilan, Dinitramin, Difenilamin, Disulfoton, Ditalimfos, Edifenfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Etrimfos, Famfur, Fenarimol, Fenklorfos, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Fluakloralin, Fludioksonil, Fluensulfon, Flumetrin, Flumioksazin, Fluotrimazol, Flupiradifuron, Fluorokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutriafol, Fluvalinat-tau (suma izomera), Folpet, Formotion, Furilazol, gama-BHC (lindan), Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izodrin, Izoksaben, Izoksadifen-etil, Izoksation, Jasmolin I, Jasmolin II, lambda-Cihalotrin, Lenacil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanol, N-(2,4-dimetilfenil) formamid, Nitralin, Nitrofen, Nitrotral-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksikloridan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Paration-metil, Pebulat, Penkonazol, Pendimetalin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Protiokonazol, Protiofos, Piraklofos, Piraflufenetil, Piridaben, Piridafention, Pirimetanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametarin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbamil, Tolklofos-metil, Tolilfluanid, Tolilfluanid metabolit, trans-Klordan,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Triciklazol, Tridemorf, Trietazin, Trifluralin, Unikonzol, Valifenalat, Vernolat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 1-Naphthylacetic acid,, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allethrin-3,4 (Bioallethrin), Allidochlor, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Benfluralin, Benzyladenine, beta-BHC, beta-Endosulfan, Bifenthrin, Binapacryl, Biphenyl, Binapacryl, Bitrex (Denatonium benzoate), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Captan, Carbophenothion, Carbophenothion methyl, Chinomethionat, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolate, Cinerin I, Cinerin II, Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Clomazone, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Demeton-S-methyl, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlofluanid, Dichlorvos, Diclobutrazol, Dicloran, Dicofol, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Dinitramine, Diphenamid, Diphenylamine, Disulfoton, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Etrimfos, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfathion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Flumioxazine, Fluotrimazol, Flupyradifurone, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutriafol, Fluvalinate (sum of isomers), Folpet, Formothion, Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Iprobenfos, Iprodione, Isobenzan, Isodrin, Isoxaben, Isoxadifen ethyl, Isoxathion, Jasmolin I, Jasmolin II, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Ofurace, Paclobutrazol, Parathion, Parathion-methyl, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phorate, Phosfolan,, Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prothioconazole, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tricyclazole, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vernolate, Vinclozolin</i> Ukupno / In total 270				
	E5-2 Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroaniline, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Aktrinatin (suma izomera), Alaklor,	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Aldrin, Aletrin -3,4 (Bioaletrin), alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Azinfos-methyl, Benfluralin, Benziladenin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Bupirimat, Karbofenotion, Karbofenotion metil, Kinometionat, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorbenzid, Klorobenzilat, Kloroneb, Klortalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Cinidon etil, cis-Klordan, cis-Permetrin, Klodinafop propargil, Klok-vintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklorvos, Diklobutrazol, Dikloran, Dikofol, Dieldrin, Dimefoks, Dimektaklor, Difenamid, Dimetipin, Dimetilan, Dinitramin, Difenilamin, Disulfoton, Ditalimfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Famfur, Fenarimol, Fenklorfos, Fenitrotion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutriafol, Fluvalinat-tau (suma izomera), Fonofos, Formotion, Fostiazat (suma izomera), Furilazol, gama-BHC (lindan), Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izofenos metil, Izoksaben, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, Nitalin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Orizalin, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Paration-metil, Pebulat, Penkonazol, Pendimetalin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Fosfolan, Fosmet, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Protiokonazol, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Spiromezifen, Sulkotrin, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflutrin,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbazil, Tolklofos-metil, Tolilflu-anid, Tolilflu-anid metabolit, trans-Klordan, Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Tridemorf, Trietazin, Trifluralin, Unikonazol, Valifenalat, Vernolat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifeno, (Z)-Pyrifeno, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-Tetrachloroaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin,, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allethrin-3,4 (Bioallethrin), alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Azinphos methyl, Benfluralin, Benzyladenine, beta- BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Binapacryl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-methyl, Bromfenvinfos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Carbophenothion, Carbophenothion methyl, Chinomethionat, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlorvos, Diclobutrazol, Dicloran, Dicofol, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Dinitramine, Diphenamid, Diphenylamine, Disulfoton, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfuralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutriafol, Fluvalinate (sum of isomers), Folpet, Formothion, Fosthiazate (sum of isomers), Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Iprodione, Isobenzan, Isocarbofos, Isodrin, Isofenphos-methyl, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, Nitralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychloran, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Ofurace, Oryzalin, Paclobutrazol, Parathion, Parathion-methyl, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phosfolan, Phosmet, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prothioconazole, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesifen, Sulcotrione, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vernolate, Vinclozolin</i> Ukupno / In total 258				
	E5-3 Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroanilin, 2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen,	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Akrinatri (suma izomera), Alaklor, Aldrin, Aletrin -3,4 (Bioaletrin), Alidoklor, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Benfluralin, Benziladenin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Kaptan, Karbofenotion, Karbofenotion metil, Karfentrazon etil, Kinometionat, Klorbenzid, Klorbufam, Klordekon, Kloretoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorobenzilat, Kloroneb, Klortalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Klordan, cis-Permetrin, Klodinafop propargil, Klokintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklofluamid, Diklorvos, Diklobutrazol, Dikloran, Dikofol, Dieldrin, Dimefoks, Dimetaklor, Difenamid, Dimetipin, Dimetilan, Dinitramin, Difenilamin, Disulfoton, Ditalimfos, Edifenfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Famfur, Fenarimol, Fenklorfos, Fenflutrin, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flupiradifuron, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutolanil, Flutriafol, Fluvalinat-tau (suma izomera), Formotion, Furilazol, gama-BHC (lindan), Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izofenfos metil, Izoksaben, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metidation, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitalin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Paration-metil, Pebulat, Penkonazol, Pendi-metalin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Procimidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Protiokonazol, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piriproksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen, Siltiofam, Simazin, Spiromezifen, Sulkotrin, Sulfentrazon, Sulfotep, Sulprofos,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Swep, Tebukonazol, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tio-nazin, THPI (Tetrahidroftalimid), Tiokarbazil, Tolklofos-metil, Tolilfluaniid, Tolilfluaniid metabolit, trans-Klordan, Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Tridemorf, Trieta-zin, Trifluralin, Unikonazol, Valifenalat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allethrin-3,4 (Bioallethrin), Allidochlor, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Atraton, Azinphos-ethyl, Benfluralin, Benzyladenine, beta-BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Binapacryl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Captan, Carbophenothion, Carbophenothion methyl, Carfentrazone ethyl, Chinomethionat, Chlorbenside, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolate, Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlofluaniid, Dichlorvos, Diclobutrazol, Dicloran, Dicozol, Dieltrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Dinitramine, Diphenamid, Diphenylamine, Disulfoton, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Famphur, Fenarimol, Fenchlorphos, Fenfluthrin, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfathion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flupyradifurone, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutolanil, Flutriafol, Fluvalinate (sum of isomers), Formothion, Furilazole, gamma-BHC				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		(Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos methyl, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nutralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychloridan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Ofurace, Paclobutrazol, Parathion, Parathion-methyl, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phorate, Phosfolan., Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prosulfocarb, Prothioconazole, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesife, Sulcotrione, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbamil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vinclozolin Ukupno / In total 260				
	E5-4 Povrće i žitarice s visokim udjelom škroba i/ili proteina te niskim udjelom vode i masti	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (Z)-Pirifenox, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroaniline, 2,4-Dimetilanilin,	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	Vegetables and cereals - high starch and/or protein content and low water and fat content	2,4'-Metoksiklor , 2-Fenilfenol, 3,4'-Dikloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Aktrinatin (suma izomera), Alaklor, Aldrin, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Azinfos-metil, Benfluralin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitreks (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Kadusafos, Karbofenotion, Karbofenotion metil, Kinometionat, Klorbenzid, Klorbufam, Klordek, Klorotoksifos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorobenzilat, Kloroneb, Kloropropilat, Klortalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Klordan, cis-Nonaklor, cis-Permetrin, Klodinafop propargil, Klokvincet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklorvos, Diklobutrazol, Dikloran, Dieldrin, Dimetaklor, Dimeptin, Dimetilan, Difenilamin, Disulfoton, Ditalimfos, Edifenfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Famfur, Fenarimol, Fenklorfos, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotion, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetarin, Fluotrimazol, Flurokloridon, Flurprimidol, Flutiacet-metil, Fluvalinat-tau (suma izomera), Formotion, Furilazol, gama-BHC (lindan), Genit, Halfenproks, Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izo-karbofos, Izodrin, Izofenfos-metil, Izoksadifen-etil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metazaklor, Metakrifos, Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitrilin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksiklordan, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Paration-metil, Pebulat, Penkonazol, Pendimetalin, Pentakloroanilin, Pentakloroanisol, Pentaklorbenzen, Pentanoklor (Solan), Forat, Fosalon, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procimidon, Prodiamin, Profenofos, Profluralin, Prometrin, Propaklor, Propanil, Propazin, Propetamfos, Propizoklor, Propizamid, Prosulfokarb, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Piridafention, Pirimetanil, Piri-proksifen, Kvinalfos, Kvintozen, Resmetrin (suma izomera), S421, Silafluofen,				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Siltiofam, Simazin, Spiromezifen, Sulfentrazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimifos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Tiokarbazil, Tolklofos-metil, Tolilfluaniid, Tolilfluaniid metabolit, trans-Klordan, Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Triciklazol, Tridemorf, Trietazin, Trifluralin, Unikonazol, Vinklozolin / (E)-Metominstrobin, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6-Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 4,4'-Dichlorobenzophenone, 4,4'-Methoxychlor olefin,, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Azinphos methyl, Benfluralin, beta-BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Cadusafos, Carbophenothion, Carbophenothion methyl, Chinomethionat, Chlorbenside, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapyr, Chlorfenprop-methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chloropropylate, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal-dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolate, Cinidon ethyl, cis-Chlordan, cis-Nonachlor, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet-mexyl, Coumaphos, Crimidine, Crotoxyphos, Cyanazine, Cyanofenphos, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlorvos, Diclobutrazol, Dicloran, Dieldrin, Dimethachlor, Dimethipin, Dimetilan, Diphenamid, Diphenylamine, Disulfoton, Ditalimfos, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfathion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flurochloridon, Flurprimidol, Fluthiacet-methyl, Fluvalinate (sum of isomers), Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos-methyl, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methacrifos, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paclobutrazol, Parathion, Parathion-methyl, Pebulate, Penconazole, Pendimethalin, Pentanochlor (Solan), Phorate, Phosfolan., Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Prodiamine, Profenofos, Profluralin, Prometryn, Propachlor, Propanil, Propazin, Propetamphos, Propisochlor, Propyzamide, Prosulfocarb, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesifen, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbamil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transluthrin, trans-Nonachlor, trans-Permethrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Vinclozolin Ukupno / In total 257				
	E5-5 Teški matriksi <i>Difficult commodities</i>	Metoda za određivanje ostataka pesticida u hrani biljnog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg (E)-Metominostrobin, (E)-Pirifenoks, (Z)-Pirifenoks, 1,2,3,4-Tetraklorobenzen, 1,2,4,5-Tetraklorobenzen, 1,4-dimetilnaftalen, 2,3,5,6-Tetrakloroaniline,	GC-MS/MS	HRN EN 15662:2018 (EN 15662:2018)	23/07/2014	09/02/2023

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		2,4-Dimetilanilin, 2,4'-Metoksiklor, 2-Fenilfenol, 3,4-Dikloroanilin, 3-Kloroanilin, 4,4'-Diklorobenzofenone, 4,4'-Metoksiklor olefin, Acetoklor, Aklonifen, Aktrin (suma izomera), Alaklor, Aldrin, Alidoklor, alfa-BHC, alfa-Endosulfan, Amitraz, Antrakvinone, Azinfos-etil, Azinfos metil, Benfluralin, beta-BHC, beta-Endosulfan, Bifentrin, Bifenil, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-metil, Bromfenvinfos, Bromociklen, Bromofos, Bromofos-etil, Bromopropilat, Kadusafos, Karbofenotion, Karbofenotion metil, Karfentrazon etil, Kinometionat, Klorbenzid, Klorbufam, Klordekon, Kloretofos, Klorfenapir, Klorfenprop-metil, Klorfenson, Klormefos, Klorobenzilat, Kloroneb, Klortalonil, Klorprofam, Klorpirifos, Klorpirifos-metil, Klortal-dimetil, Klortion, Klortiofos (suma izomera), Klozolinat, Cinidon etil, cis-Kloran, cis-Permetrin, Klodinafop propargil, Klokintocet-meksil, Krimidin, Krotoksifos, Cianazin, Cianofenfos, Cianofos, Ciflutrin (suma izomera), Cipermetrin (suma izomera), Ciprodinil, Ciprofuram, delta-BHC, Deltametrin, Dialat (cis i trans), Diazinon, Diklobenil, Diklofention, Diklorvos, Diklobutrazol, Dikloran, Dieldrin, Dimefoks, Dimetaklor, Difenamid, Dimetipin, Dimetilan, Difenamid, Difenilamin, Disulfoton, Ditalimfos, Endosulfan sulfat, Endrin, EPN, EPTC, Etalfluralin, Etoprofos, Etoksikvin, Etofenproks, Etridiazol, Famfur, Fenarimol, Fenklorfos, Fenitrotrion, Fenoksaprop etil, Fenson, Fensulfotio, Fenvalerat (suma izomera), Fipronil, Flamprop izopropil, Flamprop metil, Fluazifop-P-butil, Flukloralin, Fludioksonil, Fluensulfon, Flumetrin, Fluotrimazol, Flurokloridon, Flurprimidol, Flusilazol, Flutiacet-metil, Flutriafol, Flualinat-tau (suma izomera), Formotion, Fostiazat (suma izomera), Furilazol, gama-BHC (lindan), Genit, Heptaklor, Heptaklor-egzo-epoksid, Heksaklorbenzen, Heksazinon, Imiprotrin, Jodofenfos, Ioksonil oktanoat, Iprobenfos, Iprodion, Izobenzan, Izokarbofos, Izodrin, Izofenfos metil, Izoksaben, Izoksadifenetil, Izoksation, lambda-Cihalotrin, Lenacil, Leptofos, Linuron, Mandestrobin, Metalaksil (Mefenoksam), Metidation, Metoksiklor, Metolaklor (S-Metolaklor), Mireks, Miklobutanil, N-(2,4-dimetilfenil) formamid, Nitrilin, Nitrofen, Nitrotal-izopropil, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurak, Oksadiazon, Oksikloran, Oksifluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paklobutrazol, Paration, Pebulat, Penkonazol, Pendimetilin, Pentakloroanilin, Pentakloroanisol, Pentaklorobenzen, Pentanoklor (Solan), Forat, Fosfolan, Fosmet, Ftalimid, Piperonil butoksid, Piperofos, Pirimifos etil, Pirimifos metil, Plifenat, Prokloraz, Procidon, Profluralin, Prometrin, Propaklor, Propanil, Propetamfos, Propizoklor, Propizamid, Prosulfokarb, Protiofos, Piraklofos, Piraflufen-etil, Piridaben, Pirdafention, Pirimetanil, Piripoksifen, Kvinalfos, Kvintozen, Resmetrin (suma				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		izomera), S421, Silafluofen, Siltiofam, Simazin, Spiromesifen, Sulkotrion, Sulfentazon, Sulfotep, Sulprofos, Swep, Tebukonazol, Tebufenpirad, Tebupirimi- fos, Tebutam, Teknazen, Teflutrin, Terbacil, Terbufos, Terbutilazin, Tetradifon, Tetrametrin (suma izomera), Tetrasul, Tionazin, THPI (Tetrahidroftalimid), Ti- okarbazil, Tolklofos-metil, Tolilfluaniid, Tolilfluaniid metabolit, trans-Klordan, Transflutrin, trans-Permetrin, Triadimefon, Triadimenol, Trialat, Triamifos, Triazofos, Trikloronat, Tridemorf, Trietazin, Trifluralin, Unikonazol, Valifena- lat, Vinklozolin / (E)-Metominstrobin, (E)-Pyrifenox, (Z)-Pyrifenox, 1,2,3,4-Tetrachlorbenzene, 1,2,4,5-Tetrachlorbenzene, 1,4-dimethylnaphthalene, 2,3,5,6- Tetrachloroaniline, 2,4-Dimethylaniline, 2,4'-Methoxychlor, 2-Phenylphenol, 3,4-Dichloroaniline, 3-Chloroaniline, 4,4'-Dichlorobenzophenone, 4,4'- Methoxychlor olefin,, Acetochlor, Aclonifen, Acrinathrin (sum of isomers), Alachlor, Aldrin, Allidochlor, alpha-BHC, alpha-Endosulfan, Amitraz, Anthraquinone, Azinphos-ethyl, Azinphos methyl, Benfluralin, beta- BHC, beta-Endosulfan, Bifenthrin, Biphenyl, Bitrex (Denatonium benzoat), Bromacil, Bromfenvinfos-methyl, Bromfenvinphos, Bromocyclen, Bromophos, Bromophos-ethyl, Bromopropylate, Cadusafos, Carbophenothion, Carbophenothion methyl, Carfentrazone-ethyl, Chinomethionat, Chlorbenside, Chlorbufam, Chlordecone, Chlorethoxyfos, Chlorfenapir, Chlorfenprop- methyl, Chlorfenson, Chlormephos, Chlorobenzilate, Chloroneb, Chlorothalonil, Chlorpropham, Chlorpyrifos, Chlorpyrifos-methyl, Chlorthal- dimethyl, Chlorthion, Chlorthiophos (sum of isomers), Chlozolate, Cinidon ethyl, cis-Chlordan, Cis-Permethrin, Clodinafop-propargyl, Cloquintocet- mexyl, Crimidine, Crotoxyphos, Cyanazine, Cyanofenphos, Cyanophos, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Cyprodinil, Cyprobufam, delta-BHC, Deltamethrin, Diallylate (cis & trans), Diazinon, Dichlobenil, Dichlofenthion, Dichlorvos, Diclobutrazol, Dicloran, Dieldrin, Dimefox, Dimethachlor, Dimethipin, Dimetilan, Diphenamid, Diphenylamine, Disulfoton, Edifenphos, Endosulfan sulfate, Endrin, EPN, EPTC, Ethalfluralin, Ethoprophos, Ethoxyquin, Etofenprox, Etridiazole, Famphur, Fenarimol, Fenchlorphos, Fenitrothion, Fenoxaprop-ethyl, Fenson, Fensulfothion, Fenvalerate (sum of isomers), Fipronil, Flamprop-isopropyl, Flamprop-methyl, Fluazifop-P-butyl, Fluchloralin, Fludioxonil, Fluensulfone, Flumethrin, Fluotrimazol, Flurochloridon, Flurprimidol, Flusilazole, Fluthiacet-methyl, Flutriafol, Fluvalinate (sum of isomers), Formothion, Fosthiazate (sum of				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		isomers), Furilazole, gamma-BHC (Lindane), Genite, Halfenprox, Heptachlor, Heptachlor-exo-epoxide, Hexachlorobenzene, Hexazinone, Imiprothrin, Iodofenphos, Ioxynil octanoate, Iprobenfos, Iprodione, Isobenzan, Isocarbophos, Isodrin, Isofenphos-methyl, Isoxaben, Isoxadifen ethyl, Isoxathion, lambda-Cyhalothrin, Lenacil, Leptophos, Linuron, Mandestrobin, Metalaxyl (Mefenoxam), Metazachlor, Methidathion, Methoxychlor, Metolachlor (S-Metolachlor), Mirex, Myclobutanil, N-(2,4-dimethylphenyl) formamide, Nitralin, Nitrofen, Nitrothal-isopropyl, o,p'-DDD, o,p'-DDE, o,p'-DDT, Ofurace, Oxadiazon, Oxychlordan, Oxyfluorfen, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paclobutrazol, Parathion, Pebulate, Penconazole, Pendimethalin, Pentachloroaniline, Pentachloroanisole, Pentachlorobenzene, Pentanochlor (Solan), Phorate, Phosfolan., Phosmet, Phthalimide, Piperonyl butoxide, Piperophos, Pirimiphos ethyl, Pirimiphos-methyl, Plifenate, Prochloraz, Procymidone, Profluralin, Prometryn, Propachlor, Propanil, Propetamphos, Propisochlor, Propyzamide, Prosulfocarb, Prothiofos, Pyraclofos, Pyraflufen-ethyl, Pyridaben, Pyridaphenthion, Pyrimethanil, Pyriproxyfen, Quinalphos, Quintozene, Resmethrin (sum of isomers), S421, Silafluofen, Silthiofam, Simazine, Spiromesifen, Sulfentrazone, Sulfotep, Sulprofos, Swep, Tebuconazole, Tebufenpyrad, Tebupirimifos, Tebutam, Tecnazene, Tefluthrin, Terbacil, Terbufos, Terbutylazine, Tetradifon, Tetramethrin (sum of isomers), Tetrasul, Thionazin, THPI (Tetrahydrophthalimide), Tiocarbazil, Tolclofos-methyl, Tolyfluanid, Tolyfluanid metabolite, trans-Chlordane, Transfluthrin, trans-Permethrin, Triadimefon, Triadimenol-1, Tri-allate, Triamiphos, Triazophos, Tribufos, Trichloronat, Tridemorph, Trietazin, Trifluralin, Uniconazole, Valifenalate, Vinclozolin Ukupno / In total 254				
	E5-6 Meso i mesni proizvodi <i>Meat and meat</i> <i>products</i>	Metoda za određivanje ostataka pesticida u hrani životinjskog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of animal origin –</i> <i>gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,001 mg/kg	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-1- 2,16/03 6. izdanje / <i>edition</i> (2026-04-09)	24/01/2005	09/04/2026

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		4,4'-diklorobenzofenon, 4,4'-metoksiklor olefin, aldrin, alfa-BHC (alfa-heksaklorocikloheksan), alfa-endosulfan (Endosulfan I), azinfos etil, beta-BHC (beta-hheksaklorocikloheksan), beta-endosulfan (Endosulfan II), bifentrin, klorbenzid, klorfenskon, kjlorfenvinfos (suma izomera), klorobenzilat, kloroneb, klorpirifos, klorpirifos-metil, cis-klordan, cis-nonaklor, cis-permetrin, ciflutrin (suma izomera), cipermetrin (suma izomera), delta-BHC (delta-heksaklorocikloheksan), deltametrin, diazinon, dieldrin, endosulfan ether, endosulfan sulfate, endrin, endrin keton, famoksadon, fenson, fention, fention sulfon, fenvalerate (suma izomera), fipronil, gama-BHC (Lindan), heptaklor, heptaklor-endo-epoksid, heptachlor-egzo-epoksid, heksaklorbenzen, indoksakarb, izodrin, lambda-cihalotrin, malation, metidation, metoksiklor, mireks, o,p'-DDD, o,p'-DDE, o,p'-DDT, Oksiklordan, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paration, Paration metil, pendimetalin, pentakloroanisol, pentaklorbenzen, pentaklorotioanisol, fosmet, pirimifos metil, profenofos, pirazofos, kintozen, teknazen, tetradifon, trans-klordan, trans-nonaklor, trans-permetrin, triazofos, vinklozolin / 4,4'-Dichlorobenzophenone, 4,4'-methoxychlor olefin, Aldrin, alpha-BHC (alpha-hexachlorocyclohexane), alpha-Endosulfan (Endosulfan I), Azinphos ethyl, beta-BHC (beta-hexachlorocyclohexane), beta-Endosulfan (Endosulfan II), Bifenthrin, Chlorbenside, Chlorfenson, Chlorfenvinphos (sum of isomers), Chlorobenzilate, Chloroneb, Chlorpyrifos, Chlorpyrifos-methyl, cis-Chlordane, cis-Nonachlor, cis-Permethrin, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), delta-BHC (delta-hexachlorocyclohexane), Deltamethrin, Diazinon, Dieldrin, Endosulfan ether, Endosulfan sulfate, Endrin, Endrin ketone, Famoxadone, Fenson, Fenthion, Fenthion sulfone, Fenvalerate (sum of isomers), Fipronil, gamma-BHC (Lindane), Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Indoxacarb, Isodrin, lambda-Cyhalothrin, Malathion, Methidation, Methoxychlor, Mirex, o,p'-DDD, o,p'-DDE, o,p'-DDT, Oxychlordan, p,p'-DDD, p,p'-DDE, p,p'-DDT, Parathion, Parathion methyl, Pendimethalin, Pentachloroanisol, Pentachlorobenzene, Pentachlorothioanisol, Phosmet, Pirimiphos methyl, Profenofos, Pyrazophos, Quintozene, Tecnazene, Tetradifon, trans-Chlordane, trans-Nonachlor, trans-Permethrin, Triazophos, Vinclozolin Ukupno / In total 72				

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	E5-7 Mlijeko i mliječni proizvodi, sladoled i puding <i>Milk and dairy products, icecream and pudding</i>	Metoda za određivanje ostataka pesticida u hrani životinjskog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of animal origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,0008 mg/kg Endrin / <i>Endrin</i> 0,002 mg/kg Klorfenvinfos (suma izomera), 2,4'-Metoksiklor, 4,4'-Diklorobenzofenon, 4,4'-Metoksiklor olefin, Aldrin, alfa-BHC (alfa-heksaklorocikloheksan), alfa-Endosulfan (Endosulfan I), Azinfos etil, Azinfos metil, beta-BHC (beta-heksaklorocikloheksan), beta-Endosulfan (Endosulfan II), Bifentrin, Klorbenzid, Klorfenson, Klorobenzilat, Kloroneb, Klorpirifos, Klorpirifos-metil, cis-Klordan, cis-Nonaklor, cis-Permetrin, Ciflutrin (suma izomera), Cypermethrin (suma izomera), delta-BHC (delta-heksaklorocikloheksan), Deltametrin, Diazinon, Dieldrin, Endosulfan eter, Endosulfan sulfat, Endrin aldehid, Endrin keton, Famoksadon, Fenson, Fention, Fention okson, Fention sulfon, Fention sulfoksid, Fenvalerat (suma izomera), Fipronil, gama-BHC (Lindan), Heptaklor, Heptaklor-endo-epoksid, Heptaklor-egzo-epoksid, Heksaklorbenzen, Indoksakarb, Izodrin, lambda-Cihalotrin, Malation, Metidation, Metoksiklor, Mireks, o,p'-DDD, o,p'-DDE, o,p'-DDT, Oksiklordan, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paration, Paration metil, Pendimetalin, Pentakloroanisol, Pentaklorbenzen, Pentaklorotioanisol, Fosmet, Pirimifos metil, Profenofos, Pirazofos, Kintozen, Remetrin (suma izomera), Teknazen, Tetradifon, trans-Klordan, trans-Nonaklor, trans-Permetrin, Triazofos, Vinklozolin / <i>Chlorfenvinphos (sum of isomers), 2,4'-Methoxychlor, 4,4'-Dichlorobenzophenone, 4,4'-methoxychlor olefin, Aldrin, alpha-BHC (alpha-hexachlorocyclohexane), alpha-Endosulfan (Endosulfan I), Azinphos ethyl,</i>	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-1- 2,16/03 6. izdanje / <i>edition</i> (2026-04-09)	01/12/2004	09/04/2026

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		Azinphos methyl, beta-BHC (beta-hexachlorocyclohexane), beta-Endosulfan (Endosulfan II), Bifenthrin, Chlorbenseide, Chlorfenson, Chlorobenzilate, Chloroneb, Chlorpyrifos, Chlorpyrifos-methyl, cis-Chlordane, cis-Nonachlor, cis-Permethrin, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), delta-BHC (delta-hexachlorocyclohexane), Deltamethrin, Diazinon, Dieldrin, Endosulfan ether, Endosulfan sulfate, Endrin aldehyde, Endrin ketone, Famoxadone, Fenson, Fenthion, Fenthion oxon, Fenthion sulfone, Fenthion sulfoxide, Fenvalerate (sum of isomers, Fipronil, gamma-BHC (Lindane), Heptachlor, Heptachlor-endo-epoxide, Heptachlor-exo-epoxide, Hexachlorobenzene, Indoxacarb, Isodrin, lambda-Cyhalothrin, Malathion, Methidation, Methoxychlor, Mirex, o,p'-DDD, o,p'-DDE, o,p'-DDT, Oxychlordan, p,p'-DDD, p,p'-DDE, p,p'-DDT, Parathion, Parathion methyl, Pendimethalin, Pentachloroanisole, Pentachlorobenzene, Pentachlorothioanisole, Phosmet, Pirimiphos methyl, Profenofos, Pyrazophos, Quintozene, Resmethrin (sum of isomers), Tecnazene, Tetradifon, trans-Chlordane, trans-Nonachlor, trans-Permethrin, Triazophos, Vinclozolin Ukupno / In total 78				
	E5-8 Jaja Eggs	Metoda za određivanje ostataka pesticida u hrani životinjskog podrijetla metodom plinske kromatografije s masenom spektrometrijom (GC-MS/MS) <i>Method for the determination of pesticide residues in foods of animal origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,005 mg/kg 2,4'-Metoksiklor, 4,4'-Diklorobenzofenon, 4,4'-metoksiklor olefin, Aldrin, alfa-BHC (alfa-heksaklorocikloheksan), alfa-Endosulfan (Endosulfan I), Azinfos etil, Azinfos metil, beta-BHC (beta-heksaklorocikloheksan), beta-Endosulfan (Endosulfan II), Bifenthrin, Ciflutrin (suma izomera), Cipermetrin (suma izomera), cis-Klordan, cis-Nonaklor, cis-Permetrin, delta-BHC (delta-heksaklorocikloheksan), Deltametrin, Dijazinon, Dieldrin, Endosulfan eter, Endosulfan sulfat, Endrin, Endrin aldehid, Endrin keton, Famoksadon, Fenson, Fention, Fenthion oxon, Fention sulfon, Fenvalerat (suma izomera), Fipronil, Fosmet, gama-BHC (Lindan), Heksaklorobenzen, Heptaklor, Heptaklor-egzo-	GC-MS/MS	Vlastita metoda In-house method SOP-LKH-1- 2,16/03 6. izdanje / edition (2026-04-09)	15/06/2021	09/04/2026

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		epoksd (izomer B), Heptaklor-endo-epoksid (izomer A), Indoksakarb, Izodrin, Klorbenzid, Klorfenson, Klorfenvinfos (suma izomera), Klorobenzilat, Kloroneb, Klorpirifos, Klorpirifos-metil, Kintozen, lambda-Cihalotrin, Malation, Metidation, Metoksiklor, Mireks, o,p'-DDD, o,p'-DDE, o,p'-DDT, Oksiklordan, p,p'-DDD, p,p'-DDE, p,p'-DDT, Paration, Paration metil, Pendimetalin, Pentakloroanisol, Pentaklorobenzen,, Pentaklorotioanisol, Pirazofos, Pirimifos-metil, Profenofos, Resmetrin (suma izomera), Teknazen, Tetradifon, trans-Klordan, trans-Nonaklor, trans-Permetrin, Triazofos, Vinklozolin / 2,4'-Methoxychlor, 4,4'-Dichlorobenzophenone, 4,4'-methoxychlor olefin, Aldrin, alpha-BHC (alpha-hexachlorocyclohexane), alpha-Endosulfan (Endosulfan I), Azinphos ethyl, Azinphos methyl, beta-BHC (beta-hexachlorocyclohexane), beta-Endosulfan (Endosulfan II), Bifenthrin, Cyfluthrin (sum of isomers), Cypermethrin (sum of isomers), Chlorbenside, Chlorfenson, Chlorfenvinfos (sum of isomers), Chlorobenzilate, Chloroneb, Chlorpyrifos, Chlorpyrifos-methyl, cis-Chlordane, cis-Nonachlor, cis-Permethrin, delta-BHC (delta-hexachlorocyclohexane), Deltamethrin, Dieldrin, Endosulfan ether, Endosulfan sulfate, Endrin aldehyde, Endrin ketone, Endrin, Famoxadone, Fenson, Fenthion, Fenthion oxone, Fenthion sulfone, Fenvalerate (sum of isomers), Fipronil, gamma-BHC (Lindane), Heptachlor, Heptachlor-exo-epoxide (isomer B), Heptachlor-endo-epoxide (A), Hexachlorobenzene, Indoxacarb, Isodrin, lambda-Cyhalothrin, Malathion, Methidation, Methoxychlor, Mirex, o,p'-DDD, o,p'-DDE, o,p'-DDT, Oxychlordan, p,p'-DDD, p,p'-DDE, p,p'-DDT, Parathion, Parathion methyl, Pendimethalin, Pentachloroanisolet, Pentachlorobenzene, Pentachlorothioanisole, Phosmet, Pirimiphos methyl, Profenofos, Pyrazophos, Quintozene, Diazinon, Resmethrin (sum of isomers), Tecnazene, Tetradifon, trans-Chlordane, trans-Nonachlor, trans-Permethrin, Triazophos, Vinclozolin Ukupno / In total 77				
E6	Hrana biljnog podrijetla <i>Food of plant origin</i>	Određivanje ostataka pesticida u hrani pojedinačnim metodama (Single Residue Method - SRM)	LC-MS/MS			

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
		<i>Determination of pesticide residues in food - Single Residue Methods (SRM)</i>				
	E6-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje matrina u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Matrine in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,005 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/229 II. izdanje / <i>edition</i> (2021-06-10)	15/09/2020	10/06/2021
	E6-2 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water</i>	Određivanje maleinskog hidrazida u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Maleic-hydrazide in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,05 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/234 I. izdanje / <i>edition</i> (2021-04-20)	20/04/2021	20/04/2021
	E6-3 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom	Određivanje ditianona u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Dithianon by liquid chromatography-mass spectrometry (LC-MS/MS)</i>	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/220 II. izdanje / <i>edition</i> (2021-11-03)	11/09/2019	03/11/2021

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i> Povrće i žitarice s visokim udjelom škroba i/ili proteina te niskim udjelom vode i masti <i>Vegetables and cereals - high starch and/or protein content and low water and fat content</i>	0,01 mg/kg				
	E6-4 Voće i povrće s visokim udjelom	Određivanje klormekvata i mepikvata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom	LC-MS/MS	Vlastita metoda <i>In-house method</i>	21/11/2017	21/11/2017

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	vode <i>Fruits and vegetables – high water content</i>	spektrometrijom (LC-MS/MS) <i>Determination of Chlormequat and Mepiquat by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg		SOP-LKH-41/189 I. izdanje / edition (2017-11-21) QuPPE-PO, V12		
	E6-5 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Određivanje etefona u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Ethephon in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/210 II. izdanje / edition (2020-04-15) QuPPE-PO, V12	15/04/2020	15/04/2020
	E6-6 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje glifosata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Glyphosate by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/221 I. izdanje / edition (2019-07-16) QuPPE-PO, V12	16/07/2019	16/07/2019

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	E6-7 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>	Određivanje fosetila i fosfonske kiseline QuPPe metodom u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Fosetyl and Phosphonic acid following QuPPe method in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg Fosetil/Fosetyl 0,02 mg/kg Fosfonska kiselina/Phosphonic acid	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/236 I. izdanje / <i>edition</i> (2021-12-02) QuPPe-PO, V12	02/12/2021	02/12/2021
	E6-8 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje morfolina, dietanolamina i trietanolamina metodom tekućinske kromatografije s masenom spektrometrijom LC-MS/MS <i>Determination of morpholine, diethanolamine and triethanolamine by liquid chromatography with LC-MS / MS mass spectrometry</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/239 I. izdanje / <i>edition</i> (2022-02-28)	28/02/2022	28/02/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i>					
	E6-9 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Određivanje klorata i perklorata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Chlorate and Perchlorate in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/246 I. izdanje / <i>edition</i> (2022-12-02) <i>QuPPe-PO, V12</i>	02/12/2022	02/12/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	E6-10 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i>	Određivanje parakvata i dikvata u hrani biljnog podrijetla metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of Paraquat and Diquat in food of plant origin by liquid chromatography-mass spectrometry (LC-MS/MS)</i> 0,01 mg/kg	LC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/245 I. izdanje / <i>edition</i> (2022-12-02) <i>QuPPE-PO, VI2</i>	02/12/2022	02/12/2022
E7	Hrana biljnog podrijetla <i>Food of plant origin</i>	Određivanje ostataka pesticida u hrani pojedinačnim metodama (Single Residue Method - SRM) <i>Determination of pesticide residues in food - Single Residue Methods (SRM)</i>	GC-MS/MS			
	E7-1 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i>	Određivanje captana, folpeta i njihovih metabolita u hrani biljnog podrijetla GC-MS/MS metodom <i>Determination of Captan, Folpet and their metabolites in food of plant origin – gas chromatographic method with mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg Captan, Folpet, THPI, Ftalimid / Captan, Folpet, THPI, Phthalimide	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/248 I. izdanje / <i>edition</i> (2022-11-16)	16/11/2022	16/11/2022
	E7-2 Voće i povrće s visokim udjelom					

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i> Povrće i žitarice s visokim udjelom škroba i/ili proteina te niskim udjelom vode i masti <i>Vegetables and cereals - high starch and/or protein content and low water and fat content</i>	Određivanje ostataka fenbutatin oksida u hrani biljnog podrijetla GC-MS/MS tehnikom <i>Determination of Fenbutatin oxide in food of plant origin by gas chromatography - mass spectrometry (GC-MS/MS)</i> 0,002 mg/kg Fenbutatin oksid / Fenbutain oxide	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/201 III. izdanje / <i>edition</i> (2022-02-16)	13/08/2018	16/02/2022

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	E7-3 Voće i povrće s visokim udjelom vode <i>Fruits and vegetables – high water content</i> Voće i povrće s visokim udjelom kiseline i visokim udjelom vode <i>Fruits and vegetables – high acid content and high water content</i> Voće i povrće s visokim udjelom ulja i srednjim te niskim udjelom vode <i>Fruits and vegetables – high oil content and intermediate or low water content</i> Povrće i žitarice s visokim udjelom škroba i/ili proteina te niskim udjelom vode i masti <i>Vegetables and</i>	Određivanje amitraza u hrani bilnog podrijetla GC-MS/MS tehnikom <i>Determination of Amitraz in food of plant origin by gas chromatography - mass spectrometry (GC-MS/MS)</i> 0,01 mg/kg Amitraz / Amitraz	GC-MS/MS	Vlastita metoda <i>In-house method</i> SOP-LKH-41/206 III. izdanje / <i>edition</i> (2021-09/12)	30/10/2018	12/09/2021

Oznaka/ Identification	Materijali/ Proizvodi Materials/ Products	Vrsta ispitivanja/Svojstvo Type of test / Property Granica kvantifikacije / Limit of quantification	Tehnika / Technique	Metoda ispitivanja Test method	U primjeni od/ In use from	Datum posljednje izmjene/ Date of last change
	cereals - high starch and/or protein content and low water and fat content					
E8	E8-1 Hrana biljnog podrijetla (uključujući mlinske i pekarske proizvode, tjesteninu i brzo smrznuta tijesta) <i>Food of plant origin (including mill and bakery products, pasta and quick- frozen dough)</i>	Određivanje toksina metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of toxins - liquid chromatography with mass spectrometry (LC- MS/MS)</i> 1 µg/kg Aflatoksin B1, B2, G1, G2, HT2- toksin, T2-toksin, ohratoksin A/aflatoxin B1, B2, G1, G2, HT2-toxin, T2-toxin, ochratoxin 10 µg/kg Zearalenon/zearalenone 50 µg/kg Deoksinivalenol/deoxynivalenol	LC-MS/MS	Vlastita metoda In-house method SOP-LKH- 02,41/203 VI. izdanje / edition (2024-07-12)	07/09/2018	12/7/2024
	E8-2 Mlijeko i mliječni proizvodi <i>Milk</i>	Određivanje toksina metodom tekućinske kromatografije s masenom spektrometrijom (LC-MS/MS) <i>Determination of toxins - liquid chromatography with mass spectrometry (LC- MS/MS)</i> 0,02 µg/kg Aflatoksin M1/aflatoxin M1	LC-MS/MS	Vlastita metoda In-house method SOP-LKH- 02,41/203 VI. izdanje / edition (2024-07-12)	21/02/2022	12/7/2024
E9	Hrana <i>Food</i>	Određivanje odabranih elemenata metodom atomske apsorpcijske spektrometrije (AAS) <i>Determination of selected elements - atomic absorption spectrometry (AAS)</i> 0,005 mg/kg Pd 0,004 mg/kg Cd	AAS	HRN EN 14083:2005 (EN 14083:2003)	21/4/2005	27/02/2025

Priprema analitičkih uzorka / Preparation of analytical samples:

- Tlo / *Soil*: modificirana / *modified* HRN EN ISO 54321:2021 (*ISO 54321:2020; EN ISO 54321:2021*); HRN ISO 14869-2:2004 (*ISO 14869-2:2002*)
- Otpad / *Waste*: modificirana / *modified* HRN EN 13657:2008 (*EN 13657:2002*); modificirana / *modified* HRN EN 12457-4:2005 (*EN 12457-4:2002*); HRN ISO 14869-2:2004 (*ISO 14869-2:2002*)
- Eluat otpada / *Waste leachate* HRI CEN/TR 16192:2020
- Muljevi, sedimenti, biota / *Sludges, sediments, biota*: modificirana / *modified* HRN EN ISO 54321:2021 (*ISO 54321:2020; EN ISO 54321:2021*)
- Kruta oporabljena goriva / *Solid recovered fuels*: [HRN EN ISO 3884:2025](#) (*EN ISO 3884:2025*)
- Hrana / *Food*: HRN EN 13805:2014 (*EN 13805:2014*)
- Hrana biljnog podrijetla / *Food of plant origin*: HRN EN 15662:2018 (*EN 15662:2018*)
- Hrana biljnog podrijetla / *Food of plant origin*: HRN EN 15662:2018 modificirana sa alkalnom hidrolizom / *modified (EN 15662:2018) with alkaline hydrolysis*
- Hrana biljnog podrijetla / *Food of plant origin*: QuPPE-PO-Method

Napomene / Notes:

ICP-OES - optička emisijska spektrometrija induktivno vezane plazme / ICP OES - *inductively coupled plasma atomic emission spectroscopy*

AAS – atomska apsorpcijska spektrometrija / AAS – *atomic absorption spectroscopy*

AAS (FIMS) – “Flow Injecton Mercury System” atomska apsorpcijska spektrometrija / AAS (FIMS) – *flow injecton mercury system atomic absorption spectroscopy*

ICP-MS - masena spektrometrija induktivno vezane plazme / ICP - MS - *inductively coupled plasma / mass spectrometry*

UV-Vis - spektrofotometrija / *spectrophotometry*

GC-MS - plinska kromatografija / spektrometrija masa / GC-MS - *Gas chromatography - mass spectrometry*

GC-MS/MS – plinska kromatografija / tandemna spektrometrija masa / GC-MS/MS - *Gas chromatography-tandem mass spectrometry*

LC-MS/MS – tekućinska kromatografija / tandemna spektrometrija masa / LC-MS/MS – *Liquid chromatography-tandem mass spectrometry*

Početak fleksibilnog područja akreditacije 2012.godina. / The beginning of flexible scope of accreditation 2012.**Praćenje izmjena**

Datum izmjene	Opis izmjene
22/5/2025	Novo izdanje usklađeno s Potvrdom o akreditaciji (2025-05-22)
15/7/2025	Dopuna metode E7: Određivanje ostataka pesticida u hrani pojedinačnim metodama (Single Residue Method - SRM)
19/01/2026	Izmjene u metodama: A1, B1, B1-2, C4 i C8
20/01/2026	Izmjena u metodi: B1
10/6/2026	Izmjene u metodama: A3, B1, B2-1, C5, E2, E5-6, E5-7, E5-8
06/7/2026	Izmjene u metodi: E1
24/8/2026	Izmjene u metodi: A1, A2, A3, B1, B3-1, B3-2